To be used only under the direction of a physician - Non-Sterile .
Not made with natural rubber latex - For Single Patient Use Only

Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.

Hand For use only with Flowtron® DVT-prevention pumps manufactured by ArjoHuntleigh. Must not be used with Flowtron® Hydroven 3 or Flowtron® Hydroven 12 Intermittent Pneumatic Compression (IPC) pumps.

Hand More comprehensive information on the garment can be found in the relevant Flowtron DVT pump Instructions For Use document.

Instructions For Use

- Plug the pump into a suitable electrical outlet. Do not turn the pump on at this time.
- Remove the garments from the sealed bag. The garments may be used on either leg. Unfold the garment and position the inflatable bladder directly beneath the patient's skin.
- Snugly wrap the garment around the patient's leg and secure the fastener tabs. Repeat for the other leg. Garments should be positioned in such a way that they do not create any potential for constriction pressure points on the patient's limb. Use appropriate wrist straps or securing devices to prevent inflation lines from being pulled or secured inside the garment next to patient skin, and regularly check the patient's skin for signs of redness or pressure points. The garment is most effective in preventing venous stasis when the garment air bladder is located in the posterior position. If the garment cannot be placed at the posterior position, which will still help to prevent venous stasis.
- Attach the garments to the pump tubing set ensuring a "click" is heard from each snap-lock connector.
- If your Flowtron DVT pump has an adjustable pressure regulator, turn the regulator dial to the recommended pressure unless otherwise directed by the physician. For further details, refer to the relevant Flowtron DVT pump Instructions For Use document.
- Turn the pump on. The green power indicator lights should illuminate.

Indications
Intermittent Pneumatic Compression (IPC) is indicated for use for the prevention of deep vein thrombosis due to the presence of risk factors for thrombus formation during orthopaedic, trauma, urologic and general surgery. It is also suitable for other patient groups including neurology, critical care, general medical and obstetrics.

© and ™ are trademarks belonging to the ArjoHuntleigh group of companies.
As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

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