# Flowtron

## DVT

Calf Garment • Attele Surale • Wadenmanschette • Kuit Manchét • Fascia per Polpaccio  
Prenda Para Pantorrilla • Vadmanschett • Læg Forbindingen • Ενδύμα γαστροκνημίας

**EN** To be used only under the direction of a physician - Non-Sterile -  
Not made with natural rubber latex - For Single Patient Only

**Caution**  
Federal law restricts this device to sale by or on the order of a licensed practitioner.

**ES** For use only with Flowtron® DVT-prevention pumps manufactured by ArjoHuntleigh. Must not be used with Flowtron® Hydropren 3 or Flowtron® Hydroyen 12 Intermittent Pneumatic Compression (IPC) pumps.

**ES** More comprehensive information on the garment can be found in the relevant Flowtron DVT pump Instructions For Use document.

#### Instructions for Use

1. Plug the pump into a suitable electrical outlet. Do not turn the pump on at this time.
2. Remove the garments from the sealed bag. The garments may be used on either leg. Unfold the garment and position the inflatable bladder directly behind the patient's calf.
3. Snugly wrap the garment around the patient's leg and secure the fastener tabs. Repeat for the other leg. Garments should be positioned in such a way that they do not create any potential for constant pressure points on the patient's limb. If using apparatus with straps or securing devices - for example lithotomy stirrups, ensure tubing is not placed inside the strap next to patient skin, and regularly check the patient's skin for signs of redness or pressure points. The garment is most effective in preventing venous stasis when the garment air bladder is located in the posterior position. If the garment cannot be placed at the posterior position, the pump should be used on the calf at alternative positions all of which will still help to prevent venous stasis.
4. Attach the garments to the pump tubing set ensuring a 'click' is heard from each snap-lock connector.
5. If using Flowtron DVT pump has an adjustable pressure regulator, turn the regulator dial to the recommended pressure unless otherwise directed by the physician. For further details, refer to the relevant Flowtron DVT pump Instructions For Use document.
6. Turn the pump on. The green power indicator lights should illuminate.

#### Indications

Intermittent Pneumatic Compression (IPC) is indicated for use for the prevention of deep vein thrombosis due to the presence of risk factors for thrombus formation during orthopaedic, trauma, urologic and general surgery. It is also suitable for other patient groups including neurology, critical care, general medical and obstetrics.

© and ™ are trademarks belonging to the ArjoHuntleigh group of companies  
As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

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#### Recommendations

- Garments should be removed regularly to inspect the skin for signs of redness or pressure points
- The garments should be applied to the patient pre-operatively, prior to the induction of anaesthesia
- The system should be used continuously for no less than 72 hours post-operatively or until the patient becomes fully ambulatory
- In the non-surgical patient, the system should be initiated immediately once the risk of DVT formation is identified

#### Contraindications

- IPC should not be used in the following conditions:
1. Severe arteriosclerosis or other ischaemic vascular diseases.
  2. Known or suspected acute Deep Vein Thrombosis (CDVT) or phlebitis.
  3. Severe congestive cardiac failure or any condition where an increase of fluid to the heart may be detrimental.
  4. Pulmonary embolism.
  5. Any local condition in which the garments would interfere, including gangrene, recent skin graft, dermatitis or untreated, infected leg wounds.

If you are unsure of any contraindications refer to the patient's physician before using the device.

#### Cautions

1. Garments should be removed immediately if the patient experiences itching, numbness or pain.
2. When used for DVT prophylaxis, continuous use is recommended and any interruption of therapy for a substantial length of time should be at the discretion of the physician.
3. Lower limb positioning in relation to the garment and tubing should also be considered particularly when a patient is unconscious, or has reduced sensation / ability to move their legs.
4. Additional care should be taken when placing the garments on any deformed leg or foot, or on legs with significant oedema.
5. This product cannot be adequately cleaned and / or sterilized by the user in order to facilitate safe reuse and is therefore intended for single patient use. Attempts to clean or sterilize these devices may result in a biocompatibility, infection or product failure risk to the patient.

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ArjoHuntleigh Inc.
Oy Vestek AB
ArjoHuntleigh AB