71cm (28")

EN To be used only under the direction of a physician - Non-Sterile
- 100% Latex Free - For Single Patient Use Only

Caution
Federal law restricts this device to sale by or on the order of a licensed practitioner.

ES Use for only with **Flowtron® DVT**-prevention pumps manufactured by ArjoHuntleigh or Huntleigh Healthcare. *Must not be used with Flowtron® Hydroven 3 or Flowtron® Hydroven 12 Intermittent Pneumatic Compression (IPC) pumps.*

DE More comprehensive information on the garment can be found in the relevant **Flowtron DVT pump Instructions For Use** document.

- Inструкciones para Uso**
- Plug the pump into a suitable electrical outlet. Do not turn the pump on at this time.
 - Remove the garments from the sealed bag. The garments may be used on either leg. Unfold the garment and position the inflatable bladder around the patient's calf.
 - Snugly wrap the garment around the patient's leg and secure the fastener tabs. Repeat for the other leg.
 - Attach the garments to the pump tubing set ensuring a "click" is heard from each snap-lock connector.
 - If your **Flowtron DVT** pump has an adjustable pressure regulator, turn the regulator dial to the recommended pressure unless otherwise directed by the physician. For further details, refer to the relevant **Flowtron DVT pump Instructions For Use** document.
 - Turn the pump on. The green power indicator lights should illuminate.

Indications
Intermittent Pneumatic Compression (IPC) is indicated for use for the prevention of deep vein thrombosis due to the presence of risk factors for thrombosis formation during orthopaedic, trauma, urologic and general surgery. It is also suitable for other patient groups including neurology, critical care, general medical and obstetrics.

Recommendations
General recommendations:

- Garments should be removed regularly to inspect the skin
- Where appropriate, patients should be instructed in the proper use of the system, the purpose of therapy and any problems that should be reported to the nursing staff

© and ™ are trademarks belonging to the ArjoHuntleigh group of companies.
As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

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DVT prophylaxis:

- The garments should be applied to the patient pre-operatively, prior to the induction of anaesthesia
- The system should be used continuously for no less than 72 hours post-operatively or until the patient becomes fully ambulatory
- If the garment cannot be applied to the operative limb during surgery, it may be applied to the limb once the patient reaches the recovery unit
- In the non-surgical patient, the system should be initiated immediately once the risk of DVT formation is identified

Contraindications

- IPC should not be used in the following conditions:
- Severe arteriosclerosis or other ischaemic vascular diseases.
 - Known or suspected acute Deep Vein Thrombosis (DVT) or phlebitis.
 - Severe congestive cardiac failure or any condition where an increase of fluid to the heart may be detrimental.
 - Pulmonary embolism.
 - Any local condition in which the garments would interfere, including gangrene, recent skin graft, dermatitis or untreated, infected leg wounds.
- If you are unsure of any contraindications refer to the patient's physician before using the device.

Cautions

- Garments should be removed immediately if the patient experiences tingling, numbness or pain.
- When used for DVT prophylaxis, continuous use is recommended and any interruption of therapy for a substantial length of time should be at the discretion of the physician.

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