## ENT-ACC101 Trapeze Mount Assembly

### Instructions for Use

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
#### Design Policy and Copyright

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies unless otherwise stated.

© ArjoHuntleigh 2016

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either wholly or in part without the consent of ArjoHuntleigh.

Model **ENT-ACC101** is a lifting pole intended to assist patients in positioning themselves on the bed. **ENT-ACC101** is supplied with an adjustable chain and handle, and is suitable for use with ArjoHuntleigh Citadel Plus Bariatric Care System.


*Do not exceed the Safe Working Load of 150kg (330 lbs.)*

*Do not use the lifting poles to push or pull the bed.*

*Do not suspend other equipment from the lifting poles.*

*Only use the recommended chain and handle assembly, i.e. part number 20969.*

*The lifting poles must be used as a set for proper use of ENT-ACC101.*


*Refer to the bed instructions for use to confirm compatibility before using this accessory.*

*The lifting poles must be positioned so that the handle assembly hangs over the middle of the bed.*

*Do not allow the chain and handle assembly to hang outside the perimeter of the bed.*

*To reduce the risk of patient entrapment, the head end panel must remain fitted at all times while the lifting poles are installed on the bed.*

*The bed should not be occupied by a patient while the lifting poles are being fitted.*


*Before installing the lifting poles, check that the bed's lifting pole sockets are not damaged in any way.*

*Ensure that the lifting poles and handle assembly are positioned to suit the size and medical needs of the patient.*

*The chain and handle assembly must not be used by persons under twelve years of age.*


*Inspect the chain, handle assembly and trapeze mount assembly regularly. If any sign of wear or damage is found, remove it from use immediately and replace the complete unit.*

## INSTALLATION & USE


Insert the lifting poles into each socket (1) at the head end of the bed, with the horizontal part of the pole above the center of the bed.

Twist the pole so that the notch (2) locates over the pin at the bottom of the socket.

When the lifting pole is correctly installed, the indicator line (3) will be level with the top of the socket.

Remove the threaded knobs (4) from the trapeze mount assembly (5). From the head end of the bed, slide the trapeze mount assembly (5) over the horizontal part of the lifting poles. Ensure that the holes in the top plate are located over the four projecting pegs located on top of the lifting poles. Secure threaded knobs (4) to the trapeze mount assembly (5).

To adjust the height of the handle above the bed: press desired chain loop onto the “s” hook provided.

*Do not adjust the chain while a patient is holding the handle.*

*Do not twist the chain through more than 180° as this may cause premature wear.*

After use: remove the trapeze mount assembly, lift the poles out of the sockets and store the accessory away from the bed.

## CLEANING

*Remove the trapeze mount assembly from the lifting poles before cleaning.*

*Only use approved disinfectants and cleaning methods. Other products may cause degradation of the equipment and possibly premature failure.*

*Use all cleaners and disinfectants with care. Obey the manufacturer's instructions for use.*

*Do not disinfect the trapeze mount assembly in a pressure chamber, autoclave or steam cleaner as this may damage the product.*


Lifting pole: refer to the instructions for use supplied with the bed.

### Product weight (including trapeze mount assembly)

ENT-ACC101 = 20.1kg

## SPARE PARTS

Item	Part Number	Description	QTY
1	AH1001111	Assembly, Trapeze Mount (Threaded Knobs Included)	1
2	KKX51060	Threaded Knob	3
3	AH1001055	Assembly, Fixed Position Lifting Pole	2
4	AH1001056	Assembly, Trapeze Handle	1

## EXPLANATION OF SYMBOLS USED

Warning of possible hazard to system, patient or staff

Manufacturer


Important Operational Information

Serial Number


Reference Number

Accessory Weight

Date of Manufacturer

Safe Working Load

Refer to Instruction Manual

Separate electrical and electronic components for recycling in accordance with 2012/19/EC (WEEE).

Conforms with the Medical Device Directive (93/42/EEC) and has been subject to the conformity procedures laid down in the council directive

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## NORGE

ArjoHuntleigh Norway AS  
Olaf Helsets vei 5  
N-0694 OSLO  
Tel: +47 22 08 00 50  
Faks: +47 22 08 00 51  
E-mail: no.kundeservice@arjohuntleigh.com

## ÖSTERREICH

ArjoHuntleigh GmbH  
Dörrstrasse 85  
AT-6020 INNSBRUCK  
Tel: +43 (0) 512 204 160 0  
Fax: +43 (0) 512 204 160 75

## POLSKA

ArjoHuntleigh Polska Sp. z o.o.  
ul. Ks Piotra Wawrzyniaka 2  
PL-62-052 KOMORNIKI (Poznan)  
Tel: +48 61 662 15 50  
Fax: +48 61 662 15 90  
E-mail: arjo@arjohuntleigh.com

## PORTUGAL

ArjoHuntleigh em Portugal  
MAQUET Portugal, Lda. (Distribuidor Exclusivo)  
Rua Poeta Bocage n.º 2 - 2G  
PT-1600-233 Lisboa  
Tel: +351 214 189 815  
Fax: +351 214 177 413  
E-mail: Portugal@arjohuntleigh.com

## SUISSE / SCHWEIZ

ArjoHuntleigh AG  
Fabrikstrasse 8  
Postfach  
CH-4614 HÄGENDORF  
Tél/Tel: +41 (0) 61 337 97 77  
Fax: +41 (0) 61 311 97 42

## SUOMI

Oy Vestek AB  
Martinkuja 4  
FI-02270 ESPOO  
Puh: +358 9 8870 120  
E-mail: info@vestek.fi

## SVERIGE

ARJO Scandinavia AB  
Hans Michelsensgatan 10  
SE-211 20 MALMÖ  
Tel: +46 (0) 10 494 7760  
Fax: +46 (0) 10 494 7761  
E-mail: kundservice@arjohuntleigh.com

## UNITED KINGDOM

ArjoHuntleigh UK  
ArjoHuntleigh House  
Houghton Hall Park  
Houghton Regis  
UK-DUNSTABLE LU5 5XF  
Tel: +44 (0) 1582 745 700  
Fax: +44 (0) 1582 745 745  
E-mail: sales.admin@arjohuntleigh.com

## USA

ArjoHuntleigh Inc.  
2349 W Lake Street Suite 250  
US-Addison, IL 60101  
Tel: +1 630 307 2756  
Free: +1 800 323 1245 Institutional  
Free: +1 800 868 0441 Home Care  
Fax: +1 630 307 6195  
E-mail: us.info@arjohuntleigh.com

REV 17: 03/2016

[www.arjohuntleigh.com](http://www.arjohuntleigh.com)

---

## GETINGE GROUP

**Getinge Group** is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **Getinge** and **Maquet**. **Getinge** provides solutions for infection control within healthcare and contamination prevention within life sciences. **Maquet** specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

---

**ARJOHUNTLEIGH**  
GETINGE GROUP

ArjoHuntleigh focuses on patient handling and hygiene, disinfection, DVT prevention, medical beds, therapeutic surfaces and diagnostics.

ArjoHuntleigh AB  
Hans Michelsensgatan 10  
211 20 Malmö, Sweden  
[www.arjohuntleigh.com](http://www.arjohuntleigh.com)