# Lifeguard LG20 & LG50

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Instructions for Use

**ARJOHUNTLEIGH**

GETINGE GROUP

EN

ZH

# Lifeguard LG20 & LG50

## Instructions for Use

**ARJOHUNTLEIGH**  
GETINGE GROUP

## **WARNING**

To avoid injury, always read this Instructions for Use and accompanied documents before using the product.

Mandatory to read the Instructions for Use

Design Policy and Copyright

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies.

© ArjoHuntleigh 2000.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of ArjoHuntleigh.

## **Contents**

<b>1.</b>	<b>Introduction.....</b>	<b>3</b>
<b>2.</b>	<b>Clinical Applications .....</b>	<b>4</b>
<b>3.</b>	<b>Operation.....</b>	<b>5</b>
<b>4.</b>	<b>Accessories.....</b>	<b>19</b>
<b>5.</b>	<b>Decontamination.....</b>	<b>21</b>
<b>6.</b>	<b>Care and Preventive Maintenance .....</b>	<b>23</b>
<b>7.</b>	<b>Warranty and Service .....</b>	<b>25</b>
<b>8.</b>	<b>Dimensions and Data .....</b>	<b>26</b>

## Warnings and Cautions

**WARNINGS** given in this manual identify possible hazards in procedures or conditions which, if not correctly followed, could result in death, injury or other serious adverse reactions.

**Cautions** given in this manual identify conditions or procedures which, if not correctly followed, could result in equipment failure or damage.

*Notes* given in this manual (indicated by ) are used to explain or amplify a procedure or condition.

Optional items not fitted to all trolleys are indicated by ★

### General Warnings

#### **WARNING**

**The age, size and condition of the patient should be assessed by a competent clinically qualified person before using the trolley, giving particular consideration to appropriate use of the safety sides.**

**Do not use accessories that are not designed or approved for use with the trolley.**

**Always apply the brakes when the trolley is stationary.**

**Raise the safety sides before moving the trolley with a patient.**

**Children and restless or disorientated patients should not be left unattended at any time.**

**Before turning (rolling) a patient on the trolley or performing a lateral transfer, adjust the mattress platform to a comfortable working height; be aware that the trolley may overbalance if the patient's weight is not evenly distributed across the mattress platform.**

## 1. Introduction

---

This manual gives instructions for the use and maintenance of the Lifeguard® patient trolleys, model numbers LG20 and LG50.

All models have a radio translucent mattress platform. A full-length X-ray cassette tray can be supplied as an option.

**Model LG20**



**Model LG50**

---

## 2. Clinical Applications

---

<b>Area of Use</b>	A & E departments, admissions, day surgery and general hospital areas.	
<b>Indications</b>	<p>The trolleys may be used for general and emergency admissions, examination, diagnostic X-ray and initial treatment.</p> <p>They are also appropriate for patient transfer between departments, wards and theatre.</p>	
<b>Contra-indications</b>	Maximum safe working load for Model LG50 is 250 kg.	
	Maximum safe working load for Model LG20 is 180 kg.	

### 3. Operation

---


#### **Safety**

For the safety of operator, patient and equipment:

- Always use a mattress of the correct size and type for the trolley.
- Do not expose the equipment, especially the mattress, to naked flames, cigarettes etc.
- Before operating the trolley, make sure that the patient is positioned correctly to avoid entrapment and imbalance.
- When the trolley is operated, check that no obstacles are in its path.
- Keep these instructions handy for reference.

#### **Push Bar Handles**

To adjust the height of the handles, press the spring-loaded release button (1) and raise or lower the handle as required (2). Let go of the button to lock the handle in the high or low position.

**Push Bar Handles**

## Brakes and Steering

Linked brake pedals are situated at the head end of the trolley; model LG50 has additional pedals at the foot end.


When the pedals are horizontal, the castors are free to rotate and swivel.

When the pedals are in the BRAKE position, the brakes are applied on all four castors.


When the pedals are in the STEER position, one foot end castor is prevented from swivelling and the trolley can be pushed from the head end.

### **Easitrack® Fifth Wheel Steering System ★**

When the pedals are in the STEER position, the central fifth wheel is prevented from swivelling; the trolley will keep on a straight course in any direction and can be steered by pivoting about the fifth wheel.


**LG50 Brake and Steering Control**

**LG20 Brake and Steering Control**

## Height and Tilt Adjustment Model LG50

The height and the tilt angle of the mattress platform are adjusted using a group of five foot pedals on either side of the trolley. The three pedals with a black ▲ in the bottom row are used to raise the platform. The two pedals with a red ▼ in the top row are used to lower the platform.

**LG50 Height and Tilt Control Pedals**

### ***Tilt Adjustment***

**To raise one end of the platform:** Operate the ‘raise’ pedal (black ▲ ) nearest that end of the trolley with a pumping action.

**To lower one end of the platform:** Press down on the ‘lower’ pedal (red ▼ ) nearest that end of the trolley and keep it pressed.

*There is an additional ‘lower’ pedal (red ▼ ) that enables head down tilt to be applied from the head end of the trolley.*

### ***Height Variation***


**To raise the platform:** Operate the central ‘raise’ pedal (black ▲ ) with a pumping action.

**To lower the platform:** Press down both ‘lower’ pedals (red ▼ ) and keep them pressed.

If the platform is raised or lowered when tilted, the tilt angle will be maintained until the leading end reaches the limit of travel, when the platform will return to the horizontal condition.

## Height and Tilt Adjustment Model LG20


***Tilt Adjustment*** The angle of head down tilt is adjustable by raising the lever located just underneath the head end of the mattress platform.

### **LG20 Tilt Adjustment**

Take the weight of the mattress platform and raise the lever, then lift or lower the mattress platform to adjust the amount of tilt. Release the lever when the required angle has been reached.

***Height Variation*** The height of the mattress platform can be adjusted from either side of the trolley by using the two foot pedals A and B.

### **LG20 Height Adjustment**

**To raise the platform:** Operate the foot pedal A with a pumping action.


**To lower the platform:** Raise pedal A or press pedal B down and keep it pressed. Release the pedal when the platform is at the desired height or has reached the lowest position.

## Backrest

*The amount of effort needed to raise or lower the backrest will vary depending on patient weight and mobility. With a very heavy patient, the backrest may fall quickly if it is not supported when released.*

**To adjust the backrest:** Grasp the handle firmly, press the release knob and lift the backrest assisted by the gas spring or push it down against the action of the gas spring. Let go of the release knob to lock the backrest at any angle.


**Backrest Adjustment**

## **Safety Sides**

**To raise the safety side:** Hold the top rail and lift until it locks into position at maximum height.

**To lower the safety side:** Hold the top rail with one hand, pull the release button at the foot end and lower the rail toward the head.


**Safety Side Release Button**

## **X-Ray Photography ★**

The mattress and platform are radio translucent. A full-length cassette tray can be included as a factory-fitted option, which can be accessed from either side of the trolley.

**To load or unload cassettes:** Lower the tray by turning the handles at both ends of the trolley to the vertical position. Load or unload the cassette and raise the tray again by turning the handles to the horizontal position.

**Raising and Lowering the X-Ray Cassette Tray**

## **Mattress**

The mattress is secured to the platform by means of hook and loop tapes on the mattress underside and the laminate top.

## Oxygen Therapy Equipment

### WARNING


Ensure that the oxygen line does not become entangled or trapped which could cause damage to the line or loss of supply.

Ensure that when the trolley is operated it does not collide with the installed equipment, especially at low height.

#### **Model LG50**


A recess in the base of the trolley will accommodate a size D, E or F oxygen cylinder. The stainless steel liner can be lifted out for cleaning.

In addition, facilities are provided for the installation and use of oxygen therapy systems. A short rail (Medirail) is provided at the head end of the trolley to accommodate the flow meter and injection suction unit.

**Mounting Rail (Medirail) at Head End**


### **Mounting Bracket in Base of Trolley**

A universal mounting bracket is installed in the base of the trolley to accept the suction receiver jar.

If the bracket is not already fitted, remove the plastic cover from the hole in the base cover moulding and attach the bracket with two 6mm Button Head screws and nuts using the pre-drilled holes in the base frame.

**Model LG20**


The oxygen cylinder holder is mounted on the base of the trolley. It will accommodate a size D, E or F oxygen cylinder.

**Oxygen Cylinder Holder Shown on LG20**

*There is an optional Medirail that can be attached to the head end of the mattress platform on model LG20.*

**LG20 Medirail**


## IV Drip Rod MB-ACC32

Install the IV rod in one of the four tubular sockets at the corners of the platform (1).


### WARNING

**The maximum total load on the hooks must not exceed 5kg**

**The maximum load attached to the lower tube must not exceed 10kg**

**To adjust the height:** hold the telescoping rod firmly, loosen the locking collar (2), lift or lower the inner rod as required (3) and re-tighten the locking collar firmly (4). Fold out the hooks (5).


**IV Drip Rod Location and Operation**


*Be careful not to lift the IV rod out of the socket when raising the hook.*

**Storage:** Adjust the IV rod to minimum length and tighten the locking collar.

The IV rod is stored in a recess on the base of the trolley (model LG50) or in a tube underneath the mattress platform (model LG20).

**LG50 IV Rod Storage**


**LG20 IV Rod Storage**

**Integral IV Drip  
Rod  
799.501**

The integral IV drip rod is located at the head end of the trolley.


**WARNING**

**The maximum total load on the hooks must not exceed 5kg**

Lift the IV rod to the vertical position (1) and allow it to drop into the tubular socket on the platform (2).

**To adjust the height:** hold the telescoping rod firmly, loosen the locking collar (3), lift or lower the hook as required (4) and re-tighten the locking collar firmly (5).

**Integral IV Drip Rod Operation**

*Be careful not to lift the IV rod out of the socket when raising the hook.*

**Storage:** Adjust the IV rod to minimum height and tighten the locking collar, then lift the rod out of the socket and fold it down.

---

## 4. Accessories

---

A range of accessories, including mattresses, is available from ArjoHuntleigh. Details are available on request.

Accessories currently available for *Lifeguard* LG20 and LG50 are shown in the table below.

ACCESSORY NAME	PRODUCT CODE
Loose drip rod	157003
Foot end panel & bow	157009
Foot panel with monitor table	157018
Clip-on table	157021
Arm rest	157006PE02
Safety side pads (grey)	99901PE16
DIN rail kit (head end)	LG/ACC03
Oxygen bottle holder (PD size)	LG/ACC04
DIN rail kit (foot end)	LG/ACC05
Mattress with two-way stretch cover	LGMAT/10
Two part transfer mattress	LGMAT/13
Two part transfer mattress "Pennyla"	LGMAT/13PN
Bi-Flex pressure reducing mattress	TRM2
Bi-Flex pressure redistribution mattress	TRM3

Some accessories may not be available in all markets. ArjoHuntleigh reserves the right to alter specifications and availability without notification.

The following section describes some commonly used accessories. Other accessories are supplied with instructions for use.

### **LGMAT/10**

The LGMAT/10 is a one piece mattress with two-way stretch anti-static PU cover, used to provide the patient with adequate comfort during transport and to facilitate cleaning and disinfection.

### **LGMAT/13**


The LGMAT/13 transfer mattress is designed to allow patients to be moved easily and with minimal discomfort from a *Lifeguard* trolley to an adjacent surface.

**Foot Panel with  
Monitor Table  
157018**

Install the foot panel in the square sockets at the end of the mattress platform, with the clipboard on the outside.


**WARNING**

**Maximum safe working load 12kg evenly distributed.**

To use the monitor table, lift the clipboard to the horizontal position. Use the straps to secure the monitor to the table.

To remove the foot panel, lift it out of the sockets.

**Monitor Table Installed**

---

## 5. Decontamination

---

The following processes are recommended, but should be adapted to comply with the local or national guidelines (Decontamination of Medical Devices) which may apply within the Healthcare Facility or the country of use. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The *Lifeguard* system should be routinely decontaminated between patients and at regular intervals while in use; as is good practice for all reusable medical devices.

### **WARNING**

**Protective clothing should always be worn when carrying out decontamination procedures.**

### **Caution**

**Do not use Phenol-based solutions or abrasive compounds or pads during the decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.**

#### ***To Clean***

Clean all exposed surfaces and remove any organic debris by wiping with a cloth moistened with a simple (neutral) detergent and water. Dry thoroughly.

#### ***Chemical Disinfection***

To protect the integrity of the cover we recommend a chlorine-releasing agent, such as sodium hypochlorite, at a strength of 1,000ppm available chlorine (this may vary from 250ppm to 10,000ppm depending on local policy and contamination status).

Wipe all cleaned surfaces with the solution, rinse and dry thoroughly.

Alcohol based disinfectants (maximum strength 70%) may be used as an alternative.

Ensure the product is dry before storage.

If an alternative disinfectant is selected from the wide variety

available we recommend that suitability for use is confirmed with the chemical supplier prior to use.

**To clean the mattress:** wipe all surfaces with a disposable cloth soaked in a neutral detergent and hand hot water. Rinse with clean water and dry with disposable paper towels.

**To disinfect the mattress:** wipe all surfaces with a solution of sodium dichloroisocyanurate (NaDCC) at a concentration of 1,000ppm (0.1%) available chlorine.

## 6. Care and Preventive Maintenance

This product is subject to wear and tear during use. To ensure that it continues to perform within its original specification, preventive maintenance procedures should be carried out at the intervals shown.

### WARNING

**This list indicates the minimum level of preventive maintenance recommended by the manufacturer. More frequent inspections should be carried out when the product is subjected to heavy use or aggressive environments, or where required by local regulations.**

**Failure to carry out these checks, or continuing to use the product if a fault is found, may compromise the safety of both the patient and carer/user. Preventive maintenance can help to prevent accidents.**

<i>Actions to be done by carer/user</i>	<i>Daily</i>	<i>Weekly</i>	<i>Yearly</i>
Clean and disinfect as shown in Chapter 5	✓		
Check operation of safety sides		✓	
Visually check castors		✓	
Inspect the mattress as described on the next page	✓		
Check safety side pads (if fitted)		✓	

If the result of any of these tests is unsatisfactory, contact ArjoHuntleigh or an approved service agent.

### WARNING

**The procedures below must be carried out by suitably trained and qualified personnel. Failure to do so may result in injury or an unsafe product.**

<b><i>Actions to be done by qualified personnel</i></b>	<b><i>Daily</i></b>	<b><i>Weekly</i></b>	<b><i>Yearly</i></b>
<b>Check the operation of the lifting mechanism as described below</b>			✓
<b>Check the operation of the backrest raise and lower</b>			✓
<b>Check operation of castors, paying special attention to braking and steering functions</b>			✓
<b>Check the hydraulic pump for leaks.</b>			✓
<b>Check all nuts, bolts and other fasteners are present and correctly tightened</b>			✓
<b>Check all accessories, paying particular attention to fasteners and moving parts</b>			✓

### ***Mattress***

Check the cover, especially around the seams, for cuts or splits that may allow fluid ingress. Replace the cover if necessary.

Unzip the cover and inspect the mattress core for signs of fluid ingress, e.g. staining or dampness. Replace the complete mattress if necessary.

### ***Lifting Mechanism***

Check the operation. A steady firm resistance should be felt over the full stroke of the 'raise' pedals. If the pedal pressure is variable or yielding, this may indicate the presence of air in the system and can be remedied as follows:

1. Raise the platform to full height. Continue pumping for four or five pedal strokes then lower the platform fully.
2. Repeat this action two or three times.
3. If firm and steady pedal operation is not restored, contact ArjoHuntleigh or an approved service agent.

---

## **7. Warranty and Service**

---



ArjoHuntleigh standard terms and conditions apply to all sales. A copy is available on request. These contain full details of warranty terms and do not limit the statutory rights of the consumer.

For service, maintenance and any questions regarding this product, please contact your local ArjoHuntleigh office or authorised distributor. A list of ArjoHuntleigh offices can be found at the back of this manual.

### **Product Lifetime**

The lifetime of this product is typically ten (10) years. “Lifetime” is defined as the period during which the product will maintain the specified performance and safety, provided it has been maintained and operated in conditions of normal use in accordance with the requirements in these instructions.

## 8. Dimensions and Data

	Model LG50	Model LG20
Overall length	211 cm	211 cm
Overall width	77 cm	77 cm
Height of mattress platform (max.)	93 cm	92 cm
(min.)	56 cm	60 cm
Tilt (max.) head down	12°	12°
Tilt (max.) foot down	12°	N/A
Backrest adjustment	0 to 87°	0 to 87°
Castors (diameter)	20 cm	20 cm
Product Weight (approximate)	138 kg	120 kg
With X-ray cassette tray	145 kg	127 kg
Safe working load	250 kg	180 kg
		
Mattress size	195 x 67 x 6.5 cm 195 x 67 x 7.5 cm 195 x 67 x 13 cm	195 x 67 x 6.5 cm 195 x 67 x 7.5 cm 195 x 67 x 13 cm

### STORAGE

Handle with care. Do not drop. Avoid shock or violent impact.	
This equipment should be stored in a clean, dry and well-ventilated area which meets the following conditions:	
Temperature	-10°C to +50°C
Relative Humidity	20% to 90% at 30°C, non-condensing
Atmospheric Pressure	700hPa to 1060hPa

Lectrolite® is a registered trademark of Herculite Products Inc.

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**HONG KONG**

ArjoHuntleigh (Hong Kong) Limited  
Room 411-414, 4/F, Manhattan Centre,  
8 Kwai Cheong Road, Kwai Chung, N.T.,  
HONG KONG  
Tel: +852 2960 7668  
Fax: +852 2960 1711

**ITALIA**

Arjo Italia S.p.A.  
Via Giacomo Peroni 400-402  
IT-00131 ROMA  
Tel: +39 (0) 6 87426211  
Fax: +39 (0) 6 87426222  
E-mail: Italy.promo@arjo.com

**MIDDLE EAST**

Huntleigh Middle East  
Office G005 - Nucleotide Complex,  
Biotechnology & Research Park,  
P.O.Box 214742, Dubai,  
United Arab Emirates  
Tel: +971 (0)4 447 0942  
E-mail: Info.ME@arjo.com

**NEDERLAND**

Arjo BV  
Biezenwei 21  
4004 MB TIEL  
Postbus 6116  
4000 HC TIEL  
Tel: +31 (0) 344 64 08 00  
Fax: +31 (0) 344 64 08 85  
E-mail: info.nl@arjo.com

**NEW ZEALAND**

Arjo Ltd  
34 Vestey Drive  
Mount Wellington  
NZ-AUCKLAND 1060  
Tel: +64 (0) 9 573 5344  
Free Call: 0800 000 151  
Fax: +64 (0) 9 573 5384  
E-mail: nz.info@Arjo.com

**NORGE**

Arjo Norway AS  
Olaf Helset's vei 5  
N-0694 OSLO  
Tel: +47 22 08 00 50  
Faks: +47 22 08 00 51  
E-mail: no.kundeservice@arjo.com

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---

ArjoHuntleigh is a world-leading provider of integrated products and solutions that improve the lives of patients and residents with reduced mobility. We help healthcare facilities deliver wellness and effective everyday care, early mobilisation, safe patient handling, venous thromboembolism prevention, pressure injury prevention, hygiene routines, bariatric care and diagnostics. With extensive knowledge and experience, we strive to improve efficiency and ensure a safer and dignified environment for caregivers and their patients.

---

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ArjoHuntleigh AB