INSTRUCTIONS FOR USE

# AtmosAir 9000A

Mattress Replacement System

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Bedienungsanleitung · Gebruiksaanwijzing · Manuel d'utilisation · Istruzioni per l'uso  
Instrucciones de uso · Brugsvejledning · Instruktioner för användning  
Instruções de utilização · Kullanım Talimatları · Οδηγίες χρήσης · Käyttöohjeet  
Bruksanvisning

## **WARNING**

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.

**Mandatory to read the *Instructions for Use*.**

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

The content of this publication may not be copied either whole or in part without the consent of Arjo.

## DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

ARJO HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJO PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY ARJO SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL ARJO BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJO TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties except as set forth in the written limited warranty included in this publication or with this product. Information in this publication may be subject to change at any time. Contact Arjo for updates.

### Important Information For Users

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.
- If applicable, ensure the electrical installation of the room complies with the appropriate national / local electrical wiring standards.

**Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.**

### Notice

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the product information label for specific voltage.

Intentionally left blank

# Table Of Contents

<b>Introduction.....</b>	<b>2</b>
Indications .....	2
Contraindications.....	2
Risks and Precautions.....	2
Safety Information .....	3
Serious Incident.....	4
<b>Preparation for Use .....</b>	<b>5</b>
<b>Mattress Installation.....</b>	<b>5</b>
AtmosAir 9000A MRS.....	5
AtmosAir Pump Installation .....	6
LED Indicators .....	6
Patient Placement and Nursing Care .....	7
Alternating Pressure .....	7
CPR.....	8
Transport Mode .....	8
Skin Care.....	8
Incontinence / Drainage .....	8
General Operation.....	8
<b>Care and Cleaning.....</b>	<b>9</b>
Fire Barrier .....	9
Cover Cleaning .....	9
Cleaning <i>AtmosAir</i> Pump and Tubes.....	10
<b>Preventive Maintenance Schedule.....</b>	<b>11</b>
Daily Cleaning .....	11
Weekly Cleaning.....	11
Inspection / System Checks .....	11
<b>Troubleshooting .....</b>	<b>12</b>
<b>Parts Diagram - <i>AtmosAir</i> with SAT9000A.....</b>	<b>13</b>
<b>Parts Diagram - <i>AtmosAir</i> Pump 9000A .....</b>	<b>14</b>
<b>Replacement Parts .....</b>	<b>15</b>
<b>Specifications .....</b>	<b>16</b>
<b>Electromagnetic Compatibility.....</b>	<b>17</b>
<b>Symbols Used.....</b>	<b>20</b>
<b>Customer Contact Information.....</b>	<b>21</b>

# Introduction

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications, Safety Information and Risks and Precautions** sections prior to placing a patient on any AtmosAir™ with SAT™ Mattress Replacement System (MRS).

Caregivers should review this information with the patient and the patient's family and / or legal guardian. Save this User Guide in an easily accessible location for quick reference.

The *AtmosAir* 9000A MRS comes in a wide variety of lengths and widths to support various customer requirements. Contact an Arjo representative for additional product information.

## Indications

AtmosAir™ 9000A Mattress Replacement System (MRS) is intended to treat and prevent all categories of pressure ulcers.

## Contraindications

- unstable vertebral fracture
- cervical and skeletal traction

## Risks and Precautions

**Transfer** – Standard precautions should be taken during patient transfer.

**Side Rails and Restraints – WARNING:** Use or non-use of restraints, including side rails, can be critical to patient safety. Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. **See related Safety Information.**

- **The AtmosAir 9000A is a mattress replacement system - remove all other mattress components.**
- **To ensure compliance with IEC 60601-2-52:2010, an approved Arjo mattress should be used.**
- **The user is responsible for validating other combinations of bed and mattress to ensure compliance with IEC 60601-2-52:2010.**
- **For more information on suitable beds, contact your local Arjo office or approved distributor.**

**Patient Migration** – Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. **Monitor patients frequently to guard against patient entrapment.**

**Oxygen Use – DANGER:** Risk of explosion if the *AtmosAir* Pump is used in the presence of flammable anesthetics. Use of this product's pump in an oxygen-enriched environment may produce potential of fire hazard. This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. Unplug and do not use pump when using oxygen-administering equipment other than the nasal mask or half-bed-length tent type.

**Shock Hazard** – Electrical shock hazard; do not remove pump case covers. Refer to qualified service personnel.

**Transport** – To avoid deflation of the mattress SAT cells, either disconnect all hoses from mattress at the mattress connector or disconnect from the pump and fit the transport cap to the hose.

## Safety Information

**Patient Entrance / Exit** – Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

**Brakes** – Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

**Bed Height** – To minimize risk of falls or injury, the bed should always be in the lowest practical position when the patient is unattended.

**Bed Frame** – Always use a standard healthcare bed frame with this mattress, with any safeguards or protocols that may be appropriate. Bed frame and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. It is recommended that bed and side rails (if used) comply with all applicable regulations and protocols.

**CPR** – Level the bed and disconnect pump (if present). Lower side rails and initiate CPR per facility protocols. Consider use of backboard if indicated. After CPR, remove backboard (if used) reconnect pump, if present, raise siderails and reconfigure bed and accessories as in initial placement.

**Head of Bed Elevation** – Keep head of bed as low as possible to help prevent patient migration.

**Side Rails / Patient Restraints** - Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. Consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patient frequently to guard against patient entrapment.

**CAUTION:** When selecting a standard mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

**Skin Care** – Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

**I.V. and Drainage Tubes** - I.V. and drainage tubes should always have slack for alternating pressure or rotation and other patient movements.

**Fluids** – Avoid spilling fluids on pump controls. If spills do occur, clean fluid from pump wearing rubber gloves or while unit is unplugged to avoid any possibility of shock. Once fluid is removed, check operation of components in area of spill.

**Fluids remaining on controls can cause corrosion, which may cause components to fail or operate erratically, possibly producing potential hazards for patient and staff.**

**Avoid Fire Hazards** – To minimize risk of fire, connect the bed's power cord directly into a wall-mounted outlet. Do not use extension cords or multiple outlet strips. In the US, review and follow FDA's Safety Tips for Preventing Hospital Bed Fires.

**No Smoking in Bed** – Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

**Power Cord** – Position power cord to avoid a tripping hazard and / or damage to the cord. Ensure power cord is kept free from all pinch points and moving parts and is not trapped under casters. Improper handling of the power cord can cause damage to the cord, which could produce risk of fire or electric shock.

**General Protocols** – Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

**Fire Barrier Contact** – Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) be used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. If the fire barrier becomes damaged or soiled, it should be replaced.

#### **End of Life Disposal –**

- Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.
- Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.
- Pump units have electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

## **Serious Incident**

**If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.**

## Preparation for Use

For information concerning the bed frame, refer to the manufacturers' User Guide.

1. Open shipping container(s).

***Do not use sharp instruments to open boxes. Damage to mattress could result.***

2. Remove *AtmosAir* 9000A MRS from plastic protective cover.

**The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to accommodate; see Troubleshooting for more information. Wrinkles will not affect inflation or function, so mattress may be used immediately if needed.**

3. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing mattress onto a new frame or for a new patient, check mattress surface for staining and soiling; clean and / or disinfect as required (see **Care and Cleaning**).
5. Level bed and lock brakes.
6. Remove existing mattress from bed frame.

## Mattress Installation

### AtmosAir 9000A MRS

1. Position mattress on bed frame with logo facing up and product information tags at foot end of bed.
2. Ensure there are no gaps between mattress and bed frame or side rails.

***Always use a standard healthcare bed frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.***

## AtmosAir Pump Installation

### Pump installation

1. Place pump on a solid stationary surface or suspend on end of bed frame with built-in hanger.
2. Attach pump hoses to mattress, as shown with (Figure 1) below:
  - Attach hoses to the mattress connectors on the **patient left side** of the mattress for **alternating pressure**.
  - Ensure pump hoses are attached to pump.
3. Ensure air hoses are not kinked and will not be pinched by any articulated bed mechanisms.
4. Attach hose connector to the pump.
5. Plug pump unit into a properly-grounded wall outlet.
6. Verify power to this outlet is not controlled by a wall switch.
7. Switch pump on.


Figure 1 Mattress Hose Connector

### LED Indicators

1. Green LED is lit when the pressure is at the set level.
2. Yellow LED is lit when the pressure is below the set level. Check connections.


Figure 2 Control Panel

## Patient Placement and Nursing Care

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications**, **Safety Information** and **Risks and Precautions** sections prior to placing a patient on the *AtmosAir* 9000A MRS.

1. Transfer patient following all applicable safety rules and institution protocols.
2. Center patient side-to-side and head-to-foot on *AtmosAir* MRS surface.
3. Ensure all sections of the mattress fully support the patient.

**Mattress handles are for mattress transport only.**

## Alternating Pressure

The *AtmosAir* with SAT 9000A, model provides alternating pressure. The level of alternating pressure can be adjusted for patient comfort. For maximum effect, it is recommended to begin at maximum intensity and adjust as needed.


***Prior to engaging alternating pressure, ensure that bed frame has side rails and that all side rails are fully engaged in full upright and locked position. Observe patient during initial turning cycles to ensure proper patient placement***

***Ensure pump hoses are properly connected to the mattress (see *AtmosAir* Pump Installation).***

- Rotate the control knob on the pump clockwise to increase the intensity of alternating pressure.
- Rotate the control knob on the pump counter-clockwise to reduce the intensity of alternating pressure.


**Figure 3 Control Knob**

**The *AtmosAir* MRS DOES NOT provide a Nurse Hold function.**

## CPR

1. Level bed.
2. Disconnect hoses from pump to level mattress, as shown at right (Figure 4).
3. Lower or remove side rail on caregiver's side if necessary.
4. Begin CPR as per facility protocol. Consider use of backboard if indicated.
5. After CPR is performed:
  - Remove backboard if used.
  - Reconnect pump hoses to mattress.
  - Raise or install side rail as necessary.
  - Reconfigure bed and accessories as in initial placement.

**Figure 4 Disconnect Hoses**

## Transport Mode

When disconnected from the pump, the AtmosAir 9000A MRS will operate as a reactive surface with the full operation of the SAT system. For transport, there are two options:

1. Disconnect the tube assembly from the mattress. Valves will close off as soon as the tube assembly is disconnected.
2. Disconnect the tube assembly from the pump and securely place the transport cap on the tube assembly to maintain a seal.

## Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under patient are not wrinkled.

## Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace bed linens as required (see **Care and Cleaning**).

## General Operation

Avoid contact of sharp instruments with the *AtmosAir* MRS. Punctures, cuts and tears may prevent proper inflation and air pressure maintenance.

**Figure 5 Avoid Sharp Instruments**

# Care and Cleaning

The following processes are recommended, but should be adapted to comply with local institution protocols. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The *AtmosAir* MRS should be routinely decontaminated between patients and at regular intervals while in use.

***If a powered AtmosAir MRS, ensure power cord is unplugged from wall outlet before performing care and cleaning.***

***Do not use Phenol-based solutions or abrasive compounds or pads on the cover during decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.***

***Wipe down the cover immediately after it has been exposed to fluids or spills.***

## Fire Barrier

Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) be used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. Fire barrier should be replaced if soiled or damaged.

## Cover Cleaning

### Sewn AtmosAir MRS (Top cover can not be removed from base)

1. If a powered MRS, unplug pump before cleaning.
2. Remove or push bed linens to center of mattress.

***The sewn AtmosAir MRS cover is cleaned by wipe down methods only. Do not launder as damage to the base may occur.***

3. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
4. Rinse with a wet towel to remove chemicals.
5. After wipe down, dry surface with towel.
6. Ensure bed linens are refitted and not wrinkled under patient.
7. Clean pump and tubes (if needed) by wiping with a damp cloth.

## RF Welded *AtmosAir* MRS (Detachable top cover)

1. If a powered MRS, unplug pump before cleaning.
2. Push bed linens to center of mattress to wipe down. Remove to launder top cover or wipe down as described below.

***The RF Welded *AtmosAir* MRS has a detachable top cover that can be laundered. However the base is cleaned by wipe down methods only. Do not launder base or fire barrier as damage may occur.***

3. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
4. Rinse with a wet towel to remove chemicals.
5. After wipe down, dry surface with towel.
6. Ensure bed linens are refitted and not wrinkled under patient.
7. Clean pump and tubes (if needed) by wiping with a damp cloth.

### Laundering Detachable Top Cover

1. Unzip the top cover from the base for laundering. Do not launder base as damage may occur.
2. Recommended wash temperature for top cover is 60° C (140° F) for 15 minutes.
3. Maximum wash temperature is 95° C (203° F) for 15 minutes.
4. Tumble dry at 60° C (140° F) or air dry.
5. Maximum drying temperature 80° C (176° F).

## Cleaning *AtmosAir* Pump and Tubes

1. Ensure pump is unplugged.
2. Wipe surface of pump and tubes with coarse cloth using an approved disinfectant germicide mixed to manufacturer's instructions. **Do not flood any part of the pump with cleaning solution.**

***Avoid spilling fluids on pump controls. If spills do occur, clean fluid from pump wearing rubber gloves or while unit is unplugged to avoid any possibility of shock. Once fluid is removed, check operation of components in area of spill. Fluids remaining on controls can cause corrosion, which may cause components to fail or operate erratically, possibly producing potential hazards for patient and staff.***

3. Wipe off excess solution.
4. Disinfect pump and tubes with a chlorine solution (mixed according to the instructions in the **Care and Cleaning** section). Using a clean cloth, wring out excess solution until cloth is damp. Wipe pump surface and tubes with damp cloth.
5. Allow to air dry.

# Preventive Maintenance Schedule

Preventive maintenance for the *AtmosAir* MRS consists of regular cleaning (see **Care and Cleaning**) and an overall system check-out to be performed at the intervals described below.

**All components must be cleaned, disinfected and inspected after each patient's use and before use by a new patient.** Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

## Daily Cleaning

The cover should be wiped daily with a mild soap and water solution. After wipe down, dry surface with a towel.

## Weekly Cleaning

The pump and hoses should be cleaned weekly.

## Inspection / System Checks

Check each of the following before placing the *AtmosAir* MRS with a new patient:

1. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
2. Ensure mattress is free of stains and is not overly faded.
3. Ensure air inlet hoses and connectors on mattress and pump are clean and undamaged.
4. Ensure pump and power cord are clean and undamaged.
5. Ensure pump hanger brackets are secure and operate correctly.
6. Ensure power switch and comfort control knob operate correctly.
7. Attach pump to the hoses and power on to ensure there are no air leaks.

# Troubleshooting

Do not attempt troubleshooting outside this guide or where the solution recommends to contact Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

SYMPTOM	POSSIBLE CAUSE	SOLUTION
Green lamp on pump does not light when power switch is set to the ON position.	Power cord may be unplugged.	Plug in power cord.
	Wall outlet may be controlled by a wall switch.	Plug power cord into different wall outlet.
	Circuit breaker for wall switch may be tripped or blown.	Reset circuit breaker or replace fuse.
	Power cord or power switch may be damaged.	Contact Arjo for assistance.
Mattress too firm upon arrival.	Difference in altitude not sufficient to open valves.	Apply weight to mattress to open valves.
Mattress cover too wrinkled upon removal from shipping container.	Internal components have not accommodated to environment. This does not affect inflation or function.	Let mattress accommodate for 24 hours. If problem continues, contact Arjo for assistance.
Mattress is not firm. In use.	Tubing not connected properly.	Check tubing inside mattress for loose connectors.
	Tubing kinked.	Check tubing inside mattress for possible kinks.
	Tubing disconnected.	Check tubing inside mattress for possible disconnect.
	Holes in or damaged to SAT system.	Check SAT system for holes or damage, or contact Arjo for assistance.
	Pump switched off.	Switch pump on.
	Pump leak low pressure indicator active.	Remove and replace pump.
Mattress deflates during transport	Hoses not disconnected at mattress.	Disconnect hoses at mattress or reconnect pump.
	Hose cover not fitted.	Connect hose cap or reconnect pump.


# Parts Diagram - AtmosAir with SAT 9000A

All hoses were removed to improve diagram readability.

The fire barrier (standard on all mattresses) is a sleeve that fits over the foam and SAT assemblies.


The SAT Cover, SAT Cells, SAT Interliner and rotation bladders are all part of the integrated SAT System and cannot be ordered separately. See Replacement Parts for a complete list of SAT Systems.

# Parts Diagram - *AtmosAir* Pump 9000A

See Replacement Parts for a complete list of parts that may be ordered for the *AtmosAir* Pump, as well as entire pump systems.

# Replacement Parts

Do not attempt troubleshooting, maintenance or parts replacement outside this manual or where the solution recommends contacting Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Replaceable *AtmosAir* components are listed below. For more information such as pricing, additional options, or spare parts that are not on this list, please contact your local Arjo representative. Contact Arjo for additional assistance.

## ATMOSAIR MATTRESS COMPONENTS

### *AtmosAir* Foam Shells

PART TYPE.....	CATALOG NUMBER
9000 Foam Shell .....	KA4/9FS(SIZE)

### *AtmosAir* with SAT System

PART TYPE.....	CATALOG NUMBER
9000A SAT System.....	KA9ASAT

### *AtmosAir* Fire Barrier Lining

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> Fire Barrier Lining (standard on all sewn models) .....	KA4/9-FBL

Welded *AtmosAir* Fire Barrier Lining (Standard on all *AtmosAir*

Welded mattresses (all models).....*AtmosAir*-FBL

## ATMOSAIR MATTRESS COVERS

### *AtmosAir* Model 9000A Mattress with Alternating Pressure

PART TYPE.....	CATALOG NUMBER
Complete Cover, <i>Reliant IS</i> <sup>2</sup> , Welded Top, Non-Skid PU Bottom, No Handles.....	AXCRFPN(WWWLL)TXS
Complete Cover, <i>Reliant IS</i> <sup>2</sup> , Welded Top, Non-Skid PU Bottom, With Handles .....	AXCRFPU(WWWLL)TXS
Top Cover, <i>Reliant IS</i> <sup>2</sup> , Welded .....	AXCRFXX(WWWLL)TXS
Base Cover, PU Non-skid Bottom, No Handles.....	AXCXXPN(WWWLL)XXS
Base Cover, PU Non-skid Bottom, With Handles .....	AXCXXPU(WWWLL)XXS
Complete Cover, <i>Reliant IS</i> <sup>2</sup> Top, Non-Skid Bottom .....	AXCRLNB(WWWLL)TXS
Complete Cover, <i>Reliant IS</i> <sup>2</sup> Top, Vinyl Gray Bottom .....	AXCRLVG(WWWLL)TXS
Complete Cover, <i>Soflux</i> Top, Non-Skid Bottom, No Handles ....	AXCSONH(WWWLL)NXS
Complete Cover, <i>Soflux</i> Top, Non-Skid Bottom, With Handles ..	AXCSONB(WWWLL)NXS
Complete Cover, <i>Reliant IS</i> <sup>2</sup> , RF Welded Top, Non-Skid PU Bottom, No Handles.....	AXCRFPN(WWWLL)TXS
Complete Cover, <i>Reliant IS</i> <sup>2</sup> , RF Welded Top, Non-Skid PU Bottom, With Handles .....	AXCRFPU(WWWLL)TXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Bottom, No Handles.....	AXCPRPN(WWWLL)NXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Bottom, with Handles .....	AXCPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS</i> <sup>2</sup> , RF Welded .....	AXCRFXX(WWWLL)TXS
Top Cover, Premium RF Welded .....	AXCPRXX(WWWLL)NXS
Base Cover, PU Non-Skid Bottom, No Handles .....	AXCXXPN(WWWLL)XXS
Base Cover, PU Non-Skid Bottom, With Handles .....	AXCXXPU(WWWLL)XXS

## ATMOSAIR PUMPS

### AtmosAir Standard Pumps

PART TYPE.....	CATALOG NUMBER
----------------	----------------

AtmosAir 9000A Pump only

Europe .....	KAPUMP-EU
--------------	-----------

UK .....	KA9PUMP-UK
----------	------------

AtmosAir Harness / Tubing / Connectors (no pump) .....	KA9PUMPASMBLY
--	---------------

AtmosAir 9000A Pump Filter .....	KA9EZPMPFILTER
----------------------------------	----------------

AtmosAir 9000A Pump Filter Cap .....	KA9EZPMPFILCAP
--------------------------------------	----------------

## ATMOSAIR LITERATURE

PART TYPE.....	CATALOG NUMBER
----------------	----------------

AtmosAir 9000A MRS User Guide .....	416412
-------------------------------------	--------

## Specifications

Specifications subject to change without notice.

### Mattresses

Maximum Recommended Patient Weight* .....	500 lb (227 kg)
---	-----------------

Weight (based on 35 in x 80 in x 7 in MRS, other sizes will vary)

AtmosAir with SAT 9000A MRS .....	40 lb (18 kg)
-----------------------------------	---------------

### Dimensions (varies depending on model)\*\*

Minimum Length .....	75 in (190.5 cm)
----------------------	------------------

Maximum Length .....	84 in (213 cm)
----------------------	----------------

Minimum Width .....	32 in (81.5 cm)
---------------------	-----------------

Maximum Width .....	42 in (106.7 cm)
---------------------	------------------

Height .....	7 in (17.5 cm)
--------------	----------------

### Pump 230V

Weight .....	5 lb (2.5 kg)
--------------	---------------

Dimensions .....	10.25 in (26 cm) L x 4.5 in (11.5 cm) W x 3.75 in (9.5 cm) H
------------------	--

### Electrical:

Volts .....	230 VAC
-------------	---------

Frequency .....	50 Hz
-----------------	-------

Amps .....	1 Amps
------------	--------

Power Cord .....	12 ft
------------------	-------

Output .....	6 litres per minute
--------------	---------------------

Cycle Time .....	12 minutes
------------------	------------

\* Patient weight capacity may vary depending on frame use. Please consult frame manufacturer.

\*\* Additional sizes are available. Contact Arjo for more information. AtmosAir 9000 A Series up to 42 in (106.7 cm) with *Reliant IS<sup>2</sup>* cover and AtmosAir 9000A with Premium cover.

# Electromagnetic Compatibility

Although this equipment conforms with the intent of the directive 89/336/EEC in relation to electromagnetic compatibility (EMC), all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Radios, cell phones and similar devices may affect this equipment and should be kept at least 10.8 ft (3.3m) away from the equipment.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the following tables.

The following tables document compliance levels and guidance from the IEC 60601-1-2; 2007 Standard, for the electromagnetic environment in which the AtmosAir MRS should be used in a clinical environment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The AtmosAir MRS is intended for use in the electromagnetic environment specified below. The customer or user of the AtmosAir MRS should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	The AtmosAir MRS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Yes	

## Guidance and manufacturer's declaration - electromagnetic immunity

The AtmosAir MRS is intended for use in the electromagnetic environment specified below. The customer or user of the AtmosAir MRS should assure that it is used in such an environment.

Immunity Test		IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2		±6 kV Contact ±8 kV Air	±6 kV Contact ±8 kV Air	In accordance with IEC 60601-1-2:2007, floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst  IEC 61000-4-4	230VAC Pump	±2 kV	±2 kV	N/A
Surge IEC 61000-4-5		1 kV line(s) to Line(s) 2 kV line(s) to Earth	1 kV line(s) to Line(s) 2 kV line(s) to Earth	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11		5% half cycle 40% for 5 cycle 70% for 25 cycle  5% for 5 seconds	5% half cycle 40% 5 cycles 70% 25 cycles  5% for 5 seconds	N/A
Power frequency (50Hz/60Hz)magnetic field IEC 61000-4-8		3 A/M	3 A/M	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U<sub>i</sub> is the a.c. mains voltage prior to application of the test level.

## Recommended separation distances between portable and mobile RF communications equipment and the AtmosAir MRS

The AtmosAir MRS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AtmosAir MRS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AtmosAir MRS as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter  W	Separation distance according to frequency of transmitter  meters		
	150 kHz to 80 MHz  not applicable	80 MHz to 800 MHz  $d = [ \frac{P}{f} ] \sqrt{P}$	800 MHz to 2.5 GHz  $d = [ \frac{P}{f} ] \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.37	0.74
1	N/A	1.2	2.3
10	N/A	3.7	7.4
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separate distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1, At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects and people.

### Guidance and manufacturer's declaration - electromagnetic immunity

The AtmosAir MRS is intended for use in an electromagnetic environment specified below. The customer or user of the AtmosAir MRS should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3Vrms 150k - 80 MHz	3 Vrms 150k - 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the AtmosAir MRS, including cables, than the recommended separation distance calculated from the equation application to the frequency of the transmitter.  Recommended Separation Distance Battery Operated Device 3.5 $d = \left[ \frac{E_1}{7} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{E_1}{7} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 Vrms 80 MHz - 2.5 GHz	3 Vrms 80 MHz - 2.5 GHz	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (see note a) should be less than the compliance level in each frequency range (see note b). Interference may occur in the vicinity of equipment marked with the following symbol:  


NOTE 1, At 80 MHz and 800 MHz, the higher frequency range applies.


NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AtmosAir MRS is used exceeds the applicable RF compliance level above, the AtmosAir MRS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the AtmosAir MRS.

b) Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

# Symbols Used

	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745
	Important Operational Information
	Foot End
	Warning of possible hazard to system, patient or staff
	Consult User Guide
	Power
	Manufacturer
	Safe Working Load
	Wash at 60°C, 95°C Max for 15 minutes
	Do Not Iron

<b>IPX0</b>	Protection from ingress of fluids
	Type BF Applied Part
	Possible electrical shock hazard
	Attention - See User Guide
	On
	Off
	CE marking indicating conformity with European Community harmonised legislation. Figures indicate Notified Body supervision.
	Wipe down Only
	Tumble Dry at 60°C, 80°C Max
	Do not use Phenol-based cleaning solutions

Use Solution diluted to 1000ppm of available chlorine

Double Insulated

Serial Number

Do Not flip mattress

This product or its parts are designated for separate collection at an appropriate collection point. At the end of useful life, dispose of all waste according to local requirements, or contact your local Arjo representative for advice.

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

The operator must read this document (this instruction for use) before use. Note: the symbol is blue on the product label.

Model number

Do Not rotate mattress

## Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative, or visit [www.Arjo.com](http://www.Arjo.com)

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Intentionally left blank

AUSTRALIA  
Arjo Australia Pty Ltd  
78, Forsyth Street  
O'Connor  
AU-6163 Western Australia  
Tel: +61 89337 4111  
Free: +1 800 072 040  
Fax: +61 89337 9077

BELGIQUE / BELGIË  
Arjo NV/SA  
Evenbroekveld 16  
BE-9420 ERPE-MERE  
Tél/Tel: +32 (0) 53 60 73 80  
Fax: +32 (0) 53 60 73 81  
E-mail: info.belgium@arjo.be

BRASIL  
Arjo Brasil Equipamentos Médicos Ltda  
Rua Marina Ciufulli Zanfelice, 329 PB02  
Galpão - Lapa  
São Paulo - SP - Brasil  
CEP: 05040-000  
Phone: 55-11-3588-5088  
E-mail: vendas.latam@arjo.com  
E-mail: servicios.latam@arjo.com

CANADA  
Arjo Canada Inc.  
90 Matheson Boulevard West  
Suite 300  
CA-MISSISSAUGA, ON, L5R 3R3  
Tel/Tél: +1 905 238 7880  
Free: +1 800 665 4831 Institutional  
Free: +1 800 868 0441 Home Care  
Fax: +1 905 238 7881  
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA  
Arjo Czech Republic s.r.o.  
Na Strži 1702/65  
140 00 Praha  
Czech Republic  
Phone No: +420225092307  
e-mail: info.cz@arjo.com

DANMARK  
Arjo A/S  
Vassingerødvej 52  
DK-3540 LYNGBE  
Tel: +45 49 13 84 86  
Fax: +45 49 13 84 87  
E-mail:  
dk\_kundeservice@arjo.com

DEUTSCHLAND  
Arjo GmbH  
Peter-Sander-Strasse 10  
DE-55252 MAINZ-KASTEL  
Tel: +49 (0) 6134 186 0  
Fax: +49 (0) 6134 186 160  
E-mail: info-de@arjo.com

ESPAÑA  
Arjo Ibérica S.L.  
Parque Empresarial Rivas Futura, C/  
Marie Curie 5  
Edificio Alfa Planta 6 oficina 6.1.-62  
ES-28521 Rivas Vacia, MADRID  
Tel: +34 93 583 11 20  
Fax: +34 93 583 11 22  
E-mail: info.es@arjo.com  
FRANCE  
Arjo SAS

2 Avenue Alcide de Gasperi  
CS 70133  
FR-59436 RONCQ CEDEX  
Tél: +33 (0) 3 20 28 13 13  
Fax: +33 (0) 3 20 28 13 14  
E-mail: info.france@arjo.com

HONG KONG  
Arjo Hong Kong Limited  
Room 411-414, 4/F,  
Manhattan Centre,  
8 Kwai Cheong Road,  
Kwai Chung, N.T.,  
HONG KONG  
Tel: +852 2960 7600  
Fax: +852 2960 1711

ITALIA  
Arjo Italia S.p.A.  
Via Giacomo Peroni 400-402  
IT-00131 ROMA  
Tel: +39 (0) 6 87426211  
Fax: +39 (0) 6 87426222  
E-mail: Italy.promo@arjo.com

MIDDLE EAST  
Arjo Middle East FZ-LLC  
Office 908, 9th Floor,  
HQ Building, North Tower,  
Dubai Science Park,  
Dubai Barsha South  
P.O Box 11488, Dubai,  
United Arab Emirates  
Direct +971 487 48053  
Fax +971 487 48072  
Email: Info.ME@arjo.com

NEDERLAND  
Arjo BV  
Blezenwei 21  
4004 MB TIEL  
Postbus 6116  
4000 HC TIEL  
Tel: +31 (0) 344 64 08 00  
Fax: +31 (0) 344 64 08 85  
E-mail: info.nl@arjo.com

NEW ZEALAND  
Arjo Ltd  
34 Vestey Drive  
Mount Wellington  
NZ-AUCKLAND 1060  
Tel: +64 (0) 9 573 5344  
Free Call: 0800 000 151  
Fax: +64 (0) 9 573 5384  
E-mail: nz.info@Arjo.com

NORGE  
Arjo Norway AS  
Olaf Helsetts vei 5  
N-0694 OSLO  
Tel: +47 22 08 00 50  
Faks: +47 22 08 00 51  
E-mail: no.kundeservice@arjo.com

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Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB  
Hans Michelsensgatan 10  
211 20 Malmö, Sweden  
[www.arjo.com](http://www.arjo.com)

**arjo**

**CE**  
2797

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Arjo Inc.