
USER MANUAL

First Step Select

Therapy System

EN · FR · ES

Manuel D'utilisation • Manual Del Usuario

WARNING

To avoid injury, always read this Instructions For Use and accompanied documents before using the product.

Mandatory to read the Instructions For Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

IMPORTANT INFORMATION FOR USERS

In order for the product to perform properly, the manufacturer recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- **WARNING:** Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by the manufacturer. Contact the manufacturer for information regarding maintenance and repair.
- Ensure the electrical installation of the room complies with the appropriate national electrical wiring standards.

Specific indications, contraindications, warnings, precautions and safety information exist for therapeutic surface products. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

NOTICE

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the product information label for specific voltages.

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Table Of Contents

Introduction	8
• Indications	8
• Contraindications.....	8
• Risks and Precautions.....	8
• Safety Information	9
• Maximum Patient Weights.....	10
• Connecting the System to Other Devices	10
Preparation For Use and Patient Placement.....	11
• Unboxing and Inspection.....	11
• Unit Setup and Preparation for Patient Placement	11
Patient Transfer to the First Step Select Therapy System.....	13
• Ambulatory Patient Transfer.....	13
• Non-Ambulatory Patient Transfer.....	13
Patient Transfer from the First Step Select Therapy System	13
• Ambulatory Patient Transfer.....	13
• Non-Ambulatory Patient Transfer.....	13
• Air Pressure Adjustment.....	14
• Pressure Adjust HT/WT Method.....	14
• Pressure Adjust Manual Method	14
Completion of Patient Placement	15
• Home Display	15
InstaFlate Inflation Function (IF).....	15
Seat Deflate.....	16
Alerts	16
Warmer Adjust	16
Lock-Out Procedures	17
• All Button Lock-Out	17
• Pressure Adjust Button Lock-Out	18
Nursing Care	19
Cardio Pulmonary Resuscitation (CPR).....	19
Patient Transfer to the First Step Select Therapy System	19
Patient Transfer While on First Step Select Therapy System	19
Patient Bathing	19
Bedpan Placement.....	20
Bedpan Removal.....	20
Repositioning Patient using Dri-Flo™ Pads	20
Incontinence / Drainage	21
Skin Care.....	21

Care and Cleaning	21
Daily Care and Cleaning of First Step Select Therapy System While in Use.....	21
Weekly Care and Cleaning of First Step Select Therapy System While in Use.....	21
• Disassembly of First Step Select Therapeutic Surface and Cover Sheet for Laundering and Cleaning.....	21
Infection Control Procedures.....	21
• Cleaning (Single Patient Use).....	22
• Infection Control Protocol.....	22
• Laundering Procedure for Cover Sheet.....	22
• Laundering Procedure for Overlays, Mattress and Tote Bag	22
• Cleaning Procedure for Therapy Control Unit and Hoses.....	23
• Cleaning Procedure for Foam Base	23
• Cleaning and Disinfection for V-Series Adaptor	24
• Cleaning Procedure for Air Filter	24
• Institutional Cleaning (DMEs and Extended Care Facilities).....	24
• Assembly of the First Step Select Therapy System after Laundering and Cleaning.....	24
Maintenance	25
Preventative Maintenance.....	25
Unscheduled Maintenance.....	25
Troubleshooting	26
• Troubleshooting Chart for the First Step Select Therapy System.	26
Replacement Parts and Inventory List	28
Product Components.....	28
Product Literature.....	28
Accessories	28
Miscellaneous.....	28
Specifications	30
Dimensions:.....	30
Weight:	30
Electrical:	30
Symbols Used	31
Miscellaneous.....	32
• EMC	32
• Classification	32
• Environmental Storage / Transit Conditions	32
• Environmental Operating Conditions.....	32
• Leakage Testing and Ground Continuity Testing.....	32
• Leakage Testing	32
• Ground Continuity Testing	32
Questions and Information.....	33

Introduction

The First Step® Select Therapy System is a pressure relieving overlay (when used on an existing mattress) or a mattress replacement system (MRS) (when used with the optional foam base) designed to suspend the patient on a low air-loss support surface, providing effective pressure redistribution and for prevention and treatment of pressure injuries (Stages 1-4).

Each therapeutic surface is constructed with a moisture vapor permeable surface which is divided into three sections: HEAD, BODY and LEG. The pressure in each section can be adjusted to meet individual patient needs.

The therapy control unit consists of a single variable speed blower with an air warmer to provide air temperatures that can be adjusted for patient comfort.

The *First Step* Select Therapeutic Surface is a Type BF applied part.

The expected service life for the *First Step* Select Therapy Unit is five years.

Indications

The *First Step* Select Therapy System is indicated for patients who would benefit from a pressure redistribution surface assisting in the prevention and treatment of pressure injuries (Stages 1-4).

Contraindications

Patient conditions for which the application of pressure redistributing therapy on the *First Step* Select Therapeutic Surface are contraindicated include:

- Unstable spinal cord injury.
- Unstable vertebral fracture.
- Cervical traction.

Risks and Precautions

- Transfer - Precaution should be taken during patient transfer, including the locking of caster brakes and deflation of therapeutic surface. Refer to **Patient Transfer** sections in the **Patient Placement** and **Nursing Care** chapters of this manual.
- Side Rails and Restraints - Use or non-use of restraints, including side rails, can be critical to patient safety. **WARNING:** Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. See related Safety Information.
- Patient Migration - **WARNING:** As with all speciality bed products that are designed to reduce shear and pressure on the patient's skin, the risk of gradual movement and / or sinking into hazardous positions of entrapment and / or inadvertent bed exit may be increased.
- Skeletal Traction or Unstable Fracture (if not contraindicated) - With skeletal traction, unstable pelvic fractures or any other unstable fracture (to the extent not contraindicated), maintain physician directed angle of articulation and guard against risks of patient migration or inadvertent deflation of surface.
- Oxygen Use - **Danger:** Risk of explosion if used in the presence of flammable anesthetics. Use of this product's therapy control unit in an oxygen-enriched environment may produce potential of fire hazard. This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Position therapy control unit only at foot end of bed when using oxygen administering equipment.

- Electromagnetic Interference - Although this equipment conforms with the intent of the directive 89 / 336 / EEC in relation to electromagnetic compatibility, all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.
- Shock Hazard - Electrical shock hazard, do not remove cover from therapy control unit. Refer to qualified service personnel.

Safety Information

Bed Frame - Always use a standard health care bed frame with side rails and any other safeguards or protocols that may be appropriate. Frame and side rails (if used - see **Risks and Precautions**) must be properly sized to eliminate any gaps that might entrap a patient's head or body.

Brakes - Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

Bed Height - To minimize the risk of falls or injury the bed should always be in the lowest position when the patient is unattended.

Side Rails / Patient Restraints - Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local competent authority or government agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patients frequently to guard against patient entrapment.

Skin Care - Monitor skin conditions regularly, especially at bony prominences and areas where moisture or incontinence may occur or collect, and consider adjunct or alternative therapies for high acuity patients. Early intervention may be essential to preventing serious skin breakdown.

Fluids - Avoid spilling fluids on the therapy control unit. If spills do occur, unplug the unit, clean fluid from the bed and the therapy control unit wearing rubber gloves to avoid any possibility of shock. Once fluid is removed, check operation of components in area of spill.

Fluids remaining on controls can cause corrosion, which may cause components to fail or operate erratically, possibly producing potential hazards for patient and caregivers.

Lock-Outs - The lock-out feature on the therapy control unit should be used at the caregiver's discretion to ensure against unintentional or unauthorized tampering with unit settings.

Avoid Fire Hazards - To minimize risk of fire, connect the unit's power cord directly to the wall-mounted outlet. Do not use extension cords or multiple outlet strips.

No Smoking in Bed - Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

Tobacco Smoke - If routine laundering and cleaning procedures as described in the **Preventative Maintenance** section of the **Maintenance** chapter of the *First Step Select MRS Maintenance Manual* are not followed, tobacco smoke build-up may conceivably reduce the amount of air flow through the product. Severe air restrictions (whatever the source) may cause the therapy control unit to overheat and automatically deactivate.

Power Cord - Ensure power cord is kept free from all pinch points and moving parts and is not trapped under casters. Improper handling of the power cord can cause damage to the cord, which may possibly produce risk of fire or electrical shock.

General Protocols - Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

Maximum Patient Weights

- *First Step* Classic Overlay: 250 lb (113 kg)
- *First Step* Cirrus Overlay: 300 lb (136 kg)
- *First Step* Excel Overlay: 300 lb (136 kg)
- *First Step* Air Rail Overlay: 250 lb (113 kg)

Connecting the System to Other Devices

All *First Step* Select Therapy System components are designed to be used as a single system. The *First Step* Select therapy control unit should only be used with the *First Step* Select products. Any attempt to connect and use the therapy control unit with any other brand or model of mattress or overlay may result in improper operation of equipment, possibly leading to increased risk of patient injury.

The therapy control unit can become hot during use. Prolonged contact with the therapy control unit may result in a hazard to the patient or caregiver.

Preparation For Use and Patient Placement

Preparing the *First Step* Select Therapy System for use entails inspecting the unit for any damaged or unserviceable assemblies as well as performing the steps outlined in the **Unboxing and Inspection** procedure below. It is recommended that all sections of this manual be reviewed prior to product use. Carefully read the **Contraindications, Safety Information and Risks and Precautions** sections in the **Introduction** chapter of this manual prior to placing a patient on the *First Step* Select Therapy System.


Unboxing and Inspection

1. Open shipping boxes. **Do not use sharp instruments to open boxes.** Damage to overlay, cover sheet or hoses could result.
2. Confirm that all components are included. See page 28 for inventory list.
3. Remove and discard all packing materials.
4. Inspect all components for any visible damage occurring from transport. Contact the manufacturer immediately if damage is observed.

Unit Setup and Preparation for Patient Placement

-  Always use a standard health care bed frame with such safeguards and protocols as may be appropriate.
1. Level bed and lock caster wheel brakes and caster wheel steering.
 2. Remove existing mattress if using Mattress Foam Base or V-Series Adaptor (VSA).
-  When selecting a mattress and overlay combination or MRS, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

3. Install therapeutic surface on standard Mattress Foam Base:
 - Place tote bag at foot of mattress.
 - Open tote bag and remove therapeutic surface, air hoses, therapy control unit, sheets, etc.
 - Place therapeutic surface with hose connections and tote bag at foot of mattress.
 - Ensure overlay is placed with elastic harnesses on underside.
 - For overlay, slide elastic harness over each end of mattress.
 - Tuck tote bag under foot end of therapeutic surface.
4. Install therapeutic surface on V-Series Adaptor
 - Unfold the V-Series Adaptor
 - Place onto VersaCare Frame
 - Secure First Step Cirrus Overlay to the V-Series Adaptor using attached straps
 - Install the locking bars on head and foot end

-  When removing the V-Series Adaptor, be sure to disengage locking bar and magnets before removing from bed frame.

-  Always ensure bed frame is properly sized. Ensure against any gaps that might entrap the patient's head or body.

-  Avoid contact of sharp instruments with the *First Step* Select Therapy System. Punctures, cuts and tears may prevent proper inflation and air pressure operation.

5. Install *First Step* Select Fabric Cover Sheet (optional for *First Step* Cirrus; *First Step* Cirrus may be used with a standard hospital fitted sheet):
 - Draw sides of sheet down each side of bed and tuck under mattress.
 - Secure sheet around hose connections using individual hook and loop strips located on cover sheet.
6. Install therapy control unit on footboard:
 - Swing out adjustable hanging arms located on each side of therapy control unit.

- Pull down tilt display bracket located on lower backside of therapy control unit.
- Hang therapy control unit on footboard.

i If the arms of the therapy control unit cannot accommodate a large footboard, a separate hanging bracket is available as an accessory, or the therapy control unit may be placed on a clean floor.

i Verify air intake located on the back face of the therapy control unit is not blocked by a blanket or other obstructions. The therapy control unit has a thermal cutoff that will, after a period of time, automatically deactivate if the air intake is blocked.

7. Connect hoses to therapy control unit:

- Locate hose set coupler that holds three hose tab ends together.
- Turn hose set coupler such that hoses are stacked vertically (Fig. 1).
- Insert upper tab of coupler into horizontal slot on the right side of therapy control unit (Fig. 2).

i Ensure hose set is not twisted.

- Pivot coupler downward, placing hose ends over air outlets, until bottom of coupler clicks into locked position (Fig. 3).

Fig. 1


Fig. 2



Fig. 3


8. Attach hoses to therapeutic surface air inlets:

- Match respective hose end colors to air inlet colors.
- Attach hoses snugly by turning them slightly as you push hose end onto the air inlets.

9. Plug therapy control unit power cord into wall outlet.

i Be sure power to this outlet is not controlled by a wall switch.

10. Press the power button. Therapy control unit will activate and the therapeutic surface will inflate to provide a firm surface for patient transfer. The home display at right will appear.

- Soft / firm display is the home display.

i Motor will ramp-up, setting pressures during first 25 seconds of operation. After ramping-up and setting pressures, unit will be noticeably quieter.

! Using the therapy control unit outside the recommended operating conditions could pose a hazard to patient or caregiver.

- For non-ambulatory patients - Place non-ambulatory patients on a firm surface.
- For ambulatory patients - Place ambulatory patients on a deflated surface.

11. Power down therapy control unit.

- Press green power button to deactivate therapy control unit.
- Unplug therapy control unit from wall outlet.

i All settings will be held in memory, and restored when therapy control unit is activated.

i Inspect all components for visible damage. Verify basic functionality and inspect all hoses, cords and other components for visible damage. Discontinue use and immediately contact manufacturer if damage is observed.

Patient Transfer to the *First Step Select Therapy System*

It is recommended that all sections of this manual be reviewed prior to product use. Carefully read the **Contraindications, Safety Information and Risks and Precautions** sections in the **Introduction** chapter of this manual prior to placing a patient on the *First Step Select Therapy System*.

Ambulatory Patient Transfer

1. If in operation, press green power button to turn OFF therapy control unit. Therapeutic surface will begin to deflate.
2. Adjust height of both surfaces to comfortable levels for patient transfer.
3. Lower side rails when applicable.
4. Lock caster wheel brakes and caster wheel steering on both surfaces.
5. Assist patient in ambulatory transfer, following all applicable safety rules and institution protocols.
6. Center patient on the surface.
7. Raise and lock side rails on both sides of bed, if decision has been made to use side rails (refer to the **Safety Information and Risks and Precautions** sections in the **Introduction** chapter of this manual).
8. Press green power button to activate therapy control unit and adjust air pressures for patient comfort. See the **Air Pressure Adjustment** section of this chapter.
9. Lower bed height to lowest level.
10. Activate warmer and adjust temperature, if desired.

Non-Ambulatory Patient Transfer

1. Press green power button on therapy control unit to activate blower.
2. Press the *InstaFlate™ Inflation Function (IF)* to fully inflate therapeutic surface and create a firm surface for patient transfer.
3. Perform steps 2-7 from **Ambulatory Patient Transfer** section.
4. Deactivate the *InstaFlate* IF and adjust air pressures for patient comfort. See **Air Pressure Adjustment** section of this chapter.
5. Perform steps 9-10 of **Ambulatory Patient Transfer** section.

Patient Transfer from the *First Step Select Therapy System*

Ambulatory Patient Transfer

1. Press green power button to turn OFF therapy control unit. Therapeutic surface will begin to deflate.
2. Adjust height of bed to same level as surface to which patient is being transferred.
3. Lock caster wheel brakes and caster wheel steering on both surfaces.
4. Lower side rails when applicable.
5. Transfer patient following all applicable safety rules and institution protocols.

Non-Ambulatory Patient Transfer


1. Press *InstaFlate* IF to fully inflate therapeutic surface and create a firm surface for patient transfer.
2. Perform steps 2-5 from **Ambulatory Patient Transfer** section.

Air Pressure Adjustment


Each therapeutic surface is constructed with three sections: HEAD, BODY, and LEG which can be independently adjusted to meet patient needs.

Use Pressure Adjust button to adjust pressure in each individual section of therapeutic surface. Air pressures are set by using either HT/WT pre-set or manual method.

1. Press PRESSURE ADJUST button shown at right.


2. Display shown at right will appear.

Pressure Adjust HT/WT Method

Use HT/WT pre-set method to set an appropriate air pressure based on patient's height and weight. HT/WT pre-sets can be adjusted from 2 ft 0 in (0.6 m) to 6 ft 8 in (1.8 mm 20 cm) by 1-inch (2.54 cm) increments and from 50 lb (22 kg) to 300 lb (136 kg) by 5 lb (2 kg) increments.




1. Press left arrow button under HT/WT. The display shown at right will appear.
 - Press left arrow button under DECR to decrease height.
 - Press right arrow button under INCR to increase height.
 - Press center oval button under ENTER to enter selected height. Display shown at right will appear.
 - Press left arrow button under DECR to decrease weight.
 - Press right arrow button under INCR to increase weight.
 - Press center oval button under ENTER to enter selected weight.
 - Press center oval under EXIT to return to Home Display.

Pressure Adjust Manual Method

Use manual method to make major adjustments in existing pressures.

1. With Pressure Adjust display showing, press right arrow button under MANUAL. Display shown at right will appear.
 - Press left arrow button to decrease head section pressure.
 - Press right arrow button to increase head section pressure.
 - Press center oval button to enter selected head section pressure. Display shown at right will appear.
 - Press left arrow button to decrease body section pressure.
 - Press right arrow button to increase body section pressure
 - Press center oval button to enter selected body section pressure. Display shown at right will appear.
 - Press left arrow button to decrease leg section pressure.
 - Press right arrow button to increase leg section pressure.
 - Press center oval button to enter selected leg section pressure. Display shown at right will appear.
 - Press center oval button under EXIT to return to Home Display.


If HT/WT display is accessed after adjustments have been made using MANUAL, the MANUAL pressures will automatically be replaced with previous HT/WT pressures.

Completion of Patient Placement

Home Display

1. Verify patient comfort.
2. Use Home Display for up to 25% uniform pressure adjustment.
 - Press left arrow button under SOFT to decrease uniform pressure.
 - Press right arrow button under FIRM to increase uniform pressure.
3. If the decision has been made to use side rails, raise and lock side rails.

! Serious injury or death can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints.

! **Side Rails / Patient Restraints** - Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. **Monitor patients frequently to guard against patient entrapment.**

i Adjust bed height to lowest level.


! To minimize risks of injury from falls, the bed frame should always be in the lowest possible position when the patient is unattended.

4. Ensure all brakes are locked.
5. Activate warmer and adjust temperature if desired.


InstaFlate Inflation Function (IF)

Use *InstaFlate* to assist patient transfer and bathing by increasing pressure throughout the therapeutic surface. *InstaFlate* will automatically deactivate after 15 minutes. Screen will flash every 5 seconds in the last 15 seconds before deactivation.

1. Press the *INSTAFLATE* ON / OFF button shown at right.

Display shown at right will appear and flash until *InstaFlate* is deactivated.


2. Press the *INSTAFLATE* ON / OFF button to deactivate *InstaFlate*. Display shown at right will appear temporarily and then return to the Home Display. Pressure and warmer settings will be restored when *InstaFlate* is deactivated.


Seat Deflate

Use seat deflate to assist patient exit and bedpan placement by lowering pressure in body section. Seat deflate will automatically deactivate after 15 minutes. Screen will flash every 5 seconds in the last 15 seconds before deactivation.


1. Press the SEAT DEFFLATE ON/OFF button shown at right.

Display shown at right will appear and flash until seat deflate is deactivated.


2. Press the SEAT DEFFLATE ON/OFF button to deactivate seat deflate. Display shown at right will appear temporarily and then return to the home display. Pressure and warmer settings will be restored when seat deflate is deactivated.

Alerts


Use *InstaFlate* to assist patient transfer and bathing by increasing pressure throughout the therapeutic surface. *InstaFlate* will automatically deactivate after 15 minutes. The *Instaflate* screen will flash every 5 seconds in the last 15 seconds before deactivation.

Display shown at right will appear and flash until *InstaFlate* is deactivated.

Use seat deflate to assist patient exit and bedpan placement by lowering pressure in body section. Seat deflate will automatically deactivate after 15 minutes. The seat deflate screen will flash every 5 seconds in the last 15 seconds before deactivation.


Display shown at right will appear and flash until seat deflate is deactivated.

Warmer Adjust


Use warmer adjust to set three optional levels of comforting heat when desired. Use only under direction of patient's physician.

1. Press WARMER ADJUST button shown at right.


The display shown at right will appear.

2. Press the right arrow button to activate warmer .

3. To deactivate the warmer, press and hold the left arrow button until the display reads WARMER OFF.

4. The panel displays a sliding scale that identifies the level of selected warmth.
5. Use the right and left arrow buttons to increase or decrease warmth.
6. As warmth is increased or decreased, LED readout displays LOW, MEDIUM or HIGH accordingly.


Lock-Out Procedures

Use lock-outs to disable buttons on the control panel to ensure against unauthorized tampering after air pressures and warmth have been adjusted to individual patient needs.


Lock-outs can be activated to disable either the pressure adjust button only, or all buttons on the control panel.

All Button Lock-Out


- With home display showing, press and hold the PRESSURE ADJUST button.


- The screen shown at the right will appear temporarily.



- Release PRESSURE ADJUST button when Lock-Out Activated display appears
- Display will return to home display

- To deactivate All Button Lock-Out press and hold the PRESSURE ADJUST button while each display shown at right appears temporarily.



- Release PRESSURE ADJUST button when Lock-Out Released display appears.
- Display will return to home display.



Pressure Adjust Button Lock-Out




- With home display showing, press and hold PRESSURE ADJUST button while each display shown at right appears temporarily.


- Press *INSTAFLATE* button when Lock-Out Activated display appears. Display shown at right appears temporarily.
- Display will return to home display.

- To deactivate Pressure Adjust Button Lock-Out, press and hold PRESSURE ADJUST button while each display shown at right appears temporarily.

- Press *INSTAFLATE* button when Lock-Out Activated display appears. Display shown at right appears temporarily indicating button released.
- Display will return to home display.



Nursing Care

It is recommended that all sections of this manual be reviewed prior to product use. **Carefully read the Contraindications, Safety Information and Risks and Precautions sections in the Introduction chapter of this manual prior to performing nursing care for a patient on the *First Step Select Therapy System*.**

Cardio Pulmonary Resuscitation (CPR)

1. Lift up on quick release lever (Fig. 1).
2. Swing hoses up and away from therapy control unit (Fig. 2).
3. Bring bed to a level position.
4. Lower side rails if used.
5. Place back board under patient.
6. Begin CPR. Patient's body weight and the CPR process will deflate therapeutic surface.
7. After CPR is performed and patient is clinically stable:
 - Connect hose set to therapy control unit.
 - Return bed to original configuration, considering all applicable safety rules and protocols. If the decision has been made to use side rails, restraints or other accessories (see related **Risks and Precautions** and **Safety Information** in the **Introduction** chapter of this manual), then raise and lock the side rails on both sides of the bed and reposition other restraints and accessories.

Patient Transfer to the *First Step Select Therapy System*

Refer to the **Patient Placement** chapter of this manual if patient is being transferred to the *First Step Select Therapy System* for the first time.

Patient Transfer While on *First Step Select Therapy System*

1. Press the green power button to deactivate therapy control unit. Therapeutic surface will begin to deflate.
 - (i) Patient will be supported by existing mattress or foam base if using the *First Step Select Therapy System* as a mattress replacement system.
 2. Unplug the power cord from wall outlet.
 3. If side rails are present, ensure they are raised and locked.
 4. Unlock caster wheel brakes and caster wheel steering.
 5. Transport patient following all applicable safety rules and institution protocols.
 6. Plug the power cord into a grounded wall outlet when destination is reached.
- (i) Be sure power to this outlet is not controlled by a wall switch.
7. When power is resumed, overlay will resume previous air pressure settings.
 8. Adjust head of bed for patient comfort and return bed to original configuration, considering all applicable safety rules and protocols. If the decision has been made to use side rails, restraints or other accessories (see related **Risks and Precautions** and **Safety Information** in the **Introduction** chapter of this manual), then raise and lock side rails on both sides of the bed and reposition other restraints and accessories.

Patient Bathing

1. When side rails are available, raise the rail on the side of the bed opposite the caregiver to help reduce the risk of patient fall during bathing.
2. Level bed and adjust bed height to facilitate bathing.
3. Press the *InstaFlate* button to provide a firm surface for patient bathing.
4. Bathe and rinse patient.

- (i) Avoid spilling fluids on bed controls. If spills do occur, clean fluid from surface wearing rubber gloves to avoid any possibility of shock. Once fluid is removed, check operation of components in area of spill.

Fluids remaining on electronic controls can cause corrosion, which may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and caregiver.

To reduce risk of patient falls during bathing, caregiver must never leave patient unattended or turn away during bathing.

5. Cover sheet or therapeutic surface can be cleaned during patient bathing.

Replace and launder extensively soiled cover sheet or therapeutic surface.

6. Dry cover sheet or therapeutic surface with a towel. Remaining moisture will normally evaporate rapidly with air flow.
7. Return patient to supine position in center of surface.
8. Adjust warmer and bed height as required.
9. Adjust head of bed for patient comfort and return bed to original configuration, considering all applicable safety rules and protocols. If the decision has been made to use side rails, restraints or other accessories (see related **Risks and Precautions** and **Safety Information** in the **Introduction** chapter of this manual), then raise and lock side rails on both sides of the bed and reposition other restraints and accessories.

Bedpan Placement

1. Level bed.
2. Activate SEAT DEFLATE.
3. Lower side rail (if used) on caregiver's side.
4. Log roll patient toward opposite side rail.
5. Position bedpan parallel to patient's buttocks, next to patient and push bedpan down into overlay.
6. Log roll patient onto bedpan, keeping one hand on bedpan.
7. Adjust head of bed for patient comfort and return bed to original configuration, considering all applicable safety rules and protocols. If the decision has been made to use side rails, restraints or other accessories (see related **Risks and Precautions** and **Safety Information** in the **Introduction** chapter of this manual), then raise and lock side rails on both sides of the bed and reposition other restraints and accessories.

Bedpan Removal

1. Level bed.
2. Lower side rail (if used) on caregiver's side.
3. Grasp bedpan firmly with one hand, and turn patient toward opposite side with other hand.
4. Remove bedpan and cleanse patient's posterior.
5. Return patient to supine position.
6. Press SEAT DEFLATE button to deactivate seat deflate and return body cushion to previous setting.
7. Center patient on surface, with patient's buttocks centered on surface.
8. Adjust head of bed for patient comfort and return bed to original configuration, considering all applicable safety rules and protocols. If the decision has been made to use side rails, restraints or other accessories (see related **Risks and Precautions** and **Safety Information** in the **Introduction** chapter of this manual), then raise and lock side rails on both sides of the bed and reposition other restraints and accessories.

Repositioning Patient using Dri-Flo™ Pads

1. Adjust bed to lowest height position.
2. Press *InstaFlate* button to activate *InstaFlate* and create a firm surface for patient repositioning.
3. Lower side rails (if used).
4. Take positions on each side of patient (at least two people will be needed). Roll edges of *Dri-Flo* Pad on each side of patient to get a firm grip.
5. Use *Dri-Flo* Pad to slightly lift and properly position patient.
6. Press *InstaFlate* button to deactivate *InstaFlate* and return surface to previous air pressure settings.

7. Adjust head of bed for patient comfort and return bed to original configuration, considering all applicable safety rules and protocols. If the decision has been made to use side rails, restraints or other accessories (see related **Risks and Precautions and Safety Information** in the **Introduction** chapter of this manual), then raise and lock side rails on both sides of the bed and reposition other restraints and accessories.

Incontinence / Drainage

1. A *First Step* Select cover sheet or moisture-impermeable, breathable underpads should be used for incontinent patients.
 *Dri-Flo* Pad breathable underpads are recommended for incontinent patients. They are absorbent and do not block air flow. Do not use plastic-backed underpads. Plastic may tend to block moisture vapor transmission and air flow from mattress surface.
2. Watch for incontinence or drainage and provide appropriate skin care following each episode.

Skin Care

1. Remove excess moisture and keep skin dry and clean.
2. Maintain proper air pressures in accordance with the **Air Pressure Adjustment** section of the **Patient Placement** chapter of this manual.
3. Assess patient's skin regularly, particularly in areas where incontinence and drainage occur.

Monitor skin conditions regularly, especially at bony prominences and areas where moisture or incontinence may occur or collect, and consider adjunct or alternative therapies for high acuity patients. Early intervention may be essential to preventing serious skin breakdown.

Care and Cleaning

The following are the manufacturer's recommended daily and weekly cleaning and infection control procedures for the *First Step* Select Therapy System while in use. If this product is to be purchased, or utilized for long term rentals, it is recommended that the manufacturer be contacted for recommended infection control procedures to be utilized by the facility.

It is recommended that all sections of this guide be reviewed prior to product use. Carefully read the **Contraindications**, **Safety Information** and **Risks and Precautions** sections in the **Introduction** chapter of this manual prior to performing care and cleaning procedures on the *First Step* Select Therapeutic Surface.

Daily Care and Cleaning of *First Step* Select Therapy System While in Use

- 
- The patient does not need to be removed from the therapeutic surface when performing daily cleaning procedures. The patient may be bathed as the top surface is being cleaned.
1. For daily care and cleaning, refer to the procedures in the **Patient Bathing** section of the **Nursing Care** chapter in this manual.

Weekly Care and Cleaning of *First Step* Select Therapeutic Surface While in Use

In addition to daily care and cleaning, the cover sheet should be laundered or replaced weekly, or more often as needed. Excessively soiled therapy system and cover sheets should be replaced and cleaned as described in the **Infection Control Procedures** of this chapter.

Disassembly of *First Step* Select Therapeutic Surface and Cover Sheet for Laundering and Cleaning

1. Press power button to deactivate therapy control unit.
2. Unplug therapy control unit from wall outlet.
3. Disconnect hose set from therapy control unit and therapeutic surface.
4. Remove cover sheet (if used).

Infection Control Procedures

Follow all regulations outlined in OSHA Occupational Exposure to Blood borne Pathogens 29 CFR Part 1910.30 regarding infection control procedures to be used when cleaning the *First Step* Select Therapy System.

Cleaning (Single Patient Use)

i When blood and / or body fluids are present, wash fabric items separate from family laundry and consult a health care professional for proper use of the following throughout all cleaning procedures:

- Latex or rubber gloves.
- Plastic apron or gown.
- Protective eyewear.

i Do not use petroleum based solvents or alcohol.

i Chlorine bleach may fade fabric colors. Consult Arjo if there are specific infection control concerns or other issues regarding the use of bleach.

Infection Control Protocol

The *First Step Select Therapy System* must be disinfected after each patient use. Always use universal precautions, treating all used equipment as contaminated. When performing cleaning and disinfecting procedures, it is recommended that disposable gloves, an impervious gown or apron and protective eye wear be worn.

Laundering Procedure for Cover Sheet

Handle soiled items as little as possible, with minimum agitation to prevent gross microbial contamination of the air.

1. Wear protective clothing as described.
2. Launder cover sheet in washing machine according to the following instructions:
 - Spray heavily soiled areas with stain-removing prewash spray. Follow container directions for use.
 - Set machine temperature on warm wash, using the longest wash cycle. Select a cool or cold water rinse.
 - Add measured high-quality laundry detergent and allow machine to fill with water. (Follow directions on detergent box for specific washer size to follow washing machine manufacturer's recommendation).
 - If chlorine bleach has been recommended by a health care professional, add 8 oz. of high-quality chlorine bleach to the laundry soap and water.

i Do not allow concentrated bleach to contact fabric as it may result in discoloration.

- Add soiled cover sheet.

i Do not overload washing machine.

3. Wipe the top of the washing machine, lid surfaces and inside the rim with a disinfectant to ensure clean laundry is not recontaminated.
4. Remove protective clothing if used. Discard or clean according to instructions from health care professional.
5. Wash hands before handling clean laundry.
6. Remove cover sheet from the washer after full cycle is completed.
7. Place cover sheet into dryer.
8. Tumble dry on warm or low heat. **Do not use high heat.**

i Keep dryer temperature below 140° F (60° C).

i All items must be completely dry before removing them from the dryer.

9. Remove cover sheet from dryer and store in a clean area.

Laundering Procedure for Overlays, Mattress and Tote Bag

Overlays should be laundered monthly (or more often, if necessary). Wear recommended personal protective equipment.

i For the *First Step* Cirrus Overlay, plastic caps must be placed over the air inlets to prevent accumulation of water inside during laundering.

1. Set washing machine for warm wash - cold rinse cycle.
2. Spray heavily soiled areas with stain-removing prewash spray.
3. Add a powdered, high-quality laundry detergent measured in accordance with manufacturer's directions and approximately 1 cup bleach to wash water (1/2 oz of bleach per gallon of water).
4. Allow washing machine to agitate for two minutes before adding product.

i Do not overload washing machine. Overloading may prevent thorough cleaning. Do not wring mattress or overlays. Allow water to drain out.

i For the *First Step* Cirrus Overlay, before drying, the plastic caps must be removed from air inlets to allow any moisture to escape during tumble or air drying.

5. Tumble dry using low heat or no heat dryer setting, or allow to air dry. A high heat setting may damage mattresses, overlays and cover sheets.

Cleaning Procedure for Therapy Control Unit and Hoses

The therapy control unit and hoses should be cleaned weekly (or more often, if necessary). Wear recommended personal protective equipment.

1. Unplug therapy control unit from wall outlet.
2. Scrub all surfaces and hoses with a damp, coarse cloth or brush using a germicide solution.
- i** Do not allow excess solution to penetrate the therapy control unit.
3. Remove all visible soiling.
4. Wipe off excess germicide solution.
5. Wipe all surfaces with a damp cloth containing disinfectant to provide a final disinfection.
6. Allow to air dry completely.

Cleaning Procedure for Foam Base

The foam base should be cleaned weekly (or more often, if necessary). Wear recommended personal protective equipment.

1. Wet a clean, coarse cloth with germicide solution and squeeze out excess solution.
2. Lightly scrub foam base with wet cloth.
- i** If foam base is too heavily soiled, remove and replace with a new base.
3. Allow foam base to air dry completely.

Cleaning and Disinfection for V-Series Adaptor

The V-Series Adaptor should be cleaned regularly after each patient use. It is recommended that the V-Series Adaptor is removed prior to cleaning and disinfection (see page 11 for instructions on removal).

Cleaning: The surface of the V-Series Adaptor should be cleaned with a mild soap or organic solution and water prior to using standard disinfection solutions in order to remove gross contaminate.


Disinfection: All disinfection solutions must be properly diluted according to manufacturer's instructions. Please follow standard institutional cleaning and disinfection procedures.

1. The V-Series Adaptor can be disinfected using a wipe down approach with a 70% alcohol solution or bleach 1000 ppm.*
2. Rinse with a wet towel to remove chemicals.
3. After wipe down, dry surface with a towel.

Cleaning Procedure for Air Filter

The air filter should be cleaned weekly (or more often, if necessary). Wear recommended personal protective equipment.

1. Remove air filter from rear of therapy control unit by pinching center of filter and pulling outward.
2. Vacuum filter to remove dust before washing.
3. Soak filter in germicide solution for 10 minutes.
4. Squeeze out excess solution from filter.
5. Allow filter to air dry completely.
6. If air filter is ripped, torn or damaged, contact manufacturer for replacement part immediately.

Institutional Cleaning (DMEs and Extended Care Facilities)

1. Follow all listed procedures for In-Home Cleaning, with the following modifications:
 - Always follow Universal Precautions.
- i** Handle soiled items as little as possible, with minimum agitation to prevent gross microbial contamination of the air and personnel.
- Always add high-quality chlorine bleach (5.25% sodium hypochlorite) to laundry soap and water. Add 8 oz per normal size wash load or 12 oz per wash load for larger loads (this applies to standard household size washers using approximately 16 gallons of water).
 - Disinfectants used for scrubbing surfaces should be a disinfectant solution made of EPA-registered high-quality disinfectant detergent mixed according to manufacturer's instructions.
 - After scrubbing surfaces with disinfectant, wipe surfaces with a damp cloth containing chlorine solution mixed at 2 oz of high-quality chlorine bleach per gallon of water to provide a final disinfection.
2. After cleaning and drying completely, cover the product with clean plastic and store in clean area.

Assembly of the First Step Select Therapy System after Laundering and Cleaning

For re-assembly after laundering refer to **Unit Setup** in the **Preparation For Use** and **Patient Placement** section of this manual.

***Chlorine concentrations may vary from 250 ppm to 10,000 ppm depending on local policy and contamination status. If an alternative disinfectant is selected from the wide variety available, Arjo recommends that suitability for use is confirmed with the chemical supplier prior to use.**

Maintenance

This chapter contains preventative maintenance and unscheduled maintenance for the *First Step Select Therapy System*, as well as pertinent warnings. Information in this chapter should be used in conjunction with information provided in the **Troubleshooting** chapter of this manual. Certain functions may require the use of a trained professional depending on the degree of difficulty and level of knowledge or skill. Please exercise sound judgment when attempting certain complex maintenance functions and follow all standard safety protocols.

The following warnings should be observed during **all maintenance procedures**:

High voltage equipment to be serviced by trained, qualified personnel only. Dangerous voltages can present a significant shock or burn hazard to personnel.

Unauthorized modification of the *First Step Select therapy control unit* could result in hazard to patient / caregiver.

Preventative Maintenance

Preventative maintenance to be performed while unit is in use includes the following:

- **Daily Cleaning** - the cover sheet should be wiped daily with a mild soap and water solution. This can be done when the patient is bathed, using the bath water.
- **Weekly Cleaning** - the cover sheet should be laundered at least weekly, depending on soiling. The therapy control unit and hoses should be cleaned weekly. The air filter should be cleaned at least weekly. A clogged filter may result in damage to the therapy control unit.
- **Monthly Cleaning** - the mattress or overlay should be laundered at least monthly.

Preventative maintenance to be performed after each patient use includes the following:

- Inspect mattress or overlay, cover sheets, therapy control unit filter and hoses for wear or damage. Unserviceable items should be replaced with an identical item.
- Perform infection control procedures as defined in the **Infection Control Procedures** section of the **Care and Cleaning** chapter of this manual after each patient use and otherwise as outlined above.

Unscheduled Maintenance

Power must be switched off and the mains power cord unplugged prior to any maintenance procedure.

The *First Step Select Therapy Control Unit* is made up of several components which should not be opened, disassembled or otherwise modified by the user and should be replaced as a unit. All assemblies, operations, adjustments, modifications, maintenance and repairs must be carried out by qualified personnel authorized by the manufacturer. For additional information or for specific questions, please contact the manufacturer's Customer Service or refer to the **Questions and Information** section of this manual.

Troubleshooting

Troubleshooting procedures for the *First Step Select Therapy System* are provided in the charts on the following page. Each chart deals with a specific symptom and provides a **Problem / Possible Cause / Remedy** approach to use in identifying the solution. The **Problem** will describe the condition of the unit. The **Possible Cause** will describe several likely reasons for the symptom and the steps to take in verifying the cause. The **Remedy** will describe the final solution to the symptom and its cause.

It is recommended that trained, qualified personnel perform any modifications, repairs or parts replacement. Please contact the manufacturer or refer to the **Questions and Information** section.

If this product is not in good working order, for patient safety, discontinue use and make sure patient is promptly placed on a suitable alternative surface.

The following warnings should be observed during **all troubleshooting procedures**:

High voltage equipment to be serviced by trained, qualified personnel only. Dangerous voltages can present a significant shock or burn hazard to personnel.

Mains power must be unplugged prior to any maintenance procedure.

Carefully observe and mark locations and routing of all cables and tubing removed while performing the procedure in order to properly locate connections, cables and their tie-wraps during replacement. Wiring and tubing must be reinstalled exactly as before removal to ensure proper unit operation.

Troubleshooting Chart for the *First Step Select Therapy System*.

Problem	Cause	Remedy
Therapy control unit will not activate. Power light does not illuminate.	Power cord may be unplugged from the wall outlet. Therapy control unit may be plugged into a non-functioning outlet or an outlet controlled by a wall switch. Circuit breaker or fuse for the wall outlet may be blown. Therapy control unit may be damaged.	Verify power cord is properly plugged in, both at the wall and on the therapy control unit. Relocate to functioning wall outlet or outlet not controlled by a wall switch. Reset circuit breaker. Return therapy control unit for repair.
Therapy control unit is ON but therapeutic surface does not inflate.	Air hose may not be properly connected to therapy control unit. Air hose may not be properly connected to mattress or overlay The therapeutic surface may have rips, tears or punctures.	Check to ensure air hoses are properly connected to therapy control unit. Check to ensure air hoses are properly connected to mattress or overlay. Replace therapeutic surface.
A single section of therapeutic surface does not inflate.	The therapeutic surface may have rips, tears or punctures. Hose may be kinked, or not properly connected. A hose may be damaged.	Check for rips, tears or punctures and replace therapeutic surface if necessary. Ensure hoses are free of kinks and properly connected. Inspect and replace hose.

Troubleshooting Chart for the *First Step Select Therapy System*

Problem	Cause	Remedy
Therapy control unit is ON, therapeutic surface is inflated but not receiving warm air.	Heater not operating properly.	Return therapy control unit for repair.
Therapeutic surface feels warmer than expected.	Warmer is on too high or warmer is not operating properly.	Turn warmer OFF. Recheck after 10 minutes. If therapeutic surface is still too warm, return therapy control unit for repair.
Therapeutic surface does not respond when <i>InstaFlate</i> is activated.	Hoses kinked or not properly connected. Therapy control unit may be damaged.	Ensure hoses are free of kinks and properly connected. Return therapy control unit for repair.
Therapeutic surface does not respond when seat deflate is activated.	Hoses kinked or not properly connected. Therapy control unit may be damaged.	Ensure hoses are free of kinks and properly connected. Return therapy control unit for repair.
Therapeutic surface does not respond when comfort adjust is activated.	Hoses kinked or not properly connected. Therapy control unit may be damaged.	Ensure hoses are free of kinks and properly connected. Return therapy control unit for repair.
Therapy control unit does not respond to adjustment from control buttons.	Therapy control unit may be damaged.	Return therapy control unit for repair.

Replacement Parts and Inventory List

The following items are the replacement parts and literature available for the repair and maintenance of the *First Step Select Therapy System*. Refer to the drawings on the following page for identification and location of replacement parts.

Product Components

Description	Inventory List	Part No.
<i>First Step Select Therapy Control Unit (115V) with Accessories</i>	<i>First Step Select Therapy Control Unit (115V)</i> <i>First Step Select Hose Set</i> <i>First Step Select Coversheets (qty 2)</i> <i>First Step Select User Manual</i> <i>First Step Select Quick Reference Guide</i>	215200
<i>First Step Select Hose Set</i>	<i>First Step Select Hose Set</i>	201729
<i>First Step Classic Overlay</i>	<i>First Step Classic Overlay</i>	203604
<i>First Step Cirrus Overlay</i>	<i>First Step Cirrus Overlay</i>	415159
<i>First Step Excel MRS</i>	<i>First Step Excel MRS</i>	203660
<i>First Step Air Rail Overlay</i>	<i>First Step Air Rail Overlay</i>	203670
<i>First Step Select Foam Base (standard)</i>	<i>First Step Select Foam Base 35 in x 82 in (88.9 cm x 208.3 cm)</i>	MS2200
<i>First Step Select Cover Sheet</i>	<i>First Step Select Cover Sheets (qty 2)</i>	203605
<i>First Step Cirrus Tote Bag</i>	Tote Bag, <i>First Step Cirrus</i>	416102
Air Filter, <i>First Step Select Therapy Control Unit</i>	Air Filter, <i>First Step Select Therapy Control Unit</i>	200778
V-Series Adaptor	V-Series Adaptor	AH1017309

Product Literature


Description	Part No.	Quantity
<i>First Step Select Therapy System Quick Reference Guide (60601 3rd)</i>	416005	1
<i>First Step Select Therapy System User Manual (60601 3rd) English, French, Spanish</i>	416001	1
<i>First Step Select Therapy System Maintenance Manual (60601) 115V</i>	215239-AH	1

Accessories


Description	Part No.	Quantity
Foot Support	226540	1
Hanging Bracket, Footboard	203640	1

Miscellaneous


Description	Part No.	Quantity
0.10 Beaded Cable Tie-Wrap (attaches QRG to Therapy Control Unit)	20648	1


First Step Classic Therapeutic Overlay - a 5 in low profile overlay that fits on top of a foam pad (PN 203620) or an existing mattress.

First Step Cirrus Overlay - a 5 in low profile RF welded overlay that fits on top of a foam pad (PN 203620) V-series Adaptor (PN AH1017309) or an existing mattress.

First Step Excel MRS - a 7 in overlay that fits on top of a foam pad (PN 203620) or an existing mattress of thickness no greater than 3 in.

First Step Air Rail Overlay - a 5 in low profile overlay with air rails extending 1 to 1.5 in above the overlay and that fits on top of a foam pad (PN 203620) or an existing mattress.

Specifications

Specifications subject to change without notice.

Dimensions:

Therapeutic Surfaces (inflated):

First Step Classic and Air Rail	82 in L x 35 in W x 5 in H / 208 cm L x 89 cm W x 12.7 cm H
First Step Cirrus	84 in L x 35 in W x 5 in H / 213 cm L x 89 cm W x 12.7 cm H
First Step Excel	82 in L x 35 in W x 7 in H / 208 cm L x 89 cm W x 18 cm H
Mattress Foam Base	82 in L x 35 in W x 3 in H / 208 cm L x 89 cm W x 7.6 cm H
V-Series Adaptor	84 in L x 35 in W x 5.5 in H / 213 cm L x 89 cm W x 14 cm H

 MRS foam base is optional with all overlays.

Therapy Unit

115 VAC	12 in W x 10-3/4 in H x 7 in D / 30.5 cm L x 27.3 cm W x 18 cm H
---------	--

Weight:

Therapeutic Surfaces:

First Step Classic and Air Rail	3 lb (1.3 kg)
First Step Cirrus	7.5 lb (3.4 kg)
First Step Excel	4 lb (1.8 kg)
First Step Foam Base (standard)	5 lb (2.3 kg)
V-Series Adaptor (fits Hill-Rom VersaCare™ Med Surg) 19 lb (8.6 kg)	

Maximum Recommended Patient Weights

First Step Classic	250 lb (113 kg)
First Step Cirrus	300 lb (136 kg)
First Step Air Rail	250 lb (113 kg)
First Step Excel	300 lb (136 kg)

Therapy Unit

115 VAC	14 lb (6.4 kg)
---------	----------------

Electrical:

Voltage	115 VAC
Frequency	60 Hz
Amps	6 Amps
Max Electrical Leakage	100 Micro Amps
Fuse	F1 - 1.0A, 250 VAC, Time Lag F2 - 6.3A, 250 VAC, Time Lag
Fuse Breaking Capacity	F1 - 50 Amps @ 115 VAC F2 - 50 Amps @ 115 VAC
Power Cord Length	16 ft (4.9 m)

Symbols Used			
	Alternating Current		Type BF Applied Part
	Warning of possible electrical shock hazard		Warning of possible hazard to system, patient or staff Important operational information
	Unplugged prior to any maintenance procedure		This product is designated for separate collection at an appropriate collection point. For home health care, do not dispose of as household waste.
	Carefully observe and mark locations and routing of all cables and tubing removed		Reorder Number
	Manufacturer		Double Insulated
	Manufactured Date		ETL Listed, Conforms to AAMI ES60601-1 1st Edition, CSA C22.2#60601-1 3rd Edition and IEC 60601-1 3rd Edition
	Refer to Manual		Equipment not suitable for use in the presence of a Flammable Anesthetic Mixture with Air or Oxygen or Nitrous Oxide
	Important Operational Information	IP21	Ingress protection
	Consult Instructions for Use		Temperature Low and High Limits
	Humidity limitations		Do Not Open with Scissors
	Wipe Down Only		Keep dry
	Tumble Dry Temperature Limits 60° C		Do Not Iron
	Wash Temperature Limits 60° C		No Hooks
	Safe Working Load		No Phenol
	Use Solution diluted to 1000ppm of available chlorine 1000ppm NaOCl NaOCC		Serial Number
	Do Not rotate mattress		Do Not flip mattress

Miscellaneous

EMC

First Step Select Therapy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The *First Step Select Therapy System* was tested to the following standard.

- IEC 60601-1-2 Issued: 2007/03/30 Ed:3 Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests

Classification

Classification in accordance with UL 60601-1, CSA C22.2 No 601.1 and IEC. 60601-1

- Class II (UL 60601-1)
- Type BF Applied Part
- No Sterilization
- Not for Use with Flammable Anesthetic Mixture with Air, Oxygen or Nitrous Oxide.
- Continuous Operation
- Ordinary Equipment

Environmental Storage / Transit Conditions

Temperature Range: -20° to 60° C (-4° to 140° F)

Relative Humidity Range: 10% to 90%

Atmospheric Pressure Range: 500 hPa to 1060 hPa

Environmental Operating Conditions

Temperature Range: 10° to 40° C (50° to 104° F)

Relative Humidity Range: 10% to 95%

Atmospheric Pressure Range: 500 hPa to 1060 hPa

Leakage Testing and Ground Continuity Testing

This product is double insulated (Class II). It is constructed with two completely independent electrical insulation systems. The power cord is not used to ground the unit. Grounding is not required. The unit shall not be modified to ground the unit.

Arjo products certified as Double Insulated (Class II) meet the requirements of the Medical Electrical Equipment Standards (UL60601-1, IEC60601-1, and EN60601-1) and are labeled accordingly. Each box of the designation label (see **Symbols Used**) symbolizes the two independent insulation systems.

Leakage Testing

Double Insulated electrical devices (normally ungrounded apparatus) must be tested differently than other electrical devices. Per the Medical Electrical Equipment Standard (UL60601-1), the manufacturer, Underwriters Laboratories Inc. and other test organizations, all parties interested in testing for current leakage shall use the Methods identified for Double Insulation (Class II). Refer to electrical standards UL60601-1, IEC 60601-1 or EN 60601-1, Clause 19 and ANSI/AAMI ES1:1993.

Ground Continuity Testing

Double Insulated products (normally ungrounded apparatus) rely on two independent electrical insulation systems that are isolated from metal parts. Ground continuity is not required.

Questions and Information

For questions regarding this product, supplies, maintenance or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative, or see www.Arjo.com.

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AUSTRALIA
Arjo Australia Pty Ltd
78, Forsyth Street
O'Connor
AU-6163 Western Australia
Tel: +61 89337 4111
Free: +1 800 072 040
Fax: +61 89337 9077

BELGIQUE / BELGIË
Arjo NV/SA
Evenbroekveld 16
BE-9420 ERPE-MERE
Tél/Tel: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.be

BRASIL
Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelice, 329 PB02 Galpão - Lapa
São Paulo – SP – Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA
Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSIONAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strzi 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK
Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA
Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/Marie Curie 5
Edificio Alfa Planta 6 oficina 6.1.-62
ES-28521 Rivas Vacia, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com

FRANCE
Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG
Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA
Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST
Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building,North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND
Arjo BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND
Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@arjo.com

NORGE
Arjo Norway AS
Olaf Helssets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

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Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

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Arjo Inc.