Flowtron TRP30

	<div> <div>≤ 71</div> <div>cm</div> </div>
	<div> <div>≤ 28</div> <div>inch</div> </div>

EN To be used only under the direction of a physician • Non-Sterile • For Single Patient Use Only • Not made with natural rubber latex • For use with Flowtron® DVT-prevention pumps only

Caution: (Applicable to the USA market only) Federal law restricts this device to sale by or on the order of a licensed practitioner.

Description

The *Flowtron* Active Compression System is comprised of a Tri Pulse garment connected to a Tri Pulse enabled *Flowtron* DVT prevention pump manufactured by ArjoHuntleigh.


Instructions For Use

Note: Refer to the relevant Tri Pulse enabled *Flowtron* DVT prevention pump Instructions For Use, for complete information on the use of the system.

- Remove the garments from the sealed bag and record the Lot Number in the patient notes.
- The garments may be used on either leg. Unfold a garment and position the inflatable bladder directly behind the patient's calf as indicated on the garment (Fig. A).

Note: The arrow on the underside of the garment must point to the heel.

- Wrap the garment around the patient's leg (Fig. A). Then starting at the ankle and working upwards, secure each fastener tab in turn (Fig. B), ensuring that the entire garment fits snugly. Repeat with the other leg.

-  When placed on leg arrow must point to heel.

- Attach the garments to the Tri Pulse enabled *Flowtron* DVT prevention pump tubeset, ensuring a "click" is heard from the pump connector.
- Power up the pump.
- Check the icons on the pump display to confirm that the correct type of garment has been connected to each leg.
- Press "Start" to begin therapy; the indicators should be illuminated green.
- Press "Stop" to end therapy.

Indications

The intended use of the Tri Pulse garment is to help prevent Deep Vein Thrombosis (DVT). The risk factors of thrombus formation (DVT) are present during orthopaedic, trauma, urologic and general surgery. Intermittent Pneumatic Compression (IPC) is also suitable for other patient groups including neurology, critical care, general medical and obstetrics.

General Recommendations

- Check that there are no kinks in the pump tubeset and garment tubing.
- Regularly check that the garments remain correctly fitted to the patient.

Contraindications

IPC should not be used in the following conditions:

- Severe arteriosclerosis or other ischaemic vascular diseases.
- Severe congestive cardiac failure or any condition where an increase of fluid to the heart may be detrimental.
- Known or suspected acute Deep Vein Thrombosis (DVT) or phlebitis.
- Any local condition in which the garments would interfere, including gangrene, recent skin graft, dermatitis or untreated, infected leg wounds.

Note: If you are unsure of any contraindications refer to the patient's physician before using the device.

Cautions

- This product cannot be adequately cleaned and / or sterilized by the user in order to facilitate safe reuse and is therefore intended for single patient use.
- Attempts to clean or sterilize these devices may result in a biocompatibility, infection or product failure risk to the patient.
- Garments should be positioned to prevent sustained pressure points on the skin, paying particular attention to patients who are unconscious, cannot feel or have reduced sensation and/or ability to move their leg(s).
- The patient's skin should be inspected frequently during every shift.
- Clinical judgement is required to determine if the patient's skin condition requires additional protective measures, or if the therapy should be discontinued.
- Garments should be removed immediately if the patient experiences tingling, numbness or pain.
- When used for DVT prevention, continuous intermittent pneumatic compression is recommended until the patient is fully ambulatory. Uninterrupted use of the system is recommended.
- The system should be used **WITH CAUTION** on patients with insensitive extremities, diabetes, impaired circulation, or fragile or impaired skin.

Note: These are guidelines only and should not replace clinical judgement and experience.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor.

In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

End of Life Disposal

- Garment material or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.

© and ™ are trademarks belonging to the Arjo group of companies. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

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