Flowtron FG100

EN To be used only under the direction of a physician • Non-Sterile • For Single Patient Use Only • Not made with natural rubber latex


Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.

Description
For use with Flowtron® Deep Vein Thrombosis (DVT)-Prevention pumps manufactured by ArjoHuntleigh. Must not be used with the Flowtron® Hydroven 3 or Flowtron® Hydroven 12 Intermittent Pneumatic Compression (IPC) pumps, or the Flowtron® Excel DVT pump which is for calf/thigh garments only.

More comprehensive information on the garment can be found in the relevant Flowtron DVT pump Instructions For Use document.

Instructions For Use

1. Refer to the following Flowtron DVT pump IFU for complete information on setting up and operating the system.
2. Remove the garments from the sealed bags and apply to the feet as shown. The garments may be used on either foot.
3. The foot garments are fitted as follows:
 - (a) Place the foot in the centre of garment. Ensure the back of the garment is in line with the heel as indicated.
 - (b) Bring flap (1) over the top of the foot and hold in place.
 - (c) Bring flap (2) over the foot and secure.
 - (d) Bring flap (3) over the foot and secure. The fit should be snug but comfortable.
 - (e) Bring strap (4) around the back of the heel and fix in place as indicated.
 - (f) Tension strap so garment is secure and comfortable.

4. Attach the garment to the pump tubset, ensuring that a "click" is heard with each snap-lock connector. Ensure that the tubing is attached securely to each connector.
5. Turn the pump on. The pump will perform a short self-test cycle and will automatically adjust to the default settings of 130 mmHg pressure, 3 seconds inflation hold and 30 seconds cycle time. Check the garments after a few inflations and, if necessary, re-adjust for comfort and security of fit.

Indications

- Clinical applications for the foot garments are as follows:
- Prevention of Deep Vein Thrombosis (DVT).
 - Enhancement of venous & arterial circulation.
 - Prevention of venous stasis.
 - Assistance in the healing of cutaneous ulcers, including venous ulcers.
 - Reduction of acute & chronic edema.
 - Reduction of lower limb pain due to trauma or surgery.
 - Reduction of compartment pressures.

Recommendations

- General Recommendations:
- Garments should be removed regularly to inspect the skin.
 - Where appropriate, patients should be instructed in the proper use of the system, the purpose of therapy and any problems that should be reported to the nursing staff.

DVT prophylaxis:

- The garments should be applied to the patient pre-operatively, prior to the induction of anaesthesia.
- The system should be used continuously for no less than 72 hours postoperatively or until the patient becomes fully ambulatory.
- If the garment cannot be applied to the operative limb during surgery, it may be applied to the limb once the patient reaches the recovery unit.
- In the non-surgical patient, the system should be initiated immediately once the risk of DVT formation is identified.

Contraindications

- The foot garments should not be used in the following conditions:
- Severe arteriosclerosis or other ischaemic vascular diseases.
 - Known or suspected acute Deep Vein Thrombosis (DVT) or phlebitis.
 - Severe congestive cardiac failure or any condition where an increase of fluid to the heart may be detrimental.
 - Pulmonary embolism.
 - Any local condition in which the garments would interfere, including gangrene, recent skin graft, dermatitis or untreated, infected leg wounds.

If you are unsure of any contraindications refer to the patient's physician before using the device.

Cautions

- Garments should be removed immediately if the patient experiences tingling, numbness or pain.
- When used for DVT prophylaxis, continuous use is recommended and any interruption of therapy for a substantial length of time should be at the discretion of the physician.
- Patients should be instructed not stand or walk with the foot garments on.
- This product cannot be adequately cleaned and / or sterilized by the user in order to facilitate safe reuse and is therefore intended for single patient use. Attempts to clean or sterilize these devices may result in a biocompatibility, infection or product failure risk to the patient.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

End of Life Disposal

- Garment material or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.

® and ™ are trademarks belonging to the Arjo group of companies. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

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© e ™ sono marchi di proprietà del gruppo di società Arjo. La politica dell'azienda è volta al miglioramento costante, per tale ragione questa si riserva il diritto di modificare i prodotti senza preavviso. È vietato riprodurre integralmente o parzialmente il contenuto della presente pubblicazione senza espresso consenso di Arjo.