
INSTRUCTIONS FOR USE

Stretcher Slings

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使用说明书 · Bedienungsanleitung · Mode d'emploi · Istruzioni per l'uso · Gebruiksaanwijzing

Contents

Foreword	2	Loop Attachment & Detachment.....	7
Customer Support.....	2	Attach the Loops (5 Steps).....	7
Definitions in this IFU	2	Detach the Loops (2 steps).....	7
Intended Use	3	Applying the Sling.....	7
Patient/Resident Assessment	3	In Bed (21 steps).....	7
Expected Service Life	3	On Floor (28 steps)	8
Safety Instructions	4	Removing the Sling	10
Safety Practises	4	In Bed (11 steps)	10
Parts Designation	5	Cleaning and Disinfection.....	11
Preparations	5	Cleaning Instruction (4 steps)	11
Before First Use (6 steps)	5	Flites – Do NOT wash	11
Before Every Use (7 steps)	5	Cleaning Chemicals	11
After every use	5	Disinfection.....	11
Sling Selection.....	6	Care and Preventive Maintenance	12
Transfer/Transport/Toiling.....	6	Before and After Every Use.....	12
Bathing	6	When Soiled, Stained and Between Patients.....	12
Allowed Combinations	6	Storage.....	12
Safe Working Load (SWL).....	6	Service and Maintenance.....	12
		Troubleshooting	13
		Technical Specifications	13
		Label on the Sling.....	14

EN

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product. Mandatory to read the Instructions for Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies. © Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

The content of this publication may not be copied either whole or in part without the consent of Arjo.

Foreword

Thank you for purchasing Arjo equipment. Please read this *Instructions For Use (IFU)* thoroughly!

Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this *IFU*.

Definitions in this IFU

WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Means: This is important information for the correct use of this system or equipment.

Intended Use

Stretcher Sling is specifically designed to transfer a patient/resident in a supine position. This product is intended for assisted lateral transfer of patients/residents with limited ability to move. The *Stretcher Sling* should be used together with Arjo lift devices and Arjo stretcher frames in accordance with the Allowed Combinations specified in the Instruction for Use (*IFU*). Product MAA1152 made of Mesh material is also intended for bathing.

Stretcher Flites is intended for assisted lateral transfer of patients/residents with limited ability to move while using a patient/resident transfer lifter, under the supervision of trained nursing staff. *Stretcher Flites* are to be used for a limited period only, and, by nature of their design, must be treated as a disposable and patient/resident specific product. *Stretcher Flites* should be used together with Arjo lift devices in accordance with the Allowed Combinations specified in the *IFU*.

Stretcher Flites must not be used to support a resident while bathing or showering.

Stretcher Flites are disposable product and cannot be washed.

The *Stretcher Sling/Stretcher Flites* shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, and in accordance with the instructions outlined in the Instructions for Use.

The *Stretcher Sling/Stretcher Flites* is intended to be used in hospital environments, nursing homes, other healthcare facilities and home care.

The *Stretcher Sling/Stretcher Flites* should only be used for the purpose specified in this Instructions for Use. Any other use is prohibited.

Patient/Resident Assessment

We recommend that facilities establish regular assessment routines. Caregivers should assess each patient/resident according to the following criteria prior to use:

- Patient/resident is passive
- Might be almost completely bed ridden
- Often stiff or contracted joints
- Totally dependent on the caregiver
- Physically demanding for the caregiver
- Stimulation and activation is not a primary goal

EN

Safe Working Load (SWL) for *Stretcher Sling*:

- MAA1151, MAA1150, MAA1152: 272 kg (600lbs).

Safe Working Load (SWL) for *Stretcher Flites*:

- MFA5000: 272 kg.

The right type and size of slings should be used after proper assessment of each resident's size, condition and the type of lifting situation.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

Expected Service Life

The expected service life of *Stretcher Sling/Stretcher Flites* is the maximum period of useful life.

The expected life of the sling is dependent on the actual use conditions. Therefore, before use, always make sure that the sling does not show any sign of fraying, tearing or other damage and that there is no damage (i.e. cracking, bending, breaking). If any such damage is observed, do not use the sling.

The expected service life of the *Stretcher Sling/Stretcher Flites*:

Service Time:

- 2 years (*Stretcher Sling*)
- Flites: 2 weeks (based on 4 transfers per day) (*Stretcher Flites*)

Shelf Time:

- 5 years (*Stretcher Sling/Stretcher Flites*)

Safety Instructions

WARNING

To avoid injury, make sure that the patient is not left unattended at any time.

WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch on fire.

WARNING

To avoid injury, never wash a Flites. A Flites is only intended for “single patient use”.

EN

WARNING

To avoid injury, always assess the patient prior to use.

WARNING

To avoid injury, only store the equipment for a short period of time. If stored longer than stated in the *IFU*, it can cause weakness and breakage in the material.

WARNING

To avoid falling, ensure that the user weight is lower than the safe working load for all products or accessories being used.

WARNING

To avoid injury, keep the equipment away from sun/UV-light. Exposure to sun/UV-light can weaken the material.


Safety Practises

- At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.
- Patient with spasm can be lifted, but great care should be taken to support the patient’s legs.


Serious incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Parts Designation

1. Sling label (located on the outside of the sling)
2. Do not wash label (On Flites only, located on the outside of the sling)
3. Strap with loops

4. Commode hole (certain slings only)
5. Foot end
6. Head end

EN

Example of slings. Not all models are represented.

Preparations

Before First Use (6 steps)

1. Check all parts of the sling, see section “*Parts Designation*” on page 5. If any part is missing or damaged - Do NOT use the sling.
2. Read this *IFU* thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the *IFU* should be kept and easily accessible at all times.
5. Make sure to have a rescue plan ready in case of a patient emergency.
6. For questions, contact your local Arjo representative for support.

NOTE

Arjo recommends the name of the patient is written on the sling care label to avoid cross contamination between patients.

Before Every Use (7 steps)

1. Patient assessment must always determine the methods used in everyday care, see “*Intended Use*” on page 3.

WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the *IFU*.

2. Check the allowed combinations against the sling, spreader bar and lift, see “*Allowed Combinations*” on page 6.

WARNING

To avoid injury, always make sure to inspect the equipment prior to use.

3. Check all parts of the sling, see “*Parts Designation*” on page 5. If any part is missing or damaged - Do NOT use the sling. Check for:
 - fraying
 - loose stitching
 - thread breakage
 - tears
 - fabric holes
 - soiled fabric
 - damaged loops
 - unreadable or damaged label

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this *IFU*.

4. Check that the sling is clean. If the sling is not clean see “*Cleaning and Disinfection*” on page 11.
5. Always check the washing symbol on the Flites label. If the Flites has been washed, throw the Flites away. See “*Flites – Do NOT wash*” on page 11.
6.  If this symbol appears, do NOT use the Flites. The Flites label indicates it has been washed.
7. Before using the sling, read the lift *IFU* for transfer, transport and bathing.

After every use

If the patient is leaving the hospital, throw away the patient specific Flites.

Sling Selection

The patient's physical disabilities, weight distribution and general physique needs to be taken into consideration when selecting a sling.

Select the *Stretcher Sling* size by measuring against the patient's head and feet. Make sure to select the right size to avoid the patient's feet from falling out of the sling.

Transfer/Transport/Toileting

Article number	Safe Working Load	Product Description	Size
MAA1151	272 kg (600 lbs)	Soft Stretcher without commode hole	L (1610 mm), XL (1650 mm)
MAA1150	272 kg (600 lbs)	Soft Stretcher with commode hole for toileting	L (1610 mm), XL (1650 mm)
MFA5000	272 kg (600 lbs)	Flites Stretcher (Disposable and Patient Specific)	One size (1740 mm)

Bathing

Article number	Safe Working Load	Product Description	Size
MAA1152	272 kg (600 lbs)	Soft Stretcher Mesh with commode hole	L (1610 mm), XL (1650 mm)

Allowed Combinations

WARNING

To avoid injury, always follow the allowed combinations listed in this *IFU*. No other combinations are allowed.

Safe Working Load (SWL)

Always follow the lowest SWL of the total system. E.g. the *Maxi Move* lift has a SWL of 160 kg (352 lbs) and the MAA1151 sling has a SWL of 272 kg (599 lbs). This means that the *Maxi Move lift* has the lowest SWL. The patient is not allowed to weigh more than the lowest SWL.

Lift	Maxi Move	Maxi Sky 600	Maxi Sky 2
SWL	160 kg (352 lbs)	272 kg (600 lbs)	272 kg (600 lbs)
Spreader Bar	Stretcher Frame 700-19302	Stretcher Frame 700-19522	Stretcher Frame 700-15695
Sling	SWL	Size	Size
MAA1150	272kg (600 lbs)	L, XL	L, XL
MAA1151	272kg (600 lbs)	L, XL	L, XL
MAA1152	272kg (600 lbs)	L, XL	L, XL
MFA5000	272kg (600 lbs)	One Size	One Size

Loop Attachment & Detachment

Attach the Loops (5 Steps)

1. Place the loop over the spring loaded latch (A). (See Fig. 1)
2. Pull the loop down to force the latch to open.
3. Make sure that the spring loaded latch closes completely with the loop inside. (See Fig. 1)
4. Make sure that the latch is moving freely.
5. Make sure loops and straps are not twisted.

TIPS

Use the shorter loop lengths for less swing in the sling.

EN

Detach the Loops (2 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the loop.

Method 1

1. Press the spring loaded latch open.
2. Remove the loop.

Method 2

1. Pull one side of the loop over the hook and spring loaded latch. (See Fig. 2)
2. Pull down the loop. (See Fig. 3)

Applying the Sling

In Bed (21 steps)

1. Place the lift nearby.
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. The patient needs to lay down flat on his/her back.
4. Fold the sling lengthwise along the centre line, outside facing outside. The sling label is located on the outside.

WARNING

To avoid injury during log roll, make sure there is another caregiver present and/or use the side rails on the opposite side of the bed.

5. Log roll the patient. (Use local technique.) If log roll is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding sheet/tube IFU*
6. Place the folded sling over the patient's side.
7. Fold down the top of the sling and tuck it under the patient's body. If a stretcher commode sling is




used make sure the commode hole is positioned correctly. (See Fig. 4)

8. Log roll the patient to the other side.
9. Pull out the remaining part of the sling from under the patients body.
10. Return the patient to his/her back.
11. Make sure that the:
 - sling is centred and flat without creases,
 - sling's head end covers the patient's head area,
 - sling pieces are not twisted underneath the patient and
 - if a stretcher commode sling is used make sure the commode hole to positioned correctly

WARNING

To avoid injury, make sure the patient's arms are placed inside of the sling.

Continue with the steps on the next page.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

EN

12. Position the stretcher bar over the patient. (See Fig. 6)

13. Start by attaching the loops at each side of the patient's hip area. Make sure the mid part of the sling covers the patient's hip area for better comfort. (See Fig. 1)

TIPS

Use the shorter loop lengths for less swing in the sling.

14. Finally attach the remaining loops starting from the head area. Make sure to attach each loop at the same loop length as the one across.

15. If medical condition permits, raise the patient's head for better comfort by using a different loop length.

16. Slightly lift the patient to create tension in the sling. (See Fig. 7)

WARNING

To avoid the patient from falling, make sure that the sling attachments are securely attached before and during the lifting process.

17. Make sure that:

- All loops are securely attached
- The patient lays comfortably in the sling.

18. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the loops.

19. Adjust the balance of the stretcher frame, if needed. See respective Lift IFU.
20. If a stretcher commode sling is used, a bed pan can be placed under the patient, if needed.
21. Transfer the patient. See respective lift IFU.

On Floor (28 steps)

1. Make a clinical assessment of the patient. Do not use the sling if the patient has a head, neck, spine or hip injury. (Follow local routine.) Use a scoop stretcher instead. See the *Scoop Stretcher IFU*.
2. Place a pillow under the patient's head. (See Fig. 8)
3. Make sure the patient is in an area that is accessible for the lift. If needed, move the patient by using an Arjo sliding sheet/tube. See respective *Sliding sheet/tube IFU*.
4. Fold the sling lengthwise along the centre line, outside facing outside. The sling label is located on the outside.
5. Log roll the patient. (Use local technique.) If log roll is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding sheet/tube IFU*.
6. Place the folded sling over the patient's side. Make sure that the centre line is aligned with the patient's spine. (See Fig. 9)
7. Fold down the top of the sling and tuck it under the patient's body.
8. Log roll the patient to the other side.
9. Pull out the remaining part of the sling from under the patients body.

Continue with the steps on the next page.

Fig. 5

Fig. 6


Fig. 7


Fig. 8


Fig. 9

10. Return the patient to his/her back.
11. Make sure that the: (See Fig. 10)
 - sling is centred and flat without creases,
 - sling's head end covers the patient's head area,
 - sling pieces are not twisted underneath the patient and
 - if a stretcher commode sling is used make sure the commode hole is positioned correctly
12. Position the lift sideways with the lift legs wide open.
13. One caregiver needs to be close to the patient's head and the other caregiver needs to lift up the patient's legs.
14. Position one lift leg next to the patient's head and the other lift leg under the patient's legs. Make sure the stretcher frame is positioned over the patient's shoulders. (See Fig. 11)
15. Apply the brakes on the lift.
16. Lower the stretcher frame.
17. Start by attaching the loops at each side of the patient's hip area. Make sure the mid part of the sling covers the patient's hip area for better comfort. (See Fig. 1)
18. Finally attach the remaining loops starting from the head area. Make sure to attach each loop at the same loop length as the one across.
19. If medical condition permits, raise the patient's head for better comfort by using a different loop length.
20. Slightly lift the patient to create tension in the sling. (See Fig. 12)

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

EN

WARNING

To avoid injury, make sure the patient's arms are placed inside of the sling.

12. Position the lift sideways with the lift legs wide open.
13. One caregiver needs to be close to the patient's head and the other caregiver needs to lift up the patient's legs.
14. Position one lift leg next to the patient's head and the other lift leg under the patient's legs. Make sure the stretcher frame is positioned over the patient's shoulders. (See Fig. 11)
15. Apply the brakes on the lift.
16. Lower the stretcher frame.
17. Start by attaching the loops at each side of the patient's hip area. Make sure the mid part of the sling covers the patient's hip area for better comfort. (See Fig. 1)
18. Finally attach the remaining loops starting from the head area. Make sure to attach each loop at the same loop length as the one across.
19. If medical condition permits, raise the patient's head for better comfort by using a different loop length.
20. Slightly lift the patient to create tension in the sling. (See Fig. 12)
21. Make sure that:
 - All loops are securely attached
 - The patient lays comfortably in the sling.
22. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the loops.
23. Adjust the balance of the stretcher frame, if needed. See respective Lift *IFU*.
24. Continue to lift up the patient.
25. Remove the pillow.
26. Release the brakes on the lift.
27. Bring the lift legs together.
28. Transport the patient. See respective lift *IFU*.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

16. Lower the stretcher frame.
17. Start by attaching the loops at each side of the patient's hip area. Make sure the mid part of the sling covers the patient's hip area for better comfort. (See Fig. 1)

TIPS

Use the shorter loop lengths for less swing in the sling.

18. Finally attach the remaining loops starting from the head area. Make sure to attach each loop at the same loop length as the one across.
19. If medical condition permits, raise the patient's head for better comfort by using a different loop length.


Fig. 10

Fig. 11

Fig. 12

Removing the Sling

In Bed (11 steps)

1. Position the patient over the bed. (See Fig. 13)
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. Make sure the bed is in a flat position.

WARNING

To avoid injury during log roll, make sure there is another caregiver present and/or use the side rails on the opposite side of the bed.

EN

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

4. Lower the patient onto the bed.
5. Detach the sling, make sure all straps are detached. (See Fig. 14)
6. Move the lift away from the patient.

7. Log roll the patient. (Use local technique.) If log roll is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding sheet/tube IFU*
8. Tuck the sling under the patient's body.
9. Log roll the patient to the other side and remove the sling.
10. Return the patient to his/her back.
11. Make sure the patient is in a comfortable and secure position in the bed.

Fig. 13



Fig. 14

Cleaning and Disinfection

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

WARNING

To avoid material damage and injury, clean and disinfect according to this IFU.

- No other chemicals are allowed.
- Never clean with chlorine.
- Chlorine will deteriorate the surface of the material.

WARNING

To avoid injury, always remove the sling before disinfecting the lift.

WARNING

To avoid injury, never wash a Flites. A Flites is only intended for "single patient use".

All Stretcher slings, except the Flites, should be cleaned when they are soiled or stained and between patients.

Cleaning Instruction (4 steps)

(NOT Flites)

1. Detach the sling from the lift.
2. Check the symbols on the product label for cleaning recommendations.
3. Machine wash the sling at 70°C (158°F). In order to obtain sufficient level of disinfection, please wash according to local hygiene legislation.
If allowed according to the product label, tumble dry in low temperature, max 60°C (140°F).

4. Do NOT

- wash with rough surfaces or sharp objects
- use any mechanical pressure, pressing or rolling
- use bleach
- use gas sterilization
- use autoclave
- dry clean
- steam
- ironing

EN

Flites – Do NOT wash

The Flites is for single patient use.

Do NOT clean, wash, disinfect, wipe or sterilize the Flites. If the Flites has been subjected to any such treatment, it shall be discarded.

The Flites is marked with the "Do NOT wash" label. (See Fig. 15)

If the Flites is washed, the "Do NOT use/lift" label will appear. (See Fig. 16)

Cleaning Chemicals

(NOT Flites)

Use common commercial detergent without optical brightener.

No other chemicals are allowed e.g. chlorine, softener, iodine based disinfectants, bromine and ozone.

Disinfection

The only allowed disinfection is disinfection by washing.

This does not apply to the Flites. Do NOT clean, wash or disinfect the Flites.

Fig. 15

Fig. 16

Care and Preventive Maintenance

WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

EN

Before and After Every Use

Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching
- Tears
- Holes
- Discolouration and stains from bleaching
- Sling soiled or stained
- Unreadable or damaged label
- For Flites: Make sure the “Do NOT wash” symbol on the head support is intact. If damaged or missing and the underlying red “Do NOT use/lift” symbol appears, replace the sling immediately. See “*Label on the Sling*” on page 14.

When Soiled, Stained and Between Patients

Clean/Disinfect, All slings except the Flites

The caregiver should make sure that the sling is cleaned according to “*Cleaning and Disinfection*” on page 11, when it is soiled or stained and between patients.

Flites

Washing or disinfection is not allowed. Do not wipe, sterilize or soil the Flites. If the Flites has been subjected to any treatment, it shall be discarded.

Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or to excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 month), according to ISO 10535.

Troubleshooting

Problem	Action
The patient is not positioned correctly in the sling.	<ul style="list-style-type: none"> • Make sure that both shoulder loops are attached at the right length. • Make sure the straps are not twisted. • Make sure the patient is positioned on the inside of the sling. The outside has a label for reference. • Make sure the patient is centred in the sling.
The patient feels discomfort when using the sling.	<ul style="list-style-type: none"> • Make sure that there are no creases on the fabric. • Make sure that both shoulder loops are attached at the right length.
It is difficult to use the sling in combination with the lift.	Use the sling only in accordance with the combinations described in section “ <i>Allowed Combinations</i> ” on page 6.
The position of the patient in the sling should be more reclined or upward.	Use the loop lengths to adjust the patient’s position.
The patient wishes to have extra support for the neck.	Place a pillow behind the patient’s neck, inside the sling.
For MFA5000: The Flites is soiled, stained or wet.	Dispose Flites and replace it with a new one.
For MFA5000: The “Do NOT wash” label is damaged or missing and the underlying red symbol appears	Dispose Flites and replace it with a new one.

EN

Technical Specifications

General	
Safe working load (SWL) = Maximum total load	See “ <i>Allowed Combinations</i> ” on page 6.
Service life – Recommended period of use	Flites: 2 weeks, based on 4 transfers/day* All other slings: 2 years* *Refer to “ <i>Expected Service Life</i> ” on page 3.
Shelf life – Maximum period of storing new unpacked product	5 years
Model and Type	See “ <i>Sling Selection</i> ” on page 6

Operating, Transport and Storage Environment	
Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (-13°F to +140°F) Transport
Humidity	15–70% at +20°C (68°F) Operating and Storage 10–95% at +20°C (68°F) Transport

End of Life Disposal	
Package	The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.

Label on the Sling

EN

Care and Washing Symbols		Symbol for Loop
	Machine wash 70°C (158°F)	
	No bleaching allowed	
	No tumble drying allowed	
	Tumble drying Does not apply to Flites	
	No Ironing allowed	
	Do NOT wash the Flites Only applies to Flites	
	No dry cleaning allowed	
Certificates/Markings		Article Number
	CE marking indicating conformity with European Community harmonised legislation	REF XXXXXX-X Article number with -X refers to the sling size.
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745	REF XXXXXXX If the article number is stated without the size letter at the end, it means the sling is a One Size.
Flites Symbols		Fibre Content
	Only on Flites. (Disposable and patient specific)	PES Polyester
	"Do NOT wash" symbol. Located on the outside of the sling.	PU Polyurethane
	"Do NOT use/lift" symbol. The Flites has been washed. Located on the outside of the sling.	PP Polypropylene
Misc. Symbols		
		Safe Working Load (SWL)
		Patient's name symbol
		Record symbol
		Read the <i>IFU</i> before use
		Article number
		Manufacturing date and year
		Manufacturer name and address

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Intentionally left blank

Intentionally left blank

AUSTRALIA
Arjo Australia Pty Ltd
78, Forsyth Street
O'Connor
AU-6163 Western Australia
Tel: +61 89337 4111
Free: +1 800 072 040
Fax: +61 89337 9077

BELGIQUE / BELGIË
Arjo NV/SA
Evenbroekveld 16
BE-9420 ERPE-MERE
Tél/Tel: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.be

BRASIL
Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelice, 329 PB02
Galpão - Lapa
São Paulo – SP – Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA
Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSIONAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strzi 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK
Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA
Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/Marie Curie 5
Edificio Alfa Planta 6 oficina 6.1.-62
ES-28521 Rivas Vacia, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com

FRANCE
Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG
Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA
Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST
Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building,North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND
Arjo BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND
Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE
Arjo Norway AS
Olaf Helsets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

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At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges.

With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

04.SS.00-INT1

CE

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Arjo Inc.