
INSTRUCTIONS FOR USE

Paediatric Slings

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使用说明书 · Bedienungsanleitung · Mode d'emploi · Istruzioni per l'uso · Gebruiksaanwijzing

Contents

Foreword	2	Detach the Clips (3 steps)	7
Customer Support	2	Applying the Sling	8
Definitions in this IFU	2	In Bed (22 steps)	8
Intended Use	3	In Chair/ Wheelchair (17 steps)	8
Patient/Resident Assessment	3	On Floor (31 steps)	9
Expected Service Life	3	Removing the Sling	11
Safety Instructions	4	In Bed (12 steps)	11
Safety Practises	4	In Chair/Wheelchair (13 steps)	11
Serious incident	4	Cleaning and Disinfection	12
Parts Designation	5	Cleaning Instruction (5 steps)	12
Preparations	5	Cleaning Chemicals	12
Before First Use (6 steps)	5	Disinfection	12
Before Every Use (6 steps)	5	Care and Preventive Maintenance	12
Select Sling Size	6	Before and After Every Use	12
Using the Arjo Measuring Tape		When Soiled or Stained and Between Patients	12
(4 steps) (Accessory)	6	Storage	12
Sizing Without a Measuring Tape (2 steps)	6	Service and Maintenance	12
Sling Selection	6	Troubleshooting	13
Allowed Combinations	7	Technical Specifications	14
Safe Working Load (SWL)	7	Label on the Sling	14
Clip Attachment & Detachment	7		
Attach the Clips (5 steps)	7		

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product. Mandatory to read the Instructions for Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Foreword

Thank you for purchasing Arjo equipment. Please read this *Instructions For Use (IFU)* thoroughly!

Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this *IFU*.

Definitions in this IFU

WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Means: This is important information for the correct use of this system or equipment.

Intended Use

The *Paediatric Sling* is a product intended for assisted transfer of paediatric patients/residents with limited ability to move and require full body support during transfers. The *Paediatric Sling* shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, and in accordance with the instructions outlined in this *Instructions for Use*.

Paediatric Sling is intended for children. The *Paediatric Sling* is intended to be used in hospital environments, nursing homes, other healthcare facilities and home care. The *Paediatric Sling* should only be used for the purpose specified in this *Instructions for Use (IFU)*. Any other use is prohibited.

Patient/Resident Assessment

We recommend that facilities establish regular assessment routines. Caregivers should assess each patient/resident according to the following criteria prior to use:

- Patient/resident sits in wheelchair
- Has no capacity to support herself/himself at all
- Cannot stand unsupported and is not able to bear weight, not even partially
- Dependent on caregiver in most situations
- Physically demanding for caregiver
- Stimulation of remaining abilities is very important

Or, for a patient/resident who:

- Is passive
- Might be almost completely bed ridden
- Often stiff or contracted joints
- Totally dependent on the caregiver
- Physically demanding for caregiver
- Stimulation and activation is not a primary goal

Safe Working Load (SWL) for *Paediatric Sling*:

- MAA4020M: 125 kg (275 lbs).

The right type and size of slings should be used after proper assessment of each patient/resident's size, condition and the type of lifting situation.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

Expected Service Life

The expected service life of *Paediatric Sling* is the maximum period of useful life.

The operational life of the sling is dependent on the actual use conditions. Therefore, before use, always make sure that the sling does not show signs of fraying, tearing or other damage and that there is no damage (i.e. cracking, bending, breaking). If any such damage is observed, do not use the sling. If you have any doubts about sling safety, as a precaution and to ensure safety, do not use the sling.

The expected service life of the *Paediatric Sling*:

- Service Time: 2 years
- Shelf Time: 5 years

EN

Safety Instructions

EN

WARNING

To avoid injury, make sure that the patient is not left unattended at any time.

WARNING

To avoid injury, always assess the patient prior to use.

WARNING

To avoid falling, make sure that the user weight is lower than the safe working load for all products or accessories being used.

WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch on fire.

WARNING

To avoid injury, only store the equipment for a short period of time. If stored longer than stated in the *IFU*, it can cause weakness and breakage in the material.

WARNING

To avoid injury, keep the equipment away from sun/UV-light. Exposure to sun/UV-light can weaken the material.


Safety Practises

- At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.
- Patient with spasm can be lifted, but great care should be taken to support the patient's legs.

Serious incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Parts Designation

1. Stiffener pocket
2. Sling label (located on the outside of the sling)
3. Buckle
4. Attachment clip
5. Leg strap
6. Leg flap
7. Body section
8. Shoulder strap
9. Head section

Example of slings. Not all models are represented.

EN

Preparations

Before First Use (6 steps)

1. Check all parts of the sling, see section “Parts Designation” on page 5. If any part is missing or damaged, Do NOT use the sling.
2. Read this IFU thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the IFU should be kept and easily accessible at all times.
5. Make sure to have a rescue plan ready in case of a patient emergency.
6. For questions, contact your local Arjo representative for support.

NOTE

Arjo recommends the name of the patient is written on the sling care label to avoid cross contamination between patients.

Before Every Use (6 steps)

1. Patient assessment must always determine the methods used in everyday care, see “Intended Use” on page 3.

WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the IFU.

2. Make sure to select correct type and size of sling, see “Sling Selection” on page 6.

3. Check the allowed combinations against the sling, spreader bar and lift, see “Allowed Combinations” on page 7.

WARNING

To avoid injury, always make sure to inspect the equipment prior to use.

4. Check all parts of the sling, see “Parts Designation” on page 5. If any part is missing or damaged - Do NOT use the sling. Check for:
 - fraying
 - loose stitching
 - tears
 - fabric holes
 - soiled fabric
 - damaged clips
 - unreadable or damaged label

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

5. Check that the sling is clean. If the sling is not clean see “Cleaning and Disinfection” on page 12.
6. Before using the sling, read the lift IFU for transfer/transport.

Select Sling Size

Using the Arjo Measuring Tape (4 steps) (Accessory)

The Arjo Measuring Tape is a sizing accessory only intended as a guide and is an approximation. Patient width and body proportions should also be considered when selecting the correct size of Arjo sling.

1. If possible, the patient should be in a seated position.
2. Place the point marked **A** on the measuring tape over the patient's coccyx/against seat (position **A**). (See Fig. 1) and (See Fig. 2)

3. Measure from coccyx/seat (position **A**) to the top of the head (position **B**) (See Fig. 1).
4. The coloured area on the measuring tape that is levelled with the top of the patient's head indicates the required sling size. If the sling size falls in between two sizes, it is recommended to select the smaller size. (See Fig. 2)

Sizing Without a Measuring Tape (2 steps)

1. Place the sling over the patient's back.
2. Make sure the sling covers the patient from the top of the head (position **B**) to the coccyx (position **A**). (See Fig. 1)



Fig. 1

Fig. 2

Sling Selection

The patient's physical disabilities, weight distribution and general physique needs to be taken into consideration when selecting a sling.

Transfer

Article number	Safe Working Load	Product Description	Size
MAA4020M	125 kg (275 lbs)	4 Point clip sling Padded	XXS, XS, S

Allowed Combinations

WARNING

To avoid injury, always follow the allowed combinations listed in this IFU. No other combinations are allowed.

Safe Working Load (SWL)

The SWL of the sling is at 125kg (275 lbs). Always follow the lowest SWL of the total system.
 E.g. the *Maxi Twin Compact* lift/spreader bar has a SWL of 160 kg (352 lbs) and the MAA4020M sling has a SWL of 125kg (275 lbs). This means that the sling has the lowest SWL. The patient is not allowed to weigh more than the lowest SWL.

EN

	Lift	Maxi Twin Compact	Maxi Lite	Maxi 500	Maxi Twin	
	SWL	160 kg (352 lbs)	160 kg (352 lbs)	227 (500 lbs)	182 kg (401 lbs)	
	Spreader Bar (DPS)	Power Medium	Flat	Manual	Manual Medium	Power
Sling	SWL	Size	Size	Size	Size	Size
MAA4020M	125kg (275 lbs)	S	S	S	S	S

	Lift	Maxi Move			Maxi Sky	
	SWL	227kg (500 lbs)			272 kg (600 lbs)	
	Spreader Bar (DPS)	Manual Small	Manual Medium	Power Medium	Manual Flat	Power Medium
Sling	SWL	Size	Size	Size	Size	Size
MAA4020M	125kg (275 lbs)	XXS, XS, S	S	S	S	S

Clip Attachment & Detachment

NOTE

If available on your lift device, reposition the spreader bar using the Dynamic Position System (DPS) for easier clip attachment of the sling. See lift IFU for instructions.

Attach the Clips (5 steps)


- Place the clip on the spreader bar lug. (See Fig. 3)
- Pull the strap down.
- Make sure the lug is locked at the top end of the clip. (See Fig. 4)
- Make sure the strap is not squeezed in between the

- clip and the spreader bar.
- Make sure the straps are not twisted.

Detach the Clips (3 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the clips.

- Pull the strap up. (See Fig. 5)
- Make sure the lug is unlocked at the bottom of the clip.
- Remove the clip.

Applying the Sling

In Bed (22 steps)

For transfer/transport/bathing, see respective lift *IFU*.

1. Place the lift nearby.
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. The patient needs to lay down flat on his/her back.
4. Check that the stiffeners are completely inside the stiffener pockets, if any.
5. Fold the sling lengthwise along the centre line, outside facing outside. The sling label is located on the outside.

EN

WARNING

To avoid injury during log roll, make sure there is another caregiver present and/or use the side rails on the opposite side of the bed.

6. Log roll the patient. (Use local technique.) If log roll is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding Sheet/Tube IFU*
7. Place the folded sling over the patient's side. Make sure that the centre of the sling is aligned with the patient's spine, starting at the coccyx. Place the leg flaps towards the legs. **(See Fig. 6)**
8. Fold down the top of the sling and tuck it under the patient's body.
9. Log roll the patient to the other side.
10. Pull out the remaining part of the sling from under the patient's body.
11. Return the patient to his/her back.
12. Place the leg flaps underneath the patient's leg. Place a pillow under the patient's knees to ease the attachment process of the leg clips.
13. Make sure that the: **(See Fig. 7)**
 - sling is centred and flat without creases,
 - sling's head support covers the neck/head area,
 - sling pieces are not twisted underneath the patient and
 - all buckles are safely closed, if any

WARNING

To avoid injury, make sure the patient's arms are placed inside of the sling.

14. Make sure the patient's arms are inside the sling.
15. Place the lift next to the patient and apply the brakes.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

16. Attach the sling. **(See Fig. 4)**
17. For easier attachment, make sure the spreader bar is in a reclined position and if possible adjust the bed to seated position.
18. Slightly lift the patient to create tension in the sling. **(See Fig. 8)**

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

19. Make sure that all clips are attached securely and the patient lays comfortably in the sling.
20. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip.
21. Release the brakes on the lift.
22. Transfer and transport according to the lift *IFU*.

In Chair/ Wheelchair (17 steps)

1. Place the lift nearby.
2. Apply the brakes on the wheelchair.
3. Check that the stiffeners are completely inside the stiffener pockets, if any.
4. Stand in front of the patient and lean the patient forward.

Continue with the steps on the next page.

Fig. 6



Fig. 7

Fig. 8

5. Place the sling over the patient's back and head with the inside of the sling against the patient. The sling label is located on the outside.
6. Make sure that the centre of the sling is aligned with the patient's spine and coccyx.
7. Tuck in the end of the sling at the coccyx or against the seat. If needed, use an Arjo sliding sheet/tube for placing the sling under the patient. See respective *Sliding sheet/tube IFU*
8. Lean the patient back.
9. Place the leg flaps underneath the patient's legs. **(See Fig. 9)**
10. Make sure that the: **(See Fig. 10)**
 - sling is centred and flat without creases,
 - sling's head support covers the neck/head area,
 - sling pieces are not twisted underneath the patient
 - all buckles are safely closed, if any
 - patient's arms are inside the sling.

WARNING

To avoid injury, make sure the patient's arms are placed inside of the sling.

11. Make sure the spreader bar is in a seated position and place the lift in front of the patient with the lift legs wide open.
12. Apply the brakes on the lift.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

13. Attach the sling. If possible adjust the bed or spreader bar to a seated position (for easier attachment).

WARNING

Make sure straps are not caught by wheelchair or lift castors.

14. Slightly lift the patient to create tension in the sling. **(See Fig. 11)**

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

15. Make sure that all clips are attached securely and the patient lays comfortably in the sling. . If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip. **(See Fig. 4)**
16. Release the brakes on the lift.
17. Transfer/transport the patient according to the lift *IFU*.

On Floor (31 steps)

1. Make a clinical assessment of the patient. Do not use the sling if the patient has a head, neck, spine or hip injury. Use a scoop stretcher instead. See the Scoop Stretcher *IFU*.
2. Place a pillow under the patient's head. **(See Fig. 12)**
3. Make sure the patient is in an area that is accessible for the lift. If needed, move the patient by using an Arjo sliding sheet/tube. See respective *Sliding sheet/tube IFU*.
4. Check that the stiffeners are completely inside the stiffener pockets, if any.
5. **If the patient is able to sit up, continue to next step. If not able to sit up, proceed to step 9.**
6. **Patient able to sit up:** Stand behind the patient and help him/her to a seated position.
7. Place the sling over the patient's back and head with the inside of the sling against the patient. The sling label is located on the outside.
8. Make sure that the centre of the sling is aligned with the patient's spine, starting at the coccyx or against the seat. **Continue to step 16.**
9. **Patient not able to sit up:** Fold the sling lengthwise along the centre of the sling, outside facing outside. The sling label is located on the outside.

Continue with the steps on the next page.

Fig. 9

Fig. 10


Fig. 11


Fig. 12

10. Log roll the patient. (Use local technique)
(**See Fig. 13**) If log rolling is not possible, use an Arjo sliding sheet/tube for placing the sling. See respective *Sliding sheet/tube IFU*.
 11. Place the folded sling over the patient's side. Make sure that the centre of the sling is aligned with the patient's spine, starting at the coccyx. Place the leg flaps towards the legs.
 12. Fold down the top of the sling and tuck it under the patient's body.
 13. Log roll the patient to the other side. (Use local technique)
 14. Pull out the remaining part of the sling from under the patient's body.
 15. Return the patient to his/her back.
 16. Make sure that the: (**See Fig. 14**)
 - sling is centred and flat without creases,
 - sling's head support covers the neck/head area,
 - sling pieces are not twisted underneath the patient
 - all buckles are safely closed, if any
 - patient's arms are inside the sling.
- WARNING**
To avoid injury, make sure the patient's arms are placed inside the sling.
17. Position the lift sideways with the lift legs wide open. One caregiver needs to be close to the patient's head and the other caregiver needs to lift up the patient's legs.
 18. Position one lift leg next to the patient's head and the other lift leg under the patient's legs. Make sure the spreader bar is positioned over the patient's shoulders. (**See Fig. 15**)
- WARNING**
To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.
19. Lower the spreader bar with the Dynamic Positioning System (DPS) in a reclined position.
 20. Place the leg flaps underneath the patient's legs. Make sure the leg flaps are not twisted.
 21. Attach the shoulder clips. (**See Fig. 4**)
 22. Move the lift/spreader bar closer to the patient's legs.
 23. Apply the brakes on the lift.
 24. Attach the leg clips.
 25. Slightly lift the patient to create tension in the sling. (**See Fig. 16**)
- WARNING**
To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.
26. Make sure that all clips are attached securely and the patient lays comfortably in the sling. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clips.
 27. Make sure that:
 - the spreader bar is in a reclined position,
 - one caregiver pays attention to the patient's head and.
 - that the other caregiver holds up the patient's feet to prevent them from hitting the lift leg. (**See Fig. 16**)
 28. Remove the pillow.
 29. Release the brakes on the lift
 30. Bring the lift legs together.
 31. Transfer/transport the patient according to the lift *IFU*.

Fig. 13

Fig. 14

Fig. 15



Fig. 16

Removing the Sling

In Bed (12 steps)

1. Position the patient over the bed. **(See Fig. 17)**
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. If the patient will be lowered on a flat bed, make sure the spreader bar is in a reclined position.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

4. Lower the patient onto the bed. Make sure that the weight of the patient is taken up by the receiving surface before removing the clips.
5. Detach the sling. **(See Fig. 5)**
6. Move the lift away from the patient.
7. Pull out the leg flaps from under the patient's legs. **(See Fig. 14)**

WARNING

To avoid injury during log roll, make sure there is another caregiver present and/or use the side rails on the opposite side of the bed.

8. Log roll the patient. (Use local technique) If log rolling is not possible, use an Arjo sliding sheet/tube for removing the sling. See respective *Sliding sheet/tube IFU*.
9. Tuck the sling under the patient's body. **(See Fig. 18)**
10. Log roll the patient to the other side and remove the sling.
11. Return the patient to his/her back.
12. Make sure the patient is in a comfortable and secure position in the bed.

In Chair/Wheelchair (13 steps)

1. Apply the brakes on the wheelchair.
2. Position the patient over the chair/wheelchair. Open the lift legs if needed. **(See Fig. 19)**
3. Adjust the spreader bar to a seated position.
4. Lower the patient into the chair/wheelchair.
5. Make sure that the patient's lower back is positioned all the way in against chair/wheelchair to prevent the patient from sliding out.
6. Position the patient by using the DPS.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

7. Detach the sling. Make sure that the weight of the patient is taken up by the receiving surface before removing the clips. **(See Fig. 5)**
8. Release the brakes on the lift.
9. Move the lift away from the patient
10. Pull out the leg flaps from under the patient's legs and place them backwards along the patient's side. **(See Fig. 20).**
11. Stand in front of the patient and lean the patient forward.
12. Pull out the sling.
13. Lean the patient back up to a comfortable and secure position.

Fig. 17


Fig. 18


Fig. 19



Fig. 20

Cleaning and Disinfection

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

WARNING

To avoid material damage and injury, clean and disinfect according to this IFU.

- No other chemicals are allowed.
- Never clean with chlorine.
- Chlorine will deteriorate the surface of the material.

WARNING

To avoid injury, always remove the sling before disinfecting the lift.

All *Paediatric Slings* should be cleaned when they are soiled or stained and between patients.

Cleaning Instruction (5 steps)

1. Detach the sling from the lift.
2. Before washing close the buckles.
3. Machine wash the sling at 70°C (158°F). In order to obtain sufficient level of disinfection, please wash according to local hygiene legislation.
4. If allowed according to the product label, tumble dry in low temperature, maximum 60°C (140°F).
5. **Do NOT**
 - wash with rough surfaces or sharp objects
 - use any mechanical pressure, pressing or roiling during the washing and drying process
 - use bleach
 - use gas sterilization
 - use autoclave
 - dry clean
 - Steam
 - ironing

Cleaning Chemicals

For all sling material use common commercial detergent without optical brightener.

No other chemicals are allowed e.g. chlorine, softener, iodine based disinfectants, bromine and ozone.

Disinfection

The only allowed disinfection is disinfection by washing.

Care and Preventive Maintenance

WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

Before and After Every Use

Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching
- Tears
- Holes
- Discolouration and stains from bleaching
- Sling soiled or stained
- Unreadable or damaged label

When Soiled or Stained and Between Patients

Clean/Disinfect

The caregiver should make sure that the sling is cleaned according to “*Cleaning and Disinfection*” on page 12, when it is soiled or stained and between patients.

Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or to excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 month), according to ISO 10535.

Troubleshooting

Problem	Action
The patient is not positioned correctly in the sling (e.g. one strap is shorter than the other).	<ul style="list-style-type: none"> • Make sure the clip straps are not twisted. • Make sure the patient is seated on the inside of the sling. The outside has a label for reference. • Make sure the patient is centred in the sling.
The patient feels discomfort in the leg area when seated in the sling.	Make sure that there are no creases on the leg flaps of the sling.
The sling is difficult to apply (e.g. during log roll).	Use an Arjo sliding sheet/tube according to instructions in respective <i>IFU</i> .
There are difficulties applying the shoulder clips.	<ul style="list-style-type: none"> • Place the patient in a more seated position in the bed or on the floor. Place a pillow behind the patient's neck/back, outside of the sling. • Raise the back rest of the bed, if available. • Make sure the clip strap is not caught between the lug on the lift and the clip attachment. • Adjust the spreader bar to a more reclined position. Pay attention to the patient's head.
There are difficulties applying the leg clips.	<ul style="list-style-type: none"> • Bend the patient's legs or place a pillow under the patient's knees. If the sling is applied on the floor, place the patient's legs on the caregiver's knee. • Make sure the clip strap is not caught between the lug on the lift and the clip attachment. • Adjust the spreader bar to a more seated position.
It is difficult to use the sling in combination with the lift and spreader bar.	Use the sling only in accordance with the combinations described in section " <i>Allowed Combinations</i> " on page 7.
The position of the patient in the sling should be more reclined or seated.	Use the spreader bar to adjust the patient's position.
The patient wishes to have extra support for the neck.	Place a pillow behind the patient's neck, inside the sling.
The head is not supported by the sling.	Make sure the head section of the sling is buckled together with the shoulder section of the sling.

EN

Technical Specifications

General

Safe working load (SWL) = Maximum patient weight	See "Allowed Combinations" on page 7
Service life – Recommended period of use	2 years* * Refer to "Expected Service Life" on page 3.
Shelf life – Maximum period of storing new unpacked product	5 years
Model and Type	See "Sling Selection" on page 6

EN

Operating, Transport and Storage Environment

Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (+13°F to +140°F) Transport
Humidity	Max 15–70% at +20°C (68°F) Operating and Storage Max 10–95% at +20°C (68°F) Transport

End of Life Disposal

Package	The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.

Parts and Accessories

Measuring tape	MTA1000
----------------	---------

Label on the Sling

Care and Washing Symbols

	Machine wash 70°C (158°F)
	No bleaching allowed
	No tumble drying allowed
	Tumble drying
	No Ironing allowed
	No dry cleaning allowed

Certificates/Markings

	CE marking indicating conformity with European Community harmonised legislation
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

Symbol for Clip

	Use a clip spreader bar.
---	--------------------------

Fibre Content

PES	Polyester
PU	Polyurethane

Article Number

REF XXXXXX-X	Article number with -X refers to the sling size.
REF XXXXXXX	If the article number is stated without the size letter at the end, it means the sling is a One Size.

Misc. Symbols

	Safe Working Load (SWL)
	Patient's name symbol
	Record symbol
	Read the IFU before use
	Article number
	Manufacturing date and year
	Manufacturer name and address

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HONG KONG

Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

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Address page - REV 24: 04/2019

www.arjo.com

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges.

With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

04.SP.00-INT1

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Arjo Inc.