INSTRUCTIONS FOR USE

AtmosAir

Mattress Replacement System

Sistema de Troca de Colchões


ATMOSAIR WITH SAT 4000 SERIES

ATMOSAIR WITH SAT 9000 SERIES

ATMOSAIR WITH SAT T-SERIES

ATMOSAIR WITH SAT V-SERIES

ATMOSAIR WITH SAT AP0D 25

EN · PTBR

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WARNING

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.

Mandatory to read the *Instructions for Use*.

® and ™ are trademarks belonging to the Arjo group of companies

© Arjo 2019. SAT™ (Self Adjusting Technology™) is a trademark of MPTI, Ltd. Total Care and VersaCare are trademarks of Hill-Rom and / or its affiliates.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

EN

ARJO HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJO PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY ARJO SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL ARJO BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJO TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties except as set forth in the written limited warranty included in this publication or with this product. Information in this publication may be subject to change at any time. Contact Arjo for updates.

Important Information For Users

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.
- If applicable, ensure the electrical installation of the room complies with the appropriate national / local electrical wiring standards.

Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

Notice

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the product information label for specific voltage.

Intentionally left blank

Table Of Contents

Introduction	2
Indications	2
Contraindications	2
Risks and Precautions	2
Safety Information	3
Serious Incident	4
Preparation for Use	5
Mattress Installation.....	5
<i>AtmosAir</i> with SAT 4000, 9000 (all models), T-Series*, V-Series** (all models), and APOD 25 (all models).....	5
<i>AtmosAir</i> with SAT V-Series	6
<i>AtmosAir</i> Pump Installation	7
LED Indicators	8
Patient Placement and Nursing Care	9
Alternating Pressure / Rotation Adjustment.....	9
CPR.....	10
Transport Mode	10
Skin Care.....	10
Incontinence / Drainage	10
General Operation.....	10
Care and Cleaning.....	11
Fire Barrier	11
Cover Cleaning Options	11
Cleaning <i>AtmosAir</i> Pump and Tubes (A and AR models)	12
Preventive Maintenance Schedule	13
Daily Cleaning	13
Weekly Cleaning.....	13
Inspection / System Check-Out	13
Parts Diagram - <i>AtmosAir</i> with SAT 4000.....	15
Parts Diagram - <i>AtmosAir</i> with SAT 9000 and T-Series	16
Parts Diagram - <i>AtmosAir</i> with SAT V-Series	17
Parts Diagram - <i>AtmosAir</i> with SAT APod 25	18
Parts Diagram - <i>AtmosAir</i> Pump (A and AR models)	19
Replacement Parts	20
Specifications	24
Electromagnetic Compatibility	25
Symbols Used	28
Customer Contact Information.....	29

Introduction

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications, Safety Information and Risks and Precautions** sections prior to placing a patient on any AtmosAir® with SAT™ Mattress Replacement System (MRS).

Caregivers should review this information with the patient and the patient's family and / or legal guardian. Save this User Guide in an easily accessible location for quick reference.

This User Guide is comprised of the following AtmosAir MRS models / series:

- *AtmosAir* with SAT 4000 Series (non powered, indicates pump not available)
- *AtmosAir* with SAT 9000 Series (non powered)
- *AtmosAir* with SAT 9000A Series (powered with alternating pressure)
- *AtmosAir* with SAT 9000AR Series (powered with alternating pressure and rotation)
- *AtmosAir* with SAT T-Series (non powered)
- *AtmosAir* with SAT T-Series A (powered with alternating pressure)
- *AtmosAir* with SAT T-Series AR (powered with alternating pressure and rotation)
- *AtmosAir* with SAT V-Series (non powered)
- *AtmosAir* with SAT APOD 25 (non powered)
- *AtmosAir* with SAT APOD 25 A (powered with alternating pressure)
- *AtmosAir* with SAT APOD 25 AR (powered with alternating pressure and rotation)

The *AtmosAir* MRS comes in a wide variety of lengths and widths to support various customer requirements. Contact an Arjo representative for additional product information.

Indications

- *AtmosAir* MRS are intended to treat and prevent all categories of pressure injuries.

Contraindications

- unstable vertebral fracture
- cervical and skeletal traction

Risks and Precautions

Transfer – Standard precautions should be taken during patient transfer.

Rotation / Alternating Pressure (for A and AR models) – Prior to engaging rotation or alternating pressure, ensure that bed frame has side rails and that all side rails are fully engaged in full upright and locked position. Ensure pump hoses are properly connected to appropriate side of mattress for the desired function (see ***AtmosAir* Pump Installation**).

WARNING: Incorrect connection of pump hoses can increase risk of potential patient falls.

Side Rails and Restraints – WARNING: Use or non-use of restraints, including side rails, can be critical to patient safety. Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. **See related Safety Information.**

Patient Migration – Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. **Monitor patients frequently to guard against patient entrapment.**

Oxygen Use – DANGER: Risk of explosion if the *AtmosAir* Pump is used in the presence of flammable anesthetics. Use of this product's pump in an oxygen-enriched environment may produce potential of fire hazard. This equipment is not suitable for use in the presence

of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. Unplug and do not use pump when using oxygen-administering equipment other than the nasal mask or half-bed-length tent type.

Shock Hazard – Electrical shock hazard; do not remove pump case covers. Refer to qualified service personnel.

Safety Information

Patient Entrance / Exit – Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Turning – CAUTION: Prior to engaging turn feature, ensure that bed frame has side rails and that all side rails are in their full up and locked position.

Brakes – Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

Bed Height – To minimize risk of falls or injury, the bed should always be in the lowest practical position when the patient is unattended.

Bed Frame – Always use a standard healthcare bed frame with this mattress, with any safeguards or protocols that may be appropriate. Bed frame and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. It is recommended that bed and side rails (if used) comply with all applicable regulations and protocols.

CPR – Level the bed and disconnect pump (if present). Lower side rails and initiate CPR per facility protocols. Consider use of backboard if indicated. After CPR remove backboard, if used, reconnect pump, if present, raise siderails and reconfigure bed and accessories as in initial placement.

Head of Bed Elevation – Keep head of bed as low as possible to help prevent patient migration.

Side Rails / Patient Restraints – Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patients frequently to guard against patient entrapment.

CAUTION: When selecting a standard mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

I.V. and Drainage Tubes – I.V. and drainage tubes should always have slack for alternating pressure or rotation and other patient movements.

Skin Care – Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

Fluids – Avoid spilling fluids on pump controls. If spills do occur, clean fluid from pump wearing rubber gloves or while unit is unplugged to avoid any possibility of shock. Once fluid is removed, check operation of components in area of spill.

Fluids remaining on controls can cause corrosion, which may cause components to fail or operate erratically, possibly producing potential hazards for patient and staff.

Avoid Fire Hazards – To minimize risk of fire, connect the bed's power cord directly into a wall-mounted outlet. Do not use extension cords or multiple outlet strips. In the US, review and follow FDA's Safety Tips for Preventing Hospital Bed Fires.

No Smoking in Bed – Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

Power Cord – Position power cord to avoid a tripping hazard and / or damage to the cord. Ensure power cord is kept free from all pinch points and moving parts and is not trapped under casters. Improper handling of the power cord can cause damage to the cord, which may possibly produce risk of fire or electric shock.

General Protocols – Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

Fire Barrier Contact – Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. If the fire barrier becomes damaged or soiled it should be replaced.

End of Life Disposal –

- Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.
- Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.
- Pump units have electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Preparation for Use

For information concerning the bed frame refer to the manufacturers' User Guide.

1. Open shipping container(s).

Do not use sharp instruments to open boxes. Damage to mattress could result.

2. Remove *AtmosAir* with SAT MRS from plastic protective cover.

The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to accommodate; see Troubleshooting for more information. Wrinkles will not affect inflation or function, so mattress may be used immediately if needed.

3. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing mattress onto a new frame or for a new patient, check mattress surface for staining and soiling; clean and / or disinfect as required (see **Care and Cleaning**).
5. Level bed and lock brakes.
6. Remove existing mattress from bed frame.

Mattress Installation

***AtmosAir* with SAT 4000, 9000 (all models), T-Series*, V-Series** (all models), and APOD 25 (all models)**

1. Position mattress on bed frame with logo facing up and product information tags at foot end of bed.

Magnets are integrated on the *AtmosAir* T-Series bottom cover to help keep the mattress in place during foot section extension / retraction.

2. Ensure there are no gaps between mattress and bed frame or side rails.

Always use a standard healthcare bed frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.

* T-Series is compatible for use with the Hill-Rom Total Care® bed frame.

** V-Series is compatible for use with the Hill-Rom VersaCare® frame.

AtmosAir with SAT V-Series

V-Series mattress on VersaCare® Frame

1. Position mattress on bed frame with logo facing up and product information tags at foot end of bed.
2. Slide locking bar located in pocket on underside of head-end of mattress to one side, as shown at right (Figure 1).


Figure 1

3. Slip end of locking bar into slot on head-end of bed frame, as shown below (Figures 2 and 3).
4. Slip other end of locking bar into other slot on head-end of bed frame, as shown below (Figure 4).


Figure 2


Figure 3

5. Repeat steps 2 through 4 to secure foot-end of mattress.

Foot-end locking bar may need to be flexed slightly to slip into foot-end slots.

Magnets are integrated on the AtmosAir V-Series bottom cover to help keep the mattress in place during foot section extension / retraction.


Figure 4


Figure 5

When removing the AtmosAir with SAT V-Series, be sure to disengage locking bar and magnets before removing from bed frame.

AtmosAir Pump Installation

Pump installation steps for A and AR models:

1. Place pump on a solid stationary surface or suspend on end of bed frame with built-in hanger.
2. Attach pump hoses to mattress, as shown at right (Figure 6):

A Models (alternating pressure):

- Attach **blue** hoses to the mattress connectors on the **patient left side** (marked with **BLUE** tag) of the mattress for **alternating pressure**.

AR Models (alternating pressure / rotation):

- Attach **blue** hoses to the mattress connectors on the **patient left side** (marked with **BLUE** tag) of the mattress for **alternating pressure**.
- Attach **red** hoses to the mattress connectors on the **patient right side** (marked with **RED** tag) of the mattress for **rotation**.


Figure 6

3. Ensure air hoses are not kinked and will not be pinched by any articulated bed mechanisms.
4. Attach hose connectors to the pump.
5. Plug pump unit into a properly-grounded wall outlet.
6. Verify power to this outlet is not controlled by a wall switch.
7. Switch pump on.

LED Indicators

1. Green LED is lit when the pressure is at the set level.
2. Yellow LED is lit when the pressure is below the set level. Check connections.


Figure 7

Patient Placement and Nursing Care

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications**, **Safety Information** and **Risks and Precautions** sections prior to placing a patient on any *AtmosAir* MRS.

1. Transfer patient following all applicable safety rules and institution protocols.
2. Center patient side-to-side and head-to-foot on *AtmosAir* MRS surface.
3. Ensure all sections of the mattress fully support the patient.

Mattress handles are for mattress transport only.

Alternating Pressure / Rotation Adjustment

The *AtmosAir* with SAT 9000 and T-Series models provide alternating pressure and / or rotation. The level of alternating pressure or rotation can be adjusted for patient comfort.

Prior to engaging alternating pressure or rotation, ensure that bed frame has side rails and that all side rails are fully engaged in full upright and locked position.

Ensure pump hoses are properly connected to appropriate side of mattress for the desired function (see *AtmosAir Pump Installation*).

- Rotate the control knob on the pump clockwise to increase the intensity of alternating pressure relief or the level of rotation.
- Rotate the control knob on the pump counter-clockwise to reduce the intensity of alternating pressure relief or the level of rotation.

The *AtmosAir* MRS DOES NOT provide a Nurse Hold function.


Figure 8

CPR

1. Level bed.
2. Disconnect hoses from pump to level mattress (for A and AR models), as shown at right (Figure 8).
3. Lower or remove side rail on caregiver's side if necessary.
4. Begin CPR as per facility protocol. Consider use of backboard if indicated.
5. After CPR is performed:
 - Remove backboard if used
 - Reconnect pump hoses to mattress (for A and AR models).
 - Raise or install side rail as necessary.
 - Reconfigure bed and accessories as in initial placement.


Figure 9

Transport Mode

When disconnected from the pump, the *AtmosAir* 9000A MRS will operate as a reactive surface with the full operation of the SAT system. For transport, there are two options:

1. Disconnect the tube assembly from the mattress. Valves will close off as soon as the tube assembly is disconnected.
2. Disconnect the tube assembly from the pump and securely place the transport cap on the tube assembly to maintain a seal.

Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under patient are not wrinkled.

Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace bed linens as required (see **Care and Cleaning.**)

General Operation

Avoid contact of sharp instruments with the *AtmosAir* MRS. Punctures, cuts and tears may prevent proper inflation and air pressure maintenance.

Figure 10

Care and Cleaning

The following processes is recommended, but should be adapted to comply with local institution protocols. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The *AtmosAir* MRS should be routinely decontaminated between patients and at regular intervals while in use.

If a powered *AtmosAir* MRS, ensure power cord is unplugged from wall outlet before performing care and cleaning.

Do not use Phenol-based solutions or abrasive compounds or pads on the cover during decontamination process as these will damage the surface coating. Do not boil or autoclave the cover. Avoid immersing electrical parts in water during the cleaning process. Do not spray cleaning solution directly onto the pump.

Covers should be wiped down immediately after it has been exposed to fluids or spills.

Fire Barrier

Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. Fire barrier should be replaced if soiled or damaged.

Cover Cleaning Options

Sewn *AtmosAir* MRS (Top cover can not be removed from base)

1. If a powered MRS, unplug pump before cleaning.
2. Remove or push bed linens to center of mattress.

The sewn *AtmosAir* MRS cover is cleaned by wipe down methods only. Do not launder as damage to the base may occur.

3. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
4. Rinse with a wet towel to remove chemicals.
5. After wipe down, dry surface with towel.
6. Ensure bed linens are refitted and not wrinkled under patient.
7. Clean pump and tubes (if needed) by wiping with a damp cloth.

RF Welded *AtmosAir* MRS (Detachable top cover)

1. If a powered MRS, unplug pump before cleaning.
2. Push bed linens to center of mattress to wipe down, remove to launder top cover or wipedown as described below.

The RF Welded *AtmosAir* MRS has a detachable top cover that can be laundered. However the base is cleaned by wipe down methods only. Do not launder base as damage may occur.

3. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.

4. Rinse with a wet towel to remove chemicals.
5. After wipe down, dry surface with towel.
6. Ensure bed linens are refitted and not wrinkled under patient.
7. Clean pump and tubes (if needed) by wiping with a damp cloth.

Laundering Detachable Top Cover

1. Unzip the top cover from the base for laundering.
Do not launder base as damage may occur.
2. Recommended wash temperature for top cover is 60° C (140° F) for 15 minutes.
3. Maximum wash temperature is 95° C (203° F) for 15 minutes.
4. Tumble dry at 60° C (140° F) or air dry.
5. Maximum drying temperature 80° C (176° F).

Cleaning *AtmosAir* Pump and Tubes (A and AR models)

1. Ensure pump is unplugged.
2. Wipe surface of pump and tubes with coarse cloth using an approved disinfectant germicide mixed to manufacturer's instructions. **Do not flood any part of the pump with cleaning solution.**

Avoid spilling fluids on pump controls. If spills do occur, clean fluid from pump wearing rubber gloves or while unit is unplugged to avoid any possibility of shock. Once fluid is removed, check operation of components in area of spill. Fluids remaining on controls can cause corrosion, which may cause components to fail or operate erratically, possibly producing potential hazards for patient and staff.

3. Wipe off excess solution.
4. Disinfect pump and tubes with a chlorine solution (mixed according to the instructions in the **Care and Cleaning** section). Using a clean cloth, wring out excess solution until cloth is damp. Wipe pump surface and tubes with damp cloth.
5. Allow to air dry.

Preventive Maintenance Schedule

Preventive maintenance for the *AtmosAir* MRS consists of regular cleaning (see **Care and Cleaning**) and an overall system check-out to be performed at the intervals described below.

All components must be cleaned, disinfected and inspected after each patient's use and before use by a new patient. Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

Daily Cleaning

The cover should be wiped daily with a mild soap and water solution. After wipe down, dry surface with a towel.

Weekly Cleaning

The pump and hoses of A and AR models should be cleaned weekly.

Inspection / System Check-Out

Check each of the following before placing the *AtmosAir* MRS with a new patient:

1. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
2. Ensure mattress is free of stains and is not overly faded.

For A and AR models:

1. Ensure air inlet hoses and connectors on mattress and pump are clean and undamaged.
2. Ensure pump and power cord are clean and undamaged.
3. Ensure pump hanger brackets are secure and operate correctly.
4. Ensure power switch and comfort control knob both operate correctly.
5. Attach pump to the red **rotation** hoses and power on to ensure that the mattress surface tilts and there are no air leaks.
6. Attach pump to the blue **alternating pressure** hoses and power on to ensure there are no air leaks.

Troubleshooting

Do not attempt troubleshooting outside this guide or where the solution recommends to contact Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

SYMPTOM	POSSIBLE CAUSE	SOLUTION
Green lamp on pump does not light when power switch is set to the ON position. Mattress inflates only to one side when in rotation mode.	Power cord may be unplugged.	Plug in power cord.
	Wall outlet may be controlled by a wall switch.	Plug power cord into different wall outlet.
	Circuit breaker for wall switch may be tripped or blown.	Reset circuit breaker or replace fuse.
	Power cord or power switch may be damaged.	Contact Arjo for assistance.
	Tubing not connected properly.	Ensure red hoses are connected to hoses marked with red tag for rotation.
	Tubing kinked.	Check tubing inside mattress for kinks.
	Tubing disconnected.	Check tubing inside mattress for possible disconnect.
Mattress too firm upon arrival.	Hole in rotation bladder.	Check deflated rotation bladder for holes.
Mattress cover too wrinkled upon removal from shipping container.	Difference in altitude not sufficient to open valves.	Apply weight to mattress to open valves.
	Internal components have not accommodated to environment. This does not affect inflation or function.	Let mattress accommodate for 24 hours. If problem continues, contact Arjo for assistance.
Mattress is not firm.	Tubing not connected properly.	Check tubing inside mattress for loose connectors.
	Tubing kinked.	Check tubing inside mattress for possible kinks.
	Tubing disconnected.	Check tubing inside mattress for possible disconnect.
	Holes in or damaged to SAT system.	Check SAT system for holes or damage, or contact Arjo for assistance.


Parts Diagram - *AtmosAir* with SAT 4000

All hoses were removed to improve diagram readability.

The fire barrier is a sleeve that fits over the foam and SAT assemblies.

The **SAT Cover**, **SAT Cells**, **SAT Interliner** and rotation bladders (for A and AR models) are all part of the integrated **SAT System** and cannot be ordered separately. See Replacement Parts for a complete list of **SAT Systems**.


Parts Diagram - *AtmosAir* with SAT 9000 and T-Series

All hoses were removed to improve diagram readability.

The fire barrier is a sleeve that fits over the foam and SAT assemblies.

The **SAT Cover**, **SAT Cells**, **SAT Interliner** and rotation bladders (for A and AR models) are all part of the integrated **SAT System** and cannot be ordered separately. See Replacement Parts for a complete list of SAT Systems.


Parts Diagram - *AtmosAir* with SAT V-Series

All hoses were removed to improve diagram readability.

The fire barrier is a sleeve that fits over the foam and SAT assemblies.

The SAT Cover, SAT Cells, SAT Interliner are all part of the integrated SAT System and cannot be ordered separately. See Replacement Parts for a complete list of SAT Systems.


Parts Diagram - *AtmosAir* with *SAT APod 25*

All hoses were removed to improve diagram readability.

The Fire Barrier is a sleeve that fits over the foam and SAT assemblies.


The *SAT Cover*, *SAT Cells*, *SAT Interliner*, and *SAT Air Pods* are all part of the integrated *SAT System* and cannot be ordered separately. See *Replacement Parts* for a complete list of *SAT Systems*.

Parts Diagram - *AtmosAir* Pump (A and AR models)

See Replacement Parts for a complete list of parts that may be ordered for the *AtmosAir* Pump, as well as entire pump systems.

Replacement Parts

Do not attempt troubleshooting, maintenance or parts replacement outside this manual or where the solution recommends contacting Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Replaceable *AtmosAir* components are listed below. For more information such as pricing or additional spare parts that are not on this list, please contact your local Arjo representative.

ATMOSAIR MATTRESS COMPONENTS

AtmosAir Foam Shells

PART TYPE.....	CATALOG NUMBER
4000 / 9000 Foam Shell	KA4/9FS(SIZE)
T-Series Foam Shell	KA9FSTC
V-Series Foam Shell.....	KVFS3584
APod 25 Foam Shell.....	25PFS(SIZE)

AtmosAir with SAT System

PART TYPE.....	CATALOG NUMBER
4000 SAT System.....	KA4SAT
9000 SAT System	KA9SAT
9000A SAT System	KA9ASAT
9000AR SAT System	KA9ARSAT
V-Series SAT System	KVSAT
T-Series SAT System	KA9-T-SAT
T-Series SAT System A	KA9A-T-SAT
T-Series SAT System AR.....	KA9AR-T-SAT
APod 25 SAT (Fits mattresses 32-38 in wide, 75-84 in length).....	25PSAT
APod 25A SAT (Fits mattresses 32-38 in wide, 75-84 in length).....	25PASAT
APod 25AR SAT (Fits mattresses 32-38 in wide, 75-84 in length).....	25PARSAT

AtmosAir Fire Barrier Lining

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> Fire Barrier Lining (Standard on all models)	KA4/9-FBL

RF Welded *AtmosAir* Fire Barrier Lining (Standard on all *AtmosAir* RF Welded mattresses

[all models]) AtmosAir-FBL

ATMOSAIR MATTRESS COVERS

AtmosAir Model 4000 Mattress

PART TYPE	CATALOG NUMBER
Complete Cover, Soflux™ Top, Non-Skid Bottom	A4CSOVB(WWWLL)NXS
Complete Cover, Soflux Top, Vinyl Gray Bottom	A4CSOVG(WWWLL)NXS
Complete Cover, Reliant IS ² ™ Top, Non-Skid Bottom	A4CRLNB(WWWLL)TXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom	A4CBFNB(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom, No Handles	A4CBFNH(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom	A4CBFVG(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom, No Handles	A4CBFVH(WWWLL)PXS
Complete Cover, Reliant IS ² Top, Vinyl Gray Bottom	A4CRLVG(WWWLL)TXS
Complete Cover, Reliant IS ² Top, Vinyl Gray Bottom, No Handles	A4CRLVH(WWWLL)TXS
Complete Cover, Reliant IS ² , RF Welded Top, Non-Skid PU Bottom, No Handles	A4CRFPN(WWWLL)TXS
Complete Cover, Reliant IS ² , RF Welded Top, Non-Skid PU Bottom, With Handles	A4CRFPU(WWWLL)TXS
Top Cover, Reliant IS ² , RF Welded	A4CRFFX(WWWLL)TXS
Base Cover, PU Non-skid Bottom, No Handles	A4CXXPN(WWWLL)XXS
Base Cover, PU Non-skid Bottom, With Handles	A4CXXPU(WWWLL)XXS

*WWW = width in 0.5 in increments, LL = Length in 1.0 in increments

AtmosAir Model 9000 Mattress

PART TYPE	CATALOG NUMBER
Complete Cover, Soflux Top, Non-Skid Bottom	A9CSOVB(WWWLL)NXS
Complete Cover, Soflux Top, Non Skid Bottom, No Handles	A9CSONH(WWWLL)NXS
Complete Cover, Soflux Top, Vinyl Gray Bottom	A9CSOVG(WWWLL)NXS
Complete Cover, Soflux Top, Vinyl Gray Bottom, No Handles	A9CSOVH(WWWLL)NXS
Complete Cover, Reliant IS ² Top, Non-Skid Bottom	A9CRLNB(WWWLL)TXS
Complete Cover, Reliant IS ² Top, Non Skid Bottom, No Handles	A9CRLNH(WWWLL)TXS
Complete Cover, Reliant IS ² Top, Vinyl Gray Bottom	A9CRLVG(WWWLL)TXS
Complete Cover, Reliant IS ² Top, Vinyl Gray Bottom, No Handles	A9CRLVH(WWWLL)TXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom	A9CBFNB(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom, No Handles	A9CBFNH(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom	A9CBFVG(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom, No Handles	A9CBFVH(WWWLL)PXS
Complete Cover, Reliant IS ² , RF Welded Top, Non-Skid PU Bottom, No Handles	A9CRFPN(WWWLL)TXS
Complete Cover, Reliant IS ² , RF Welded Top, Non-Skid PU Bottom, With Handles	A9CRFPU(WWWLL)TXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, No Handles	A9CPRPN(WWWLL)NXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, with Handles	A9CPRPU(WWWLL)NXS
Top Cover, Reliant IS ² , RF Welded	A9CRFFX(WWWLL)TXS
Top Cover, Premium, RF Welded	A9CPRXX(WWWLL)NXS
Base Cover, PU Non-skid Bottom, No Handles	A9CXXPN(WWWLL)XXS
Base Cover, PU Non-skid Bottom, With Handles	A9CXXPU(WWWLL)XXS

AtmosAir Model 9000A Mattress with Alternating Pressure

PART TYPE	CATALOG NUMBER
Complete Cover, Soflux Top, Non-Skid Bottom	AXCSOVB(WWWLL)NXS
Complete Cover, Soflux Top, Vinyl Gray Bottom	AXCSOVG(WWWLL)NXS
Complete Cover, Reliant IS ² Top, Non-Skid Bottom	AXCRLNB(WWWLL)TXS
Complete Cover, Reliant IS ² Top, Vinyl Gray Bottom	AXCRLVG(WWWLL)TXS
Complete Cover, Soflux Top, Non-Skid Bottom, No Handles	AXCSONH(WWWLL)NXS
Complete Cover, Reliant IS ² Top, Non-Skid Bottom	AXCRLNH(WWWLL)TXS
Complete Cover, Soflux Top, Vinyl Gray Bottom, No Handles	AXCSOVH(WWWLL)NXS
Complete Cover, Reliant IS ² Top, Vinyl Gray Bottom	AXCRLVH(WWWLL)TXS
Complete Cover, Reliant IS ² , RF Welded Top, Non-Skid PU Bottom, No Handles	AXCRFPN(WWWLL)TXS
Complete Cover, Reliant IS ² , RF Welded Top, Non-Skid PU Bottom, With Handles	AXCRFPU(WWWLL)TXS

Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, No Handles.....	AXCPRPN(WWWLL)NXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, with Handles.....	AXCPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS²</i> , RF Welded.....	AXCRFX(XWWLL)TXS
Top Cover, Premium, RF Welded.....	AXCPRXX(WWWLL)NXS
Base Cover, PU Non-Skid Bottom, No Handles.....	AXCXXPN(WWWLL)XXS
Base Cover, PU Non-Skid Bottom, With Handles.....	AXCXXPU(WWWLL)XXS

AtmosAir Model 9000AR Mattress with Alternating Pressure and Rotation

PART TYPE.....	CATALOG NUMBER
Complete Cover, <i>Sofflux</i> Top, Non-Skid Bottom.....	ARCSOVB(WWWLL)NXS
Complete Cover, <i>Sofflux</i> Top, Vinyl Gray Bottom.....	ARCSOVG(WWWLL)NXS
Complete Cover, <i>Reliant IS²</i> Top, Non-Skid Bottom.....	ARCRLNB(WWWLL)TXS
Complete Cover, <i>Reliant IS²</i> Top, Vinyl Gray Bottom.....	ARCRLVG(WWWLL)TXS
Complete Cover, <i>Sofflux</i> Top, Non-Skid Bottom, No Handles.....	ARCSOVB(WWWLL)NXS
Complete Cover, <i>Reliant IS²</i> Top, Non-Skid Bottom.....	ARCRLNH(WWWLL)TXS
Complete Cover, <i>Sofflux</i> Top, Vinyl Gray Bottom, No Handles.....	ARCSOVG(WWWLL)NXS
Complete Cover, <i>Reliant IS²</i> Top, Vinyl Gray Bottom.....	ARCRLVH(WWWLL)TXS
Complete Cover, <i>Reliant IS²</i> , RF Welded Top, Non-Skid PU Bottom, No Handles.....	ARCRFPN(WWWLL)TXS
Complete Cover, <i>Reliant IS²</i> , RF Welded Top, Non-Skid PU Bottom, With Handles.....	ARCRFPU(WWWLL)TXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, No Handles.....	ARCPRPN(WWWLL)NXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, with Handles.....	ARCPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS²</i> , RF Welded.....	ARCRFX(XWWLL)TXS
Top Cover, Premium, RF Welded.....	ARCPRXX(WWWLL)NXS
Base Cover, PU Non-Skid Bottom, No Handles.....	ARCXXPN(WWWLL)XXS
Base Cover, PU Non-Skid Bottom, With Handles.....	ARCXXPU(WWWLL)XXS

AtmosAir Model APod25 Mattress

PART TYPE.....	CATALOG NUMBER
Complete Cover, <i>Sofflux</i> Top, Non-Skid Bottom.....	ADCSONB(WWWLL)NXS
Complete Cover, <i>Sofflux</i> Top, Vinyl Gray Bottom.....	ADCSOVG(WWWLL)NXS
Complete Cover, <i>Reliant IS²</i> Top, Non-Skid Bottom.....	ADCRLNB(WWWLL)TXS
Complete Cover, <i>Reliant IS²</i> Top, Vinyl Gray Bottom.....	ADCRLVG(WWWLL)TXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, No Handles.....	ADCPRPN(WWWLL)NXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, With Handles.....	ADCPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS²</i> , RF Welded.....	ADCRFXX(WWWLL)TXS
Top Cover, Premium, RF Welded.....	ADCPRXX(WWWLL)NXS
Base Cover, PU Non-Skid Bottom, No Handles.....	ADCXXPN(WWWLL)XXS
Base Cover, PU Non-Skid Bottom, With Handles.....	ADCXXPU(WWWLL)XXS

AtmosAir Model APod 25 A Mattress with Alternating Pressure

PART TYPE.....	CATALOG NUMBER
Complete Cover, <i>Sofflux</i> Top, Non-Skid Bottom.....	DXCSONB(WWWLL)NXS
Complete Cover, <i>Sofflux</i> Top, Vinyl Gray Bottom.....	DXCSOVG(WWWLL)NXS
Complete Cover, <i>Reliant IS²</i> Top, Non-Skid Bottom.....	DXCRLNB(WWWLL)TXS
Complete Cover, <i>Reliant IS²</i> Top, Vinyl Gray Bottom.....	DXCRLVG(WWWLL)TXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, No Handles.....	DXCPRPN(WWWLL)NXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom, With Handles.....	DXCPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS²</i> , RF Welded.....	DXCRFX(XWWLL)TXS
Top Cover, Premium, RF Welded.....	DXCPRXX(WWWLL)NXS
Base Cover, PU Non-Skid Bottom, No Handles.....	DXCXXPN(WWWLL)XXS
Base Cover, PU Non-Skid Bottom, With Handles.....	DXCXXPU(WWWLL)XXS

AtmosAir Model APod 25 AR Mattress with Alternating Pressure and Rotation

PART TYPE.....	CATALOG NUMBER
Complete Cover, <i>Sofflux</i> Top, Non-Skid Bottom	DRCSONB(WWWLL)NXS
Complete Cover, <i>Sofflux</i> Top, Vinyl Gray Bottom	DRCISOVG(WWWLL)NXS
Complete Cover, <i>Reliant IS²</i> Top, Non-Skid Bottom	DRCRLNB(WWWLL)TXS
Complete Cover, <i>Reliant IS²</i> Top, Vinyl Gray Bottom	DRCRLVG(WWWLL)TXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Bottom, No Handles	DRCPRPN(WWWLL)NXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Bottom, With Handles	DRCPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS²</i> , RF Welded	DRCRFXX(WWWLL)TXS
Top Cover, Premium, RF Welded	DRCPRXX(WWWLL)NXS
Base Cover, PU Non-Skid Bottom, No Handles	DRCXXPN(WWWLL)XXS
Base Cover, PU Non-Skid Bottom, With Handles	DRCXXPU(WWWLL)XXS

AtmosAir - T-Series

PART TYPE.....	CATALOG NUMBER
T-Series Mattress - Complete Cover, <i>Sofflux</i> top and Non-Skid base	ATCONB34079NXS
T-Series Mattress - Complete Cover, <i>Reliant IS²</i> top and Non-Skid base	ATCRLNB34079TXS
T-Series A Mattress with Alternating Pressure - Complete Cover, <i>Sofflux</i> top and Non-Skid base	TXCSONB34079NXS
T-Series A Mattress with Alternating Pressure - Complete Cover, <i>Reliant IS²</i> top and Non-Skid base	TXCRLNB34079TXS

AtmosAir - V-Series

PART TYPE.....	CATALOG NUMBER
V-Series Mattress - Complete Cover, <i>Sofflux</i> top and Vinyl Gray base	AVCSOVG35084NXS
V-Series Mattress - Complete Cover, <i>Reliant IS²</i> top and Vinyl Gray base	AVCRLVG35084TXS

ATMOSAIR PUMPS**AtmosAir Standard Pumps (A and AR Models only)**

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> EZ Care Pump only, Current Version	KA9PUMP
<i>AtmosAir</i> Harness / Tubing / Connectors (no pump)	KA9PUMPASMBLY
<i>AtmosAir</i> EZ Care Pump Filter	KA9EZPMPFILTER
<i>AtmosAir</i> EZ Care Pump Filter Cap	KA9EZPMPFILCAP
<i>AtmosAir</i> Elite Pump Filter, compatible with older pump version	KA9PUMPFILTER
<i>AtmosAir</i> Elite Pump Filter Cap, compatible with older pump version	KA9PUMPFILCAP

ATMOSAIR LITERATURE

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> MRS User Guide	407390-AH

Specifications

Specifications subject to change without notice.

Mattresses

Maximum Recommended Patient Weight*500 lb (227 kg)

Weight (based on 35 in x 80 in x 7 in MRS, other sizes will vary)

AtmosAir with SAT 4000 MRS30 lb (14 kg)

AtmosAir with SAT 9000 MRS33 lb (15 kg)

AtmosAir with SAT 9000A and 9000AR MRS40 lb (18 kg)

AtmosAir with SAT T-Series MRS33 lb (15 kg)

AtmosAir with SAT T-Series A and AR MRS40 lb (18 kg)

AtmosAir with SAT V-Series MRS33 lb (15 kg)

AtmosAir with SAT APOD 2533 lb (15 kg)

AtmosAir with SAT APOD 25A and AR MRS40 lb (18 kg)

Dimensions (varies depending on model)**

Minimum Length75 in (190.5 cm)

Maximum Length84 in (213 cm)

Minimum Width32 in (81.5 cm)

Maximum Width38 in (96.5 cm)

Height7 in (17.5 cm)

Pump (UL Listed)

Weight5 lb (2.5 kg)

Dimensions 10.25 in (26 cm) L x 4.5 in (11.5 cm) W x 3.75 in (9.5 cm) H

Electrical:

Volts110 VAC

Frequency60 Hz

Max Amps0.075 Amps

Power Cord12 ft

Output7 liters per minute

Cycle Time5 minutes

* Patient weight capacity may vary depending on frame use. Please consult frame manufacturer.

** Additional sizes are available. *AtmosAir* 4000 and 9000 series up to 42in (106.7cm) with *Reliant* IS² cover and *AtmosAir* 9000 with Premium cover. Contact Arjo for more information.

Electromagnetic Compatibility

This section pertains only to A and AR models.

Although this equipment conforms with the intent of the directive 89/336/EEC in relation to electromagnetic compatibility (EMC), all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.

Portable and mobile RF communications equipment can effect medical electrical equipment.

Radios, cell phones and similar devices may affect this equipment and should be kept at least 6.5 ft (2 m) away from the equipment.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the following tables.

The following tables document compliance levels and guidance from the IEC 60601-1-2; 2007 Standard, for the electromagnetic environment in which the *AtmosAir* MRS should be used in a clinical environment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The <i>AtmosAir</i> MRS is intended for use in the electromagnetic environment specified below. The customer or user of the <i>AtmosAir</i> MRS should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	The <i>AtmosAir</i> MRS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions IEC 61000-3-3	yes	

Guidance and manufacturer's declaration - electromagnetic immunity

The *AtmosAir* MRS is intended for use in the electromagnetic environment specified below.
The customer or user of the *AtmosAir* MRS should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV Contact ±8 kV Air	±6 kV Contact ±8 kV Air	In accordance with IEC 60601-1-2:2007, floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	±1 kV ±2 kV	±1kV ±2 kV	N/A
Surge IEC 61000-4-5	1 kV line(s) to Line(s) 2 kV line(s) to Earth	1 kV line(s) to Line(s) 2 kV line(s) to Earth	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% half cycle 40% for 5 cycle 70% for 25 cycle 5% for 5 seconds	5% half cycle 40% 5 cycles 70% 25 cycles 5% for 5 seconds	N/A
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/M	3 A/M	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the *AtmosAir* MRS

The *AtmosAir* MRS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *AtmosAir* MRS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *AtmosAir* MRS as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter meters		
	150 kHz to 80 MHz not applicable	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.37	0.74
1	N/A	1.2	2.3
10	N/A	3.7	7.4
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separate distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1, At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity			
The <i>AtmosAir</i> MRS is intended for use in an electromagnetic environment specified below. The customer or user of the <i>AtmosAir</i> MRS should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3Vrms 150K - 80 MHz	3 Vrms 150K - 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the <i>AtmosAir</i> MRS, including cables, than the recommended separation distance calculated from the equation application to the frequency of the transmitter. Recommended Separation Distance Battery Operated Device
Radiated RF IEC 61000-4-3	3 Vrms 80 MHz - 2.5 GHz	3 Vrms 80 MHz - 2.5 GHz	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (see note a) should be less than the compliance level in each frequency range (see note b). Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1, At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the <i>AtmosAir</i> MRS is used exceeds the applicable RF compliance level above, the <i>AtmosAir</i> MRS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the <i>AtmosAir</i> MRS. b) Over the frequency range 150kHz, field strengths should be less than $[V1]$ V/m.			

***AtmosAir* products listed in this User Guide may or may not be applicable to CE markings. To verify applicable CE marking for a particular model / series, refer to that product's mattress tags.**

Symbols Used

Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

Important Operational Information

Foot End

Warning of possible hazard to system, patient or staff

Consult User Guide

Power

Manufacturer

= Safe Working Load

Wash at 60°C, 95°C Max for 15 minutes

Do Not Iron

Use Solution diluted to 1000ppm of available chlorine

Double Insulated

IPX0

Protection from ingress of fluids

Type BF Applied Part

Possible electrical shock hazard

Attention - See User Guide

On

Off

CE marking indicating conformity with European Community harmonised legislation. Figures indicate Notified Body supervision.

Wipe down Only

Tumble Dry at 60°C, 80°C Max

Do not use Phenol-based cleaning solutions

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

The operator must read this document (this instruction for use) before use. Note: the symbol is blue on the product label.

Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative, or visit www.arjo.com.

Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelize, 329 PB02 Galpão - Lapa
São Paulo - SP - Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

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Intentionally left blank

Intentionally left blank

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

CE
2797