INSTRUCTIONS FOR USE

# AtmosAir


## Mattress Replacement System

AtmosAir with SAT 4000 series

AtmosAir with SAT 9000 series

AtmosAir with SAT T-series

AtmosAir with SAT A-series

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Bedienungsanleitung · Gebruikershandleiding · Mode d'emploi  
Istruzioni per l'uso · Instrucciones de uso · Brugsvejledning · Bruksanvisning  
Manual do usuário · Kullanıcı Kılavuzu · Εγχειρίδιο Χρήστη · Käyttöohje  
Bruksanvisning · 사용자 매뉴얼

## **WARNING**

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.

**Mandatory to read the *Instructions for Use*.**

® and ™ are trademarks belonging to the Arjo group of companies

© Arjo 2019. SAT™ (Self Adjusting Technology™) is a trademark of MPTI, Ltd. Total Care and VersaCare are trademarks of Hill-Rom and / or its affiliates.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

ARJO HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJO PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY ARJO SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL ARJO BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJO TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties except as set forth in the written limited warranty included in this publication or with this product. Information in this publication may be subject to change at any time. Contact Arjo for updates.

### Important Information For Users

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.

**Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.**

Intentionally left blank

# Table Of Contents

<b>Introduction.....</b>	<b>3</b>
Indications .....	3
Contraindications.....	3
Risks and Precautions.....	3
Safety Information .....	4
Serious Incident .....	5
<b>Preparation for Use .....</b>	<b>6</b>
<b>Mattress Installation.....</b>	<b>6</b>
<i>AtmosAir</i> with SAT 4000, 9000, <i>AtmosAir</i> with SAT T-Series* and <i>AtmosAir</i> with SAT A-Series .....	6
Patient Placement and Nursing Care .....	7
CPR.....	7
Skin Care.....	7
Incontinence / Drainage .....	7
General Operation.....	7
<b>Care and Cleaning .....</b>	<b>8</b>
Fire Barrier .....	8
Cover Cleaning Options .....	8
<b>Preventive Maintenance Schedule.....</b>	<b>9</b>
Daily Cleaning .....	9
Inspection / System Check-Out .....	9
<b>Troubleshooting .....</b>	<b>9</b>
<b>Parts Diagram - <i>AtmosAir</i> with SAT 4000.....</b>	<b>10</b>
<b>Parts Diagram - <i>AtmosAir</i> with SAT 9000 and T-Series .....</b>	<b>11</b>
<b>Parts Diagram - <i>AtmosAir</i> with SAT A-Series .....</b>	<b>12</b>
<b>Replacement Parts+ .....</b>	<b>13</b>
<b>Specifications .....</b>	<b>16</b>
<b>Symbols Used.....</b>	<b>17</b>
<b>Customer Contact Information.....</b>	<b>17</b>

Intentionally left blank

# Introduction

It is recommended that all sections of this User Guide be read prior to product use.

Carefully review the **Contraindications, Safety Information and Risks and Precautions** sections prior to placing a patient on any AtmosAir® with SAT™ Mattress Replacement System (MRS).

Caregivers should review this information with the patient and the patient's family and / or legal guardian. Save this User Guide in an easily accessible location for quick reference.

This User Guide is comprised of the following *AtmosAir* MRS series:

- *AtmosAir* with SAT 4000 Series
- *AtmosAir* with SAT 9000 Series
- *AtmosAir* with SAT T-Series
- *AtmosAir* with SAT A-Series

The *AtmosAir* MRS comes in a wide variety of lengths and widths to support various customer requirements. Contact an Arjo representative for additional product information.

## Indications

- The *AtmosAir* MRS is intended to treat and prevent all categories of pressure ulcers.

## Contraindications

- unstable vertebral fracture
- cervical and skeletal traction

## Risks and Precautions

**Transfer** – Standard precautions should be taken during patient transfer.

**Side Rails and Restraints – WARNING:** Use or non-use of restraints, including side rails, can be critical to patient safety. Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. **See related Safety Information.**

**Patient Migration** – Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. **Monitor patients frequently to guard against patient entrapment.**

**Oxygen Use - DANGER:** Risk of explosion if the *AtmosAir* pump is used in the presence of flammable anesthetics. Use of this products pump in an oxygen-enriched environment may produce potential of fire hazard. This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. Unplug and do not use pump when using oxygen-administering equipment other than the nasal mask or half-bed-length tent type.

**Shock Hazard** - Electrical shock hazard; do not remove pump case covers. Refer to qualified service personnel.

## Safety Information

**Patient Entrance / Exit** – Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

**Brakes** – Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

**Bed Height** – To minimize risk of falls or injury, the bed should always be in the lowest practical position when the patient is unattended.

**Bed Frame** – Always use a standard healthcare bed frame with this mattress, with any safeguards or protocols that may be appropriate. Bed frame and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. It is recommended that bed and side rails (if used) comply with all applicable regulations and protocols.

**CPR** - Level the bed. Lower side rails and initiate CPR per facility protocols. Consider use of backboard if indicated. After CPR remove backboard, if used, raise siderails and reconfigure bed and accessories as in initial placement.

**Head of Bed Elevation** – Keep head of bed as low as possible to help prevent patient migration.

**Side Rails / Patient Restraints** - Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patients frequently to guard against patient entrapment.

**CAUTION:** When selecting a standard mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

**I.V. and Drainage Tubes** – I.V. and drainage tubes should always have slack for alternating pressure or rotation and other patient movements.

**Skin Care** – Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

**No Smoking in Bed** – Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

**General Protocols** – Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

**Fire Barrier Contact** - Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. If the fire barrier becomes damaged or soiled it should be replaced.

#### **End of Life Disposal –**

- Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.
- Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.
- Pump units have electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

### **Serious Incident**

**If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.**

## Preparation for Use

For information concerning the bed frame refer to the manufacturers' User Guide.

1. Open shipping container(s).

***Do not use sharp instruments to open boxes. Damage to mattress could result.***

2. Remove *AtmosAir* with SAT MRS from plastic protective cover.

**The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to accommodate; see Troubleshooting for more information. Wrinkles will not affect inflation or function, so mattress may be used immediately if needed.**

3. Check mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing mattress onto a new frame or for a new patient, check mattress surface for staining and soiling; clean and / or disinfect as required (see **Care and Cleaning**).
5. Level bed and lock brakes.
6. Remove existing mattress from bed frame.

## Mattress Installation

### *AtmosAir* with SAT 4000, 9000, *AtmosAir* with SAT T-Series\* and *AtmosAir* with SAT A-Series

1. Position mattress on bed frame with logo facing up and product information tags at foot end of bed.

**Magnets are integrated on the *AtmosAir* T-Series bottom cover to help keep the mattress in place during foot section extension / retraction.**

2. Ensure there are no gaps between mattress and bed frame or side rails.

***Always use a standard healthcare bed frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.***

\* T-Series is compatible for use with the Hill-Rom Total Care™ bed frame.

## Patient Placement and Nursing Care

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications**, **Safety Information** and **Risks and Precautions** sections prior to placing a patient on any *AtmosAir* MRS.

1. Transfer patient following all applicable safety rules and institution protocols.
2. Center patient side-to-side and head-to-foot on *AtmosAir* MRS surface.
3. Ensure all sections of the mattress fully support the patient.

**Mattress handles are for mattress transport only.**

## CPR

1. Level bed.
2. Lower or remove side rail on caregiver's side if necessary.
3. Begin CPR as per facility protocol. Consider use of backboard if indicated.
4. After CPR is performed:
  - Remove backboard if used.
  - Raise or install side rail as necessary.
  - Reconfigure bed and accessories as in initial placement.

## Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under patient are not wrinkled.

## Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace bed linens as required (see **Care and Cleaning**).

## General Operation

Avoid contact of sharp instruments with the *AtmosAir* MRS. Punctures, cuts and tears may prevent proper inflation and air pressure maintenance.

**Figure 4**

## Care and Cleaning

The following processes is recommended, but should be adapted to comply with local institution protocols. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The *AtmosAir* MRS should be routinely decontaminated between patients and at regular intervals while in use.

***Do not use Phenol-based solutions or abrasive compounds or pads on the cover during decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.***

***Covers should be wiped down immediately after it has been exposed to fluids or spills.***

### Fire Barrier

Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. Fire barrier should be replaced if soiled or damaged.

### Cover Cleaning Options

#### Sewn *AtmosAir* MRS (Top cover can not be removed from base)

1. Remove or push bed linens to center of mattress.

***The sewn *AtmosAir* MRS cover is cleaned by wipe down methods only. Do not launder as damage to the base may occur.***

2. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry surface with towel.
5. Ensure bed linens are refitted and not wrinkled under patient.
6. Clean pump and tubes (if needed) by wiping with a damp cloth.

#### RF Welded *AtmosAir* MRS (Detachable top cover)

1. Push bed linens to center of mattress to wipe down, remove to launder top cover or wipedown as described below.

***The RF Welded *AtmosAir* MRS has a detachable top cover that can be laundered. However the base is cleaned by wipe down methods only. Do not launder base as damage may occur.***

2. Wipe and rinse any soiling from the mattress surface and base. Use 1000 ppm chlorine or 70% alcohol.
3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry surface with towel.
5. Ensure bed linens are refitted and not wrinkled under patient.

Laundering Detachable Top Cover

- 1.    Unzip the top cover from the base for laundering.  
      Do not launder base as damage may occur.
- 2.    Recommended wash temperature for top cover is 60° C (140° F) for 15 minutes.
- 3.    Maximum wash temperature is 95° C (203° F) for 15 minutes.
- 4.    Tumble dry at 60° C (140° F) or air dry.
- 5.    Maximum drying temperature 80° C (176° F).

Preventive Maintenance Schedule

Preventive maintenance for the *AtmosAir* MRS consists of regular cleaning (see **Care and Cleaning**) and an overall system check-out to be performed at the intervals described below.

**All components must be cleaned, disinfected and inspected after each patient’s use and before use by a new patient.** Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

Daily Cleaning

The cover should be wiped daily with a mild soap and water solution. After wipe down, dry surface with a towel.

Inspection / System Check-Out

Check each of the following before placing the *AtmosAir* MRS with a new patient:

- 1.    Check mattress surface for tears or cracking; do not use if tears or cracks are present.
- 2.    Ensure mattress is free of stains and is not overly faded.

Troubleshooting

Do not attempt troubleshooting outside this guide or where the solution recommends to contact Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.


SYMPTOM	POSSIBLE CAUSE	SOLUTION
Mattress too firm upon arrival.	Difference in altitude not sufficient to open valves.	Apply weight to mattress to open valves.
Mattress cover too wrinkled upon removal from shipping container.	Internal components have not accommodated to environment. This does not affect inflation or function.	Let mattress accommodate for 24 hours. If problem continues, contact Arjo for assistance.

## Parts Diagram - *AtmosAir* with SAT 4000

The fire barrier is a sleeve that fits over the foam and SAT assemblies.

The SAT Cover, SAT Cells and SAT Interliner are all part of the integrated SAT System and cannot be ordered separately. See Replacement Parts+ for a complete list of SAT Systems.


## Parts Diagram - *AtmosAir* with SAT 9000 and T-Series

The fire barrier is a sleeve that fits over the foam and SAT assemblies.

The SAT Cover, SAT Cells and SAT Interliner are all part of the integrated SAT System and cannot be ordered separately. See Replacement Parts+ for a complete list of SAT Systems..


## Parts Diagram - *AtmosAir* with SAT A-Series

All hoses were removed to improve diagram readability.

The fire barrier is a sleeve that fits over the foam and SAT assemblies.

The **SAT Cover**, **SAT Cells** and **SAT Interliner** are all part of the integrated **SAT System** and cannot be ordered separately. See **Replacement Parts+** for a complete list of **SAT Systems**..

## Replacement Parts+

Do not attempt troubleshooting, maintenance or parts replacement outside this manual or where the solution recommends contacting Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Replaceable *AtmosAir* components are listed below. For more information such as pricing or additional spare parts that are not on this list, please contact your local Arjo representative.

### ATMOSAIR MATTRESS COMPONENTS

#### *AtmosAir* Foam Shells

PART TYPE.....	CATALOG NUMBER
4000 / 9000 Foam Shell .....	KA4/9FS(SIZE)
T-Series Foam Shell .....	KA9FSTC
A-Series Foam Shell .....	KAA-FOAMSHELL

#### *AtmosAir* with SAT System

PART TYPE.....	CATALOG NUMBER
4000 SAT System (Fits mattresses 32-38 in wide, 75-84 in length) .....	KA4SAT
9000 SAT System (Fits mattresses 32-38 in wide, 75-84 in length) .....	KA9SAT
T-Series SAT System (Fits mattresses 36 in wide, 75-84 in length) .....	KA9-T-SAT
A-Series SAT System (Fits mattresses 32-38 in wide, 75-84 in length) .....	KA9SAT

#### *AtmosAir* Fire Barrier Lining

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> Fire Barrier Lining, (Standard on all models) .....	KA4/9-FBL
RF Welded <i>AtmosAir</i> Fire Barrier Lining (Standard on all <i>AtmosAir</i> RF Welded mattresses [all models]) .....	AtmosAir-FBL

### ATMOSAIR MATTRESS COVERS\*

#### *AtmosAir* Model 4000 Mattress

PART TYPE.....	CATALOG NUMBER
Complete Cover, Soflux™ Top, Non-Skid Bottom .....	A4CSONB(WWWLL)NXS
Complete Cover, Soflux Top, Vinyl Gray Bottom .....	A4CSOVG(WWWLL)NXS
Complete Cover, Reliant IS²™ Top, Non-Skid Bottom .....	A4CRLNB(WWWLL)TXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom .....	A4CBFNB(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom, No Handles .....	A4CBFNH(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom .....	A4CBFVG(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom, No Handles .....	A4CBFVH(WWWLL)PXS
Complete Cover, Reliant IS² Top, Vinyl Gray Bottom .....	A4CRLVG(WWWLL)TXS
Complete Cover, Reliant IS² Top, Vinyl Gray Bottom, No Handles .....	A4CRLVH(WWWLL)TXS
Complete Cover, Reliant IS², RF Welded Top, Non-Skid PU Bottom, No Handles .....	A4CRFPN(WWWLL)TXS
Complete Cover, Reliant IS², RF Welded Top, Non-Skid PU Bottom, With Handles .....	A4CRFPU(WWWLL)TXS
Top Cover, Reliant IS², RF Welded .....	A4CRFXX(WWWLL)TXS
Base Cover, PU Non-skid Bottom, No Handles .....	A4CXXPN(WWWLL)XXS
Base Cover, PU Non-skid Bottom, With Handles .....	A4CXXPU(WWWLL)XXS

\* WWW = width in 0.5 in increments, LL = Length in 1.0 in increments.

**AtmosAir Model 9000 Mattress**

<b>PART TYPE</b>	<b>CATALOG NUMBER</b>
Complete Cover, <i>Soflux</i> Top, Non-Skid Bottom	A9CSONB(WWWLL)NXS
Complete Cover, <i>Soflux</i> Top, Non Skid Bottom, No Handles	A9CSONH(WWWLL)NXS
Complete Cover, <i>Soflux</i> Top, Vinyl Gray Bottom	A9CSOVG(WWWLL)NXS
Complete Cover, <i>Soflux</i> Top, Vinyl Gray Bottom, No Handles	A9CSOVH(WWWLL)NXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> Top, Non-Skid Bottom	A9CRLNB(WWWLL)TXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> Top, Non Skid Bottom, No Handles	A9CRLNH(WWWLL)TXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> Top, Vinyl Gray Bottom	A9CRLVG(WWWLL)TXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> Top, Vinyl Gray Bottom, No Handles	A9CRLVH(WWWLL)TXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom	A9CBFNB(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Non-Skid Bottom, No Handles	A9CBFNH(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom	A9CBFVG(WWWLL)PXS
Complete Cover, Barrier Flex 1 Top, Vinyl Gray Bottom, No Handles	A9CBFVH(WWWLL)PXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> , RF Welded Top, Non-Skid PU Bottom, No Handles	A9CRFPN(WWWLL)TXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> , RF Welded Top, Non-Skid PU Bottom, With Handles	A9CRFPU(WWWLL)TXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Bottom, No Handles	A9CPRPN(WWWLL)NXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Bottom, with Handles	A9CPRPU(WWWLL)NXS
Top Cover, <i>Reliant IS<sup>2</sup></i> , RF Welded	A9CRFXX(WWWLL)TXS
Top Cover, Premium RF Welded	A9CPRXX(WWWLL)NXS
Base Cover, PU Non-skid Bottom, No Handles	A9CXXPN(WWWLL)XXS
Base Cover, PU Non-skid Bottom, With Handles	A9CXXPU(WWWLL)XXS

**AtmosAir - T-Series**

<b>PART TYPE</b>	<b>CATALOG NUMBER</b>
T-Series Mattress - Complete Cover, <i>Soflux</i> top and Non-Skid base	ATCSONB35084NXS
T-Series Mattress - Complete Cover, <i>Reliant IS<sup>2</sup></i> top and Non-Skid base	ATCRLNB35084NXS

**AtmosAir - A-Series**

<b>PART TYPE</b>	<b>CATALOG NUMBER</b>
A-Series Mattress - Complete Cover, <i>Reliant IS<sup>2</sup></i> top and Vinyl Gray base	AACRLVG(WWWLL)TXS
A-Series Mattress - Complete Cover, <i>Reliant IS<sup>2</sup></i> RF welded top, Non-Skid PU Base, No Handles	AACRFPN(WWWLL)TXS
A-Series Mattress - Top cover <i>Reliant IS<sup>2</sup></i> RF welded	AACRFXX(WWWLL)TXS
A-Series Mattress - Base cover, PU Non-Skid	AACXXPN(WWWLL)TXS
Complete Cover, <i>Reliant IS<sup>2</sup></i> RF Welded Top, Non-Skid PU Base, With Handles	AACRFPN(WWWLL)TXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Base, No Handles	AACPRPN(WWWLL)NXS
Complete Cover, Premium RF Welded Top, Non-Skid PU Base, With Handles	AACPRPU(WWWLL)NXS

Top Cover, Premium RF Welded .....	AACPRXX(WWWLL)NXS
Base Cover, Non-Skid PU, No Handles .....	AACXXPN(WWWLL)XXS
Base Cover, Non-Skid PU, With Handles.....	AACXXPU(WWWLL)XXS

**ATMOSAIR LITERATURE**

<b>PART TYPE.....</b>	<b>CATALOG NUMBER</b>
<i>AtmosAir</i> MRS User Guide.....	407384-AH
<i>AtmosAir</i> MRS User Guide (EMEA) .....	407391-EE-AH

# Specifications

Specifications subject to change without notice.

## Mattresses

<b>Maximum Recommended Patient Weight*</b>	<b>500 lb (227 kg)</b>
<b>Weight (based on 35 in x 80 in x 7 in MRS, other sizes will vary)</b>	
<i>AtmosAir</i> with SAT 4000 MRS	30 lb (14 kg)
<i>AtmosAir</i> with SAT 9000 MRS	33 lb (15 kg)
<i>AtmosAir</i> with SAT T-Series MRS	33 lb (15 kg)
<i>AtmosAir</i> with SAT A-Series MRS	33 lb (15 kg)

## Dimensions (varies depending on model)\*\*

Minimum Length	75 in (190.5 cm)
Maximum Length	84 in (213 cm)
Minimum Width	32 in (81.5 cm)
Maximum Width	38 in (96.5 cm)
Height	7 in (17.5 cm)

\* Patient weight capacity may vary depending on frame use. Please consult frame manufacturer.

\*\* Additional sizes are available. Contact Arjo for more information. *AtmosAir* 4000 and 9000 Series up to 42 in (106.7 cm) with **Reliant IS<sup>2</sup>** cover, and *AtmosAir* 9000 with Premium cover.

## Symbols Used

Important Operational Information

Foot End

Warning of possible hazard to system, patient or staff

Consult User Guide

Manufacturer

Safe Working Load

Max 95°C  
15 min.

Wash at 60°C, 95°C Max for 15 minutes

Serial Number

Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

Do Not Iron

Use Solution diluted to 1000ppm of available chlorine

Attention - See User Guide

CE marking indicating conformity with European Community harmonised legislation


Wipe down Only

Tumble Dry at 60°C, 80°C Max

Do not use Phenol-based cleaning solutions

The operator must read this document (this instruction for use) before use. Note: the symbol is blue on the product label.

## Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative, or visit [www.arjo.com](http://www.arjo.com).

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Intentionally left blank

AUSTRALIA  
Arjo Australia Pty Ltd  
78, Forsyth Street  
O'Connor  
AU-6163 Western Australia  
Tel: +61 89337 4111  
Free: +1 800 072 040  
Fax: +61 89337 9077

BELGIQUE / BELGIË  
Arjo NV/SA  
Evenbroekveld 16  
BE-9420 ERPE-MERE  
Tél/Tel: +32 (0) 53 60 73 80  
Fax: +32 (0) 53 60 73 81  
E-mail: info.belgium@arjo.be

BRASIL  
Arjo Brasil Equipamentos  
Médicos Ltda  
Rua Marina Ciufulli Zanfelice,  
329 PB02  
Galpão - Lapa  
São Paulo - SP - Brasil  
CEP: 05040-000  
Phone: 55-11-3588-5088  
E-mail: vendas.latam@arjo.com  
E-mail: servicios.latam@arjo.com

CANADA  
Arjo Canada Inc.  
90 Matheson Boulevard West  
Suite 300  
CA-MISSISSAUGA, ON, L5R 3R3  
Tel/Tél: +1 905 238 7880  
Free: +1 800 665 4831 Institutional  
Free: +1 800 868 0441 Home Care  
Fax: +1 905 238 7881  
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA  
Arjo Czech Republic s.r.o.  
Na Strži 1702/65  
140 00 Praha  
Czech Republic  
Phone No: +420225092307  
e-mail: info.cz@arjo.com

DANMARK  
Arjo A/S  
Vassingerødvej 52  
DK-3540 LYNGE  
Tel: +45 49 13 84 86  
Fax: +45 49 13 84 87  
E-mail: dk\_kundeservice@arjo.com

DEUTSCHLAND  
Arjo GmbH  
Peter-Sander-Strasse 10  
DE-55252 MAINZ-KASTEL  
Tel: +49 (0) 6134 186 0  
Fax: +49 (0) 6134 186 160  
E-mail: info-de@arjo.com

ESPAÑA  
Arjo Ibérica S.L.  
Parque Empresarial Rivas Futura,  
C/ Marie Curie 5  
Edificio Alfa Planta 6 oficina 6.1.-62  
ES-28521 Rivas Vacia, MADRID  
Tel: +34 93 583 11 20  
Fax: +34 93 583 11 22  
E-mail: info.es@arjo.com

FRANCE  
Arjo SAS  
2 Avenue Alcide de Gasperi  
CS 70133  
FR-59436 RONCQ CEDEX  
Tél: +33 (0) 3 20 28 13 13  
Fax: +33 (0) 3 20 28 13 14  
E-mail: info.france@arjo.com

HONG KONG  
Arjo Hong Kong Limited  
Room 411-414, 4/F,  
Manhattan Centre,  
8 Kwai Cheong Road,  
Kwai Chung, N.T.,  
HONG KONG  
Tel: +852 2960 7600  
Fax: +852 2960 1711

ITALIA  
Arjo Italia S.p.A.  
Via Giacomo Peroni 400-402  
IT-00131 ROMA  
Tel: +39 (0) 6 87426211  
Fax: +39 (0) 6 87426222  
E-mail: Italy.promo@arjo.com

MIDDLE EAST  
Arjo Middle East FZ-LLC  
Office 908, 9th Floor,  
HQ Building, North Tower,  
Dubai Science Park,  
Al Barsha South  
P.O Box 11488, Dubai,  
United Arab Emirates  
Direct +971 487 48053  
Fax +971 487 48072  
Email: Info.ME@arjo.com

NETHERLAND  
Arjo BV  
Biezenwei 21  
4004 MB TIEL  
Postbus 6116  
4000 HC TIEL  
Tel: +31 (0) 344 64 08 00  
Fax: +31 (0) 344 64 08 85  
E-mail: info.nl@arjo.com

NEW ZEALAND  
Arjo Ltd  
34 Vestey Drive  
Mount Wellington  
NZ-AUCKLAND 1060  
Tel: +64 (0) 9 573 5344  
Free Call: 0800 000 151  
Fax: +64 (0) 9 573 5384  
E-mail: nz.info@Arjo.com

NORGE  
Arjo Norway AS  
Olaf Helset's vei 5  
N-0694 OSLO  
Tel: +47 22 08 00 50  
Faks: +47 22 08 00 51  
E-mail: no.kundeservice@arjo.com

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Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB  
Hans Michelsensgatan 10  
211 20 Malmö, Sweden  
[www.arjo.com](http://www.arjo.com)

**arjo**

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Arjo Inc.