
INSTRUCTIONS FOR USE

Toilet Slings

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使用说明书 · Bedienungsanleitung · Mode d'emploi · Istruzioni per l'uso · Gebruiksaanwijzing

Contents

EN

Foreword	2	Loop Attachment & Detachment.....	8
Customer Support.....	2	Attach the Loops (5 Steps).....	8
Definitions in this IFU	2	Detach the Loops (2 steps).....	8
Intended Use	3	Applying the Sling.....	9
Patient/Resident Assessment	3	In Bed (24 steps).....	9
Expected Service Life	3	In Chair/ Wheelchair (24 steps).....	10
Safety Instructions	4	On Floor	10
Safety Practises	4	Removing the Sling	11
Serious incident.....	4	In Bed (14 steps).....	11
Select Sling Size.....	6	In Chair/Wheelchair (13 steps).....	11
Using the Arjo Measuring Tape (4 steps) (Accessory)	6	Cleaning and Disinfection.....	12
Sizing Without a Measuring Tape (2 steps).....	6	Cleaning Instruction (7 steps)	12
Sling Selection.....	6	Cleaning Chemicals	12
Transfer/Transport.....	6	Disinfection.....	12
Allowed Combinations.....	7	Care and Preventive Maintenance	12
Safe Working Load (SWL).	7	Before and After Every Use.....	12
Clip Attachment & Detachment.....	8	When Soiled or Stained and Between Patients	12
Attach the Clips (5 steps).....	8	Storage.....	12
Detach the Clips (3 steps).....	8	Service and Maintenance.....	12

WARNING

To avoid injury, always read this **Instructions for Use** and accompanied documents before using the product. Mandatory to read the **Instructions for Use**.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies. © Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

Foreword

Thank you for purchasing Arjo equipment. Please read this *Instructions For Use (IFU)* thoroughly!

Arjo will not be held responsible for any accidents, incidents or lack of performance that occur as a result of any unauthorized modification to its products.

Customer Support

If you require further information, please contact your local Arjo representative. The contact information is located at the end of this *IFU*.

Definitions in this IFU

WARNING

Means: Safety warning. Failure to understand and obey this warning may result in injury to you or to others.

CAUTION

Means: Failure to follow these instructions may cause damage to all or parts of the system or equipment.

NOTE

Means: This is important information for the correct use of this system or equipment.

Intended Use

The *Toilet Sling* has been specially designed to help support patients/residents while toileting, when used solely in conjunction with Arjo patient/resident lifts fitted with a four point tilting frame.

Do not use the *Toilet Sling* for lifting and transportation apart from toilet visits.

The *Toilet Sling* shall only be used by appropriately trained caregivers with adequate knowledge of the care environment, and in accordance with the instructions outlined in the Instructions for Use.

The *Toilet Sling* is intended to be used in hospital environments, nursing homes, other healthcare facilities and home care.

The *Toilet Sling* should only be used for the purpose specified in this Instructions for Use. Any other use is prohibited.

Patient/Resident Assessment

We recommend that facilities establish regular assessment routines. Caregivers should assess each patient/resident according to the following criteria prior to use:

- Patient/resident sits in wheelchair
- Has no capacity to support herself/himself at all
- Cannot stand unsupported and is not able to bear weight, not even partially
- Dependent on caregiver in most situations
- Physically demanding for caregiver

- Stimulation of remaining abilities is very important
- Requires some trunk stability as a must for this type of sling

Safe Working Load (SWL) for *Toilet Sling*:

- MAA4031M, MLA4031, MLA4531: 272 kg (600lbs).

The right type and size of slings should be used after proper assessment of each patient/resident's size, condition and the type of lifting situation.

If the patient/resident does not meet these criteria an alternative equipment/system shall be used.

EN

Expected Service Life

The expected service life of *Toilet Sling* is the maximum period of useful life.

The expected life of the sling is dependent on the actual use conditions. Therefore, before use, always make sure that the sling does not show any sign of fraying, tearing or other damage and that there is no damage (i.e. cracking, bending, breaking). If any such damage is observed, do not use the sling. If you have any doubts about sling safety, as a precaution and to ensure safety, do not use the sling.

The expected service life of the *Toilet Slings*:

Service Time:

- 2 years

Shelf Time:

- 5 years

Safety Instructions

EN

WARNING

To avoid injury, make sure that the patient is not left unattended at any time.

WARNING

To avoid injury, only use the Toilet sling on patients with upper body strength. Lack of upper body strength can cause the patient to slide out.

WARNING

To avoid falling, make sure that the user weight is lower than the safe working load for all products or accessories being used.

WARNING

To avoid injury, only store the equipment for a short period of time. If stored longer than stated in the IFU, it can cause weakness and breakage in the material.

WARNING

To avoid injury, always assess the patient prior to use.

WARNING

To avoid injury, never let the patient smoke when using the sling. The sling can catch on fire.

WARNING

To avoid injury, keep the equipment away from sun/UV-light. Exposure to sun/UV-light can weaken the material.


Safety Practises

- At any time, if the patient becomes agitated, stop transferring/transporting and safely lower the patient.
- Patient with spasm can be lifted, but great care should be taken to support the patient's legs.

Serious incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Parts Designation

1. Sling label
(located outside of the sling)
2. Stiffener pocket
3. Stiffener
4. Neck buckle
5. Attachment clip (shoulder)
6. Pull handle
7. Leg flap
8. Leg strap
9. Attachment clip (leg)
10. Mid section
11. Midsection buckle
12. Head section
13. Support belt
14. Attachment loop (leg)
15. Attachment loop (shoulder)

EN

Example of slings. Not all models are represented.

Preparations

Before First Use (6 steps)

1. Check all parts of the sling, see section "Parts Designation" on page 5. If any part is missing or damaged - Do NOT use the sling.
2. Read this *IFU* thoroughly.
3. Check that the sling is clean.
4. Choose a designated area where the *IFU* should be kept and easily accessible at all times.
5. Make sure to have a rescue plan ready in case of a patient emergency.
6. For questions, contact your local Arjo representative for support.

NOTE

Arjo recommends the name of the patient is written on the sling care label to avoid cross contamination between patients.

Before Every Use (5 steps)

1. Patient assessment must always determine the methods used in everyday care, see "Intended Use" on page 3.

WARNING

To avoid the patient from falling, make sure to select the correct sling size according to the *IFU*.

2. Check the allowed combinations against the sling, spreader bar and lift, see "Allowed Combinations" on page 7.

WARNING

To avoid injury, always make sure to inspect the equipment prior to use.

3. Check all parts of the sling, see "Parts Designation" on page 5. If any part is missing or damaged - Do NOT use the sling. Check for:
 - fraying
 - loose stitching
 - tears
 - fabric holes
 - soiled fabric
 - damaged clips/loops
 - damaged buckles
 - unreadable or damaged label

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this *IFU*.

4. Check that the sling is clean. If the sling is not clean see "Cleaning and Disinfection" on page 12.
5. Before using the sling, read the lift *IFU* for transfer.

Select Sling Size

Using the Arjo Measuring Tape (4 steps) (Accessory)

The Arjo Measuring Tape is a sizing accessory only intended as a guide and is an approximation. Patient width and body proportions should also be considered when selecting the correct size of Arjo sling.

1. If possible, the patient should be in a seated position.
2. Place the point marked **A** on the measuring tape over the patient's coccyx/against seat (position A). (See Fig. 1) and (See Fig. 2)
3. Measure from coccyx/seat (position A) to the top of the head (position **B**) (See Fig. 1).

EN



Fig. 1

Fig. 2

Sling Selection

The patient's physical disabilities, weight distribution and general physique needs to be taken into consideration when selecting a sling.

Transfer/Transport

Article number	Safe Working Load	Product Description	Size
MAA4031M	272 (600 lbs)	Clip toilet with head support	XS, S, M, L, LL, XL, XXL
MLA4031	272 (600 lbs)	Loop toilet with head support	S, M, L, XL
MLA4531	272 (600 lbs)	Loop toilet w/o head support	S, M, L, XL

Allowed Combinations

WARNING

To avoid injury, always follow the allowed combinations listed in this *IFU*. No other combinations are allowed.

Safe Working Load (SWL)

Always follow the lowest SWL of the total system. E.g. the *Maxi Twin Compact lift/spreader bar* has a SWL of 160 kg (352 lbs) and the MAA4031 sling has a SWL of 272 kg (600 lbs). This means that the *Maxi Twin Compact lift /spreader bar* has the lowest SWL. The patient is not allowed to weigh more than the lowest SWL.

EN

	Lift	Maxi Twin Compact	Maxi Lite	Maxi 500	Maxi Twin	
	SWL	160 kg (352 lbs)	160 kg (352 lbs)	227 kg (500 lbs)	182 kg (401 lbs)	
	Spreader Bar (DPS)	Power Medium	Flat	Manual	Manual Medium	Power
Sling MAA 4031M	SWL 272kg (600 lbs)	Size S, M, L, LL, XL	Size S, M, L	Size S, M, L, LL	Size S, M, L, LL	Size S, M, L, LL, XL

	Lift	Maxi Move				Maxi Sky
	SWL	227 kg (500 lbs)				272 kg (600 lbs)
	Spreader Bar (DPS)	Manual Small	Manual Medium	Power Medium	Power Large	Manual Flat
Sling MAA 4031M	SWL 272kg (600 lbs)	Size XS, S	Size S, M, L, LL	Size S, M, L, LL, XL	Size L, LL, XL, XXL	Size S, M, L

	Lift	Minstrel	Maxi 500	Maxi Twin		
	SWL	190 kg (418 lbs)	227 kg (500 lbs)	182 kg (401 lbs)		
	Spreader Bar	4 point	2 point Loop	2 point Small	2 point Medium	4 point Large
Sling	SWL	Size	Size	Size	Size	Size
MLA4531M	272kg (600 lbs)	M, L, XL	S, M, L, XL	S, M	S, M, L, XL	M, L, XL
MLA4031	272kg (600 lbs)	M, L, XL	S, M, L, XL	S, M	S, M, L, XL	M, L, XL

	Lift	Maxi Move				Maxi Sky 440	Maxi Sky 600
	SWL	125 kg (275 lbs)		227 kg (500 lbs)		200 kg (440 lbs)	272 kg (600 lbs)
	Spreader Bar	2 point Small	2 point Medium	2 point loop	4 point loop	Fix	
Sling	SWL	Size	Size	Size	Size	Size	Size
MLA4531	272kg (600 lbs)	S, M	S, M, L, XL	S, M, L, XL	M, L, XL	S, M, L, XL	S, M, L, XL
MLA4031	272kg (600 lbs)	S, M	S, M, L, XL	S, M, L, XL	M, L, XL	S, M, L, XL	S, M, L, XL

Clip Attachment & Detachment

NOTE

If available on your lift device, reposition the spreader bar using the Dynamic Position System (DPS) for easier clip attachment of the sling. See lift IFU for instructions.

EN

Attach the Clips (5 steps)

1. Place the clip on the spreader bar lug. (See Fig. 3)
2. Pull the strap down.
3. Make sure the lug is locked at the top end of the clip. (See Fig. 4)

4. Make sure the strap is not squeezed in between the clip and the spreader bar.
5. Make sure the straps are not twisted.

Detach the Clips (3 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the clip.

1. Pull the strap up. (See Fig. 5)
2. Make sure the lug is unlocked at the bottom of the clip.
3. Remove the clip.

Loop Attachment & Detachment

2 point spreader bar: First attach the shoulder loops, then attach the leg loops.

4 point spreader bar: First attach the shoulder loops in the front hooks. Then attach the leg loops in the back hooks

WARNING

Do not cross the shoulder straps.

Attach the Loops (5 Steps)

WARNING

Use the same loop length (loop colour) for the shoulder straps and the same length (loop colour) for the leg straps.

1. Place the loop over the spring loaded latch (A). (See Fig. 6)
2. Pull the loop down to force the latch to open.
3. Make sure that the spring loaded latch closes completely with the loop inside. (See Fig. 6)
4. Make sure that the latch is moving freely.
5. Make sure loops and straps are not twisted.

NOTE

The loops lengths on the straps are used to adjust the patient's position and comfort.

- Loops closer to the shoulders; more forward position
- Loops further away from the shoulders; more reclined position
- Loops closer to the legs; raised legs/bottom
- Loops further away from legs; lowered legs/bottom

Detach the Loops (2 steps)

Make sure that the weight of the patient is taken up by the receiving surface before removing the loop.

Method 1

1. Press the spring loaded latch open.
2. Remove the loop.

Method 2

1. Pull one side of the loop over the hook and spring loaded latch. (See Fig. 7)
2. Pull down the loop. (See Fig. 8)

Fig. 3


Fig. 4


Fig. 5


Fig. 6

Fig. 7



Fig. 8

Applying the Sling

In Bed (24 steps)

For transfer see respective lift *IFU*.

1. Place the lift nearby.
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. The patient needs to be in a seated position in bed.
4. Check that the stiffeners are completely inside the stiffener pockets, if any.
5. Prepare the patient for transfer/toileting by loosening the patient's clothing.
6. Lean the patient forward.
7. Place the inside of the sling against the patients back. The sling label is located on the outside. If needed, use an *Arjo Sliding Sheet/Tube* for placing the sling. See respective sliding sheet/tube *IFU*.
8. Make sure that the centre line of the sling is aligned with the patient's spine and coccyx.
9. Lean the patient back.
10. Attach the support belt around the patient's midsection, if any. This prevents the patient from falling forward.
11. Place the leg flaps underneath the patient's legs.
12. Make sure that the: (See Fig. 9)
 - sling is centred and flat without creases,
 - sling's head support covers the neck/head area,
 - sling pieces are not twisted underneath the patient and
 - patient's arms are outside the sling.
13. Apply the brakes on the lift.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

14. Attach the sling, start with the midsection and then the legs.
 - If using a loop sling, attach the loops in a more seated position.
 - If using a clip sling, position the DPS (Dynamic Positioning System) in a seated position.

15. Slightly lift the patient to create tension in the sling. (See Fig. 10)

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

EN

16. Make sure that:
 - All clips/loops are securely attached (See Fig. 4) and (See Fig. 6)
 - All straps are straight (not twisted)
 - The patient lays comfortably in the sling.

17. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip/loop.
18. Attach the neck buckle to midsection of sling, if any.
19. Check that the patient is not sliding out through the bottom of the sling. If the patient is sliding out, replace it with a smaller sling.
20. Continue lifting the patient.
21. Transfer the patient to the toilet according to the lift *IFU*.
22. Do NOT apply any lift brakes.
23. Push the sling/straps and not the patient towards the back of the toilet for better positioning.
24. Lower the patient on the toilet seat. See the lift *IFU*. (See Fig. 11)

Fig. 9

Fig. 10

Fig. 11

In Chair/ Wheelchair (24 steps)


1. Place the lift nearby.
2. Apply the brakes on the wheelchair.
3. Check that the stiffeners are completely inside the stiffener pockets, if any.
4. Prepare the patient for transfer by loosening the patient's clothing.
5. Stand in front of the patient and lean the patient forward.
6. Place the inside of the sling against the patients back. The sling label is located on the outside. If needed, use an Arjo sliding sheet/tube for placing the sling See respective *Sliding sheet/tube IFU*.
7. Make sure that the centre line is aligned with the patient's spine and coccyx.
8. Lean the patient back.
9. Attach the support belt around the patient's midsection, if any. This prevents the patient from falling forward.
10. Place the leg flaps underneath the legs.
11. Make sure that the:
 - sling is centred and flat without creases,
 - sling's head support covers the neck/head area,
 - sling pieces are not twisted underneath the patient and
 - patient's arms are outside the sling.
12. Place the lift in front of the patient and open the lift chassis legs, if needed. (See Fig. 12)
13. Apply the brakes on the lift.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

14. Attach the sling, start with the midsection and then the legs.
 - If using a loop sling, attach the loops in a more seated position.
 - If using a clip sling, position the DPS (Dynamic Positioning System) in a seated position.

Fig. 12

WARNING

Make sure straps are not caught by wheelchair or lift castors.

15. Slightly lift the patient to create tension in the sling. (See Fig. 12)

WARNING

To avoid the patient from falling, make sure that the sling attachments are attached securely before and during the lifting process.

16. Make sure that:
 - All clips/loops are securely attached (See Fig. 4) and (See Fig. 6)
 - All straps are straight (not twisted)
 - The patient lays comfortably in the sling.
17. If an adjustment is needed, lower the patient and make sure that the weight of the patient is taken up by the receiving surface before removing the clip/loop.
18. Attach the neck buckle to midsection of sling, if any.

WARNING

To avoid injury, always make sure the sling is attached correctly.

19. Check that the patient is not sliding out through the bottom of the sling. If the patient is sliding out, replace it with a smaller sling.
20. Continue to lift the patient.
21. Transfer the patient to the toilet according to the lift *IFU*.
22. Do NOT apply any lift brakes.
23. Slightly push the sling/straps and not the patient, towards the back of the toilet for better positioning.
24. Lower the patient on the toilet seat. See the lift *IFU*. (See Fig. 13)

On Floor

Do NOT use *Toilet sling*. Use an appropriate sling to lift from the floor.

Fig. 13

Removing the Sling

In Bed (14 steps)

1. Position the patient over the bed. (See Fig. 14)
2. Apply bed brakes and adjust the height of the bed to an ergonomic position, if possible.
3. Adjust the bed to a seated position
4. If using a clip sling, position the DPS (Dynamic Positioning System) in a seated position.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

5. Lower the patient onto the bed. Make sure that the weight of the patient is taken up by the receiving surface before removing the clip/loop.
6. Detach the sling.
7. Unbuckle the head section from the midsection of the sling.
8. Move the lift away from the patient.
9. Detach the sling belt around the patient's midsection, if any.
10. Pull out the leg flaps from under the patient's legs.
11. Lean the patient forward.
12. Remove the sling. Use an Arjo sliding sheet/tube for removing the sling. See respective *Sliding sheet/tube IFU*.
13. Lean the patient back.
14. Make sure the patient is in a comfortable and secure position in the bed.

Fig. 14

In Chair/Wheelchair (13 steps)


1. Apply the brakes on the wheelchair.
2. If using a floor lift, open the lift chassis legs.
3. Position the patient over the chair/wheelchair.
4. If using a clip toilet sling, position the DPS (Dynamic Positioning System) in a seated position.
5. Slightly push the sling/straps and not the patient, towards the back of the chair/wheelchair for better positioning.

WARNING

To avoid injury to the patient, pay close attention when lowering or adjusting the spreader bar.

6. Lower the patient into the chair/wheelchair.
7. Make sure the patients lower back is positioned all the way in against the chair/wheelchair to avoid the patient from sliding out.
8. Detach the sling.
9. Move the lift away from the patient. (See Fig. 15)
10. Pull out the leg flaps from under the patient's legs and place them backwards along the patient's side.
11. Stand in front of the patient and lean the patient forward.
12. Pull out the sling. If needed, use an Arjo sliding sheet/tube. See respective *Sliding sheet/tube IFU*.
13. Lean the patient back up to a comfortable and secure position.

Fig. 15

Cleaning and Disinfection

WARNING

To prevent cross-contamination, always follow the disinfection instructions in this IFU.

WARNING

To avoid material damage and injury, clean and disinfect according to this IFU.

- No other chemicals are allowed.
- Never clean with chlorine.
- Chlorine will deteriorate the surface of the material.

WARNING

To avoid injury, always remove the sling before disinfecting the lift.

All Toilet slings should be cleaned when they are soiled or stained and between patients.

Cleaning Instruction (7 steps)

1. Detach the sling from the lift.
2. Remove the stiffeners, if any, before washing.
3. Close all buckles and hook and loop attachment on the slings and belts prior to washing.
4. Check the symbols on the product label for cleaning recommendations.

Care and Preventive Maintenance

WARNING

To avoid injury to both patient and caregiver, never modify the equipment or use incompatible parts.

Before and After Every Use

Visually check all exposed parts

The caregiver shall inspect the sling before and after every use. The complete sling should be checked for all deviations listed below. If any of these deviations are visible, replace the sling immediately.

- Fraying
- Loose stitching
- Tears
- Holes
- Discolouration or stains from bleaching
- Sling soiled or stained
- Label unreadable or damaged

5. Machine wash the sling at 70°C (158°F). In order to obtain sufficient level of disinfection, please wash according to local hygiene legislation.

If allowed according to the product label, tumble dry in low temperature, max 60°C (140°F)

6. Do NOT

- wash together with other items that have rough surfaces or sharp objects.
- use any mechanical pressure, pressing or rolling
- use bleach
- use gas sterilization
- use autoclave
- dry clean
- steam
- ironing

7. Place the stiffeners, if any, back in to the sling before use.

Cleaning Chemicals

For all sling material use common commercial detergent without optical brightener. No other chemicals are allowed e.g. chlorine, softener, iodine based disinfectants, bromine and ozone.

Disinfection

The only allowed disinfection is disinfection by washing.

When Soiled or Stained and Between Patients

Clean/Disinfect

The caregiver should make sure that the sling is cleaned according to "Cleaning and Disinfection" on page 12, when it is soiled or stained and between patients. If the soil and stains are still visible, replace the sling immediately.

Storage

When not in use, the slings should be stored away from direct sunlight where they are not subject to unnecessary strain, stress or pressure, or to excessive heat or humidity. The slings should be kept away from sharp edges, corrosives or other things that could cause damage on the sling.

Service and Maintenance

It is recommended that qualified personnel inspect the condition of the sling twice a year (every 6 month), according to ISO 10535.

Troubleshooting

Problem	Action
The patient is not positioned correctly in the sling (e.g. one strap is shorter than the other).	<ul style="list-style-type: none">• Make sure the clip/loop straps are not twisted.• Make sure the patient is seated on the inside of the sling. The outside has a label for reference.• Make sure the patient is centred in the sling.
The patient feels discomfort in the leg area when seated in the sling.	Make sure that there are no creases on the leg flaps of the sling.
The sling is difficult to apply (e.g. during log roll).	Do NOT use the log roll method. Use an Arjo sliding sheet/tube according to instructions in respective <i>IFU</i> . <ul style="list-style-type: none">• Place the patient in a more seated position in the bed or on the floor. Place a pillow behind the patient's neck/back, outside of the sling.• Raise the back rest of the bed, if available.• If applicable, make sure the clip strap is not caught between the lug on the lift and the clip attachment.• Adjust the spreader bar to a more reclined position. Pay attention to the patient's head.
There are difficulties applying the shoulder clips/loops.	<ul style="list-style-type: none">• Bend the patient's legs or place a pillow under the patient's knees.• If applicable, make sure the clip strap is not caught between the lug on the lift and the clip attachment.• Adjust the spreader bar to a more seated position.
There are difficulties applying the leg clips/loops.	<ul style="list-style-type: none">• Bend the patient's legs or place a pillow under the patient's knees.• If applicable, make sure the clip strap is not caught between the lug on the lift and the clip attachment.• Adjust the spreader bar to a more seated position.
It is difficult to use the sling in combination with the lift and spreader bar.	Use the sling only in accordance with the combinations described in section " <i>Allowed Combinations</i> " on page 7.
The position of the patient in the sling should be more reclined or seated.	Use the spreader bar or the loop lengths to adjust the patient's position.
The patient wishes to have extra support for the neck.	Place a pillow behind the patient's neck, inside the sling.
The patient is sliding out of the opening of the sling.	<ul style="list-style-type: none">• Make sure the patient's arms are outside of the sling.• If the problem still occurs: Select/choose a smaller size sling.

EN

Technical Specifications

EN

General	
Safe working load (SWL) = Maximum total load	See "Allowed Combinations" on page 7
Service life – Recommended period of use	2 years* * Refer to "Expected Service Life" on page 3.
Shelf life – Maximum period of storing new unpacked product	5 years
Model and Type	See "Sling Selection" on page 6
Operating, Transport and Storage Environment	
Temperature	0°C to +40°C (+32°F to +104°F) Operating and Storage -25°C to +60°C (-13°F to +140°F) Transport
Humidity	Max 15–70% at +20°C (68°F) Operating and Storage Max 10–95% at +20°C (68°F) Transport
End of Life Disposal	
Package	The bag consists of soft plastic, recyclable to local legislations.
Sling	Slings including stiffeners/stabilizers, padding material, any other textiles or polymers or plastic materials etc. should be sorted as combustible waste.
Parts and Accessories	
Stiffeners	GX21270, TBS106 and TBS107
Belt size: S, M & L	TSB.100-C
Belt size: XL	TSB.200-C
Measuring tape	MTA1000

Label on the Sling

EN

Care and Washing Symbols		Article Number	
	Machine wash 70°C (158°F)	REF XXXXXX-X Article number with -X refers to the sling size.	
	No bleaching allowed	REF XXXXXXX If the article number is stated without the size letter at the end, it means the sling is a One Size.	
	No tumble drying allowed		
	Tumble drying		
	No Ironing allowed	PES Polyester	
	No dry cleaning allowed	PE Polyethylene	
Certificates/Markings		PU Polyurethane	
	CE marking indicating conformity with European Community harmonised legislation		
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745		
Symbol for Clip and Loop		Misc. Symbols	
	Use a clip spreader bar.		Safe Working Load (SWL)
	Use a loop spreader bar.		Patient's name symbol
			Record symbol
			Read the IFU before use
			Article number
			Manufacturing date and year
			Manufacturer name and address

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NL

Intentionally left blank

AUSTRALIA
Arjo Australia
Building B, Level 3
11 Talavera Road
Macquarie Park, NSW, 2113, Australia
Phone: 1800 072 040

BELGIQUE / BELGIË

Arjo Belgium nv
Evenbroekveld 16
9420 Erpe-Mere
Belgium
T: +32 (0) 53 60 73 80
F: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.com

BRASIL

Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufuli Zanfelice, 329 PB02
Galpão - Lapa
São Paulo - SP - Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA

Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strzi 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK

Arjo A/S
Vassingerødevej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-5525 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA

Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/Marie Curie 5
Edificio Alfa Planta 6 oficina 6.1.-62
ES-28521 Rivas Vacía, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com

FRANCE
Arjo SAS
2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG

Arjo Hong Kong Limited
Room 411-414, 4/F, Manhattan Centre,
8 Kwai Cheong Road, Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA

Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST

Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building,North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

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At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges.

With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden

www.arjo.com

04.ST.00-INT1

arjo

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Arjo Inc.