INSTRUCTIONS FOR USE

AtmosAir Stretcher

Mattress Replacement System

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Bedienungsanleitung · Gebruiksaanwijzing · Mode d'emploi · Istruzioni per l'uso
Instrucciones de uso · Brugsvejledning · Bruksanvisning · Instruções de Utilização
Kullanım Talimatları · Οδηγίες χρήσης · Käyttöohjeet · Bruksanvisning

WARNING

To avoid injury, always read this *Instructions for Use* and accompanied documents before using the product.

Mandatory to read the *Instructions for Use*.

® and ™ are trademarks belonging to the Arjo group of companies
© Arjo 2019. SAT™ (Self Adjusting Technology™) is a trademark of MPTI, Ltd.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

ARJO HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJO PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY ARJO SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL ARJO BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJO TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties except as set forth in the written limited warranty included in this publication or with this product. Information in this publication may be subject to change at any time. Contact Arjo for updates.

Important Information For Users

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.
- If applicable, ensure the electrical installation of the room complies with the appropriate national / local electrical wiring standards.

Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support systems. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

Intentionally left blank

Table Of Contents

Introduction.....	3
Indications	3
Contraindications.....	3
Risks and Precautions.....	3
Safety Information	4
Serious Incident	5
Preparation for Use	6
Mattress Installation.....	6
<i>AtmosAir Stretcher MRS</i>	<i>6</i>
Patient Placement and Nursing Care	7
CPR	7
Skin Care.....	7
Incontinence / Drainage	7
General Operation.....	7
Care and Cleaning	8
Fire Barrier	8
Cover Cleaning Options	8
Preventive Maintenance Schedule.....	9
Daily Cleaning	9
Inspection / System Check-Out	9
Troubleshooting	10
Parts Diagram - <i>AtmosAir Stretcher MRS</i>	11
Replacement Parts	12
Specifications	13
Symbols Used.....	14
Customer Contact Information.....	14

Intentionally left blank

Introduction

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications, Safety Information and Risks and Precautions** sections prior to placing a patient on any AtmosAir™ Stretcher Mattress Replacement System (MRS).

Caregivers should review this information with the patient and the patient's family and / or legal guardian. Save this User Guide in an easily accessible location for quick reference.

Indications

- The AtmosAir Stretcher MRS is intended to treat and prevent all categories of skin breakdown.

Contraindications

- unstable vertebral fracture
- cervical and skeletal traction

Risks and Precautions

Transfer – Standard precautions should be taken during patient transfer. Refer to **Patient Placement and Nursing Care** section of these instructions

Side Rails and Restraints – WARNING: Use or non-use of restraints, including side rails, can be critical to patient safety. Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. **See related Safety Information.**

Patient Migration – Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent gurney exit. **Monitor patients frequently to guard against patient entrapment.**

Safety Information

Patient Entrance / Exit – Caregiver should always aid patient in exiting the gurney. Make sure a capable patient knows how to get out of gurney safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Brakes – Caster brakes should always be locked once the gurney is in position. Verify wheels are locked before any patient transfer to or from the gurney.

Gurney Height – To minimize risk of falls or injury, the gurney should always be in the lowest practical position when the patient is unattended.

Gurney Frame – Always use a standard gurney with this stretcher mattress, with any safeguards or protocols that may be appropriate. Gurney and side rails (if used) must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.

CPR - Level the gurney. Lower side rails and initiate CPR per facility protocols. Consider use of backboard if indicated. After CPR remove backboard, if used, raise siderails and reconfigure gurney and accessories as in initial placement.

Head of Gurney Elevation – Keep head of gurney as low as possible to help prevent patient migration.

Side Rails / Patient Restraints - Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from gurney) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. This includes assessment of the gurney occupant, and the combination of gurney frame, side rail and stretcher mattress (or mattresses where overlays are used). Risk assessment should be repeated if the gurney frame, stretcher mattress, side rail or condition of patient changes. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling off gurney and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, a description of at-risk patients and guidance to further entrapment risks, refer to the FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get off gurney safely (and, if necessary, how to release the side rails) in case of fire or other emergency. **Monitor patients frequently to guard against patient entrapment.**

CAUTION: When selecting a stretcher mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66" (220 mm) to help prevent inadvertent trolley exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

Skin Care – Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

No Smoking in Bed – Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

General Protocols – Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

Fire Barrier Contact - Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the stretcher mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. If the fire barrier becomes damaged or soiled it should be replaced.

End of Life Disposal –

- Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.
- Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.
- Pump units have electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

Serious Incident

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

Preparation for Use

For information concerning the gurney frame refer to the manufacturers' User Guide.

1. Open shipping container(s).

Do not use sharp instruments to open boxes. Damage to mattress could result.

2. Remove AtmosAir Stretcher MRS from plastic protective cover.

The mattress cover may appear wrinkled when unpacked. To remove wrinkles, allow mattress up to 24 hours to accommodate; see Troubleshooting for more information. Wrinkles will not affect inflation or function, so mattress may be used immediately if needed. Apply weight to mattress to open valves.

3. Check stretcher mattress surface for tears or cracking; do not use if tears or cracks are present.
4. If re-installing stretcher mattress onto a new frame or for a new patient, check stretcher mattress surface for staining and soiling; clean and / or disinfect as required (see Care and Cleaning).
5. Level gurney and lock brakes.
6. Remove existing mattress from gurney frame.
7. Remove any hook-and-loop material from base frame of patient gurney and / or clean any remaining residue, adhesive, etc. from base frame surface.

Mattress Installation

AtmosAir Stretcher MRS

1. Remove clear backing from hook-and-loop strips located on the underside of the AtmosAir Stretcher MRS.
2. Position stretcher mattress on gurney with printed text facing up, hook-and-loop strips facing down (against gurney base) and product information label at foot end of gurney.

Use caution when positioning stretcher mattress. Once the sticky side of the hook-and-loop strips adheres to the gurney base, the strips cannot be pulled up and repositioned.

3. Press stretcher mattress down to adhere sticky side of hook-and-loop strips to gurney base.

Always use a standard gurney frame with safeguards or protocols that may be appropriate. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body.

Patient Placement and Nursing Care

It is recommended that all sections of this User Guide be read prior to product use. Carefully review the **Contraindications, Safety Information and Risks and Precautions** sections prior to placing a patient on any *AtmosAir Stretcher MRS*.

1. Inspect, clean and disinfect stretcher mattress before each new patient is placed on the product.
2. Verify all caster brakes are locked for transfer.
3. Transfer patient following all applicable safety rules and institution protocols.
4. Center patient side-to-side and head-to-foot on stretcher mattress surface.
5. Ensure all sections of the mattress fully support the patient.

Mattress handles are for mattress transport only.

CPR

1. Level gurney.
2. Lower or remove side rail on caregiver's side if necessary.
3. Begin CPR as per facility protocol. Consider use of backboard if indicated.
4. After CPR is performed:
 - Remove backboard if used.
 - Raise or install side rail as necessary.
 - Reconfigure gurney and accessories as in initial placement.

Skin Care

- Remove excess moisture and keep skin dry and clean.
- Check patient's skin regularly, particularly in areas where incontinence and drainage occur.
- Ensure linens under patient are not wrinkled.

Incontinence / Drainage

- Use moisture-impermeable underpads for incontinent patients.
- Wipe surface clean and replace bed linens as required (see **Care and Cleaning**).

General Operation

Avoid contact of sharp instruments with the *AtmosAir Stretcher MRS*. Punctures, cuts and tears may prevent proper inflation and air pressure maintenance.

Care and Cleaning

The following processes is recommended, but should be adapted to comply with local institution protocols. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The AtmosAir Stretcher MRS should be routinely decontaminated between patients and at regular intervals while in use.

Do not use Phenol-based solutions or abrasive compounds or pads on the cover during decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.

Covers should be wiped down immediately after it has been exposed to fluids or spills.

Fire Barrier

Care should be taken with regards to the fire barrier when the cover is removed. It is recommended that Personal Protective Equipment (PPE) is used during handling of the stretcher mattress with the fire barrier exposed to avoid contact with loose particles. To ensure optimal product performance, avoid over-handling of the fire barrier and handle with care. Fire barrier should be replaced if soiled or damaged.

Cover Cleaning Options

Sewn AtmosAir Stretcher MRS (Top cover can not be removed from base)

1. Remove or push bed linens to center of mattress.

The sewn AtmosAir Stretcher MRS cover is cleaned by wipe down methods only. Do not launder as damage to the base may occur.

2. Wipe and rinse any soiling from the mattress surface and base.
Use 1000 ppm chlorine or 70% alcohol.
3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry surface with towel.
5. Ensure bed linens are refitted and not wrinkled under patient.
6. Clean pump and tubes (if needed) by wiping with damp cloth.

RF Welded AtmosAir Stretcher MRS (Detachable top cover)

1. Push bed linens to center of mattress to wipe down, remove to launder top cover or wipe down as described below.

The RF Welded AtmosAir Stretcher MRS has a detachable top cover that can be laundered. However, the base is cleaned by wipe down methods only. Do not launder base as damage may occur.

2. Wipe and rinse any soiling from the mattress surface and base.
Use 1000 ppm chlorine or 70% alcohol.
3. Rinse with a wet towel to remove chemicals.
4. After wipe down, dry surface with towel.
5. Ensure bed linens are refitted and not wrinkled under patient.
6. Clean pump and tubes (if needed) by wiping with a damp cloth.

Laundering Detachable Top Cover

1. Unzip the top cover from the base for laundering.
Do not launder base as damage may occur.
2. Recommended wash temperature for top cover is 60° C (140° F) for 15 minutes.
3. Maximum wash temperature is 95° C (203° F) for 15 minutes.
4. Tumble dry at 60° C (140° F) or air dry.
5. Maximum drying temperature 80° C (176° F).

Preventive Maintenance Schedule

Preventive maintenance for the *AtmosAir Stretcher MRS* consists of regular cleaning (see **Care and Cleaning**) and an overall system check-out to be performed at the intervals described below.

All components must be cleaned, disinfected and inspected after each patient's use and before use by a new patient. Always use standard precautions, treating all used equipment as contaminated. Institutions should follow local protocols for cleaning and disinfection.

Daily Cleaning

The cover should be wiped daily with a mild soap and water solution. After wipe down, dry surface with a towel.

Inspection / System Check-Out

Check each of the following before placing the *AtmosAir Stretcher MRS* with a new patient:

1. Check stretcher mattress surface for tears or cracking; do not use if tears or cracks are present.
2. Ensure stretcher mattress is free of stains and is not overly faded.

Troubleshooting

It is recommended that all sections of this manual be reviewed before troubleshooting any *AtmosAir* Stretcher MRS.

Do not attempt troubleshooting outside this guide or where the solution recommends to contact Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

SYMPTOM	POSSIBLE CAUSE	SOLUTION
Mattress too firm upon arrival.	Difference in altitude not sufficient to open valves.	Apply weight to mattress to open valves.
Mattress cover too wrinkled upon removal from shipping container.	Internal components have not accommodated to environment. This does not affect inflation or function.	Let mattress accommodate for 24 hours. If problem continues, contact Arjo for assistance.
Mattress does not inflate or is not firm	Tubing not connected properly. Tubing kinked. Tubing disconnected. Holes in or damage to SAT™ (Self Adjusting Technology).	Check tubing inside mattress for loose connectors. Check tubing inside mattress for kinks. Check tubing inside mattress for possible disconnect. Check SAT for holes or damage or contact Arjo for assistance.


Parts Diagram - AtmosAir Stretcher MRS

All hoses were removed to improve diagram readability.

The optional fire barrier is a sleeve that fits over the foam.

The SAT Cover, SAT Cells and SAT Interliner are all part of the integrated SAT System and cannot be ordered separately. See Replacement Parts for a complete list of SAT Systems.

Replacement Parts

Do not attempt troubleshooting, maintenance or parts replacement outside this manual or where the solution recommends contacting Arjo. Any unauthorized service, modification, alteration or misuse may lead to serious injury and / or product damage and will void all applicable warranties.

Replaceable *AtmosAir* components are listed below. For more information such as pricing or additional spare parts that are not on this list, please contact your local Arjo representative.

ATMOSAIR MATTRESS COMPONENTS

AtmosAir Foam Shells

PART TYPE.....	CATALOG NUMBER
Stretcher	KSFS(SIZE)

AtmosAir with SAT System

PART TYPE.....	CATALOG NUMBER
Stretcher SAT System (Fits mattresses 25-30 in wide, 75-77 in length)	KA9SAT

AtmosAir Fire Barrier Lining

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> Fire Barrier Lining	KA4/9-FBL

AtmosAir Black Pressure Sensitive Loop

PART TYPE.....	CATALOG NUMBER
Stretcher 29" Black Pressure Sensitve Loop.....	KS29PSLOOP

ATMOSAIR MATTRESS COVERS*

AtmosAir Model Stretcher

PART TYPE.....	CATALOG NUMBER
Complete Cover, Reliant IS ² ™ Top, Non-Skid Bottom	
No Handles.....	ASCRLNH(WWWLL)TXS
Complete Cover, <i>Reliant</i> IS ² , RF Welded Top, Non-Skid PU Bottom,	
No Handles.....	ASCRFPN(WWWLL)TXS
Complete Cover, Premium, RF Welded Top, Non-Skid PU Bottom,	
No Handles.....	ASCPRPN(WWWLL)NXS
Top Cover, <i>Reliant</i> IS ² , RF Welded	ASCRFXX(WWWLL)TXS
Top Cover, Premium, RF Welded	ASCPRXX(WWWLL)NXS
Base Cover, PU Non-skid Bottom, No Handles.....	ASCXXPN(WWWLL)XXS

*WWW = width in 0.5 in increments, LL = Length in 1.0 in increments

ATMOSAIR LITTERATURE

PART TYPE.....	CATALOG NUMBER
<i>AtmosAir</i> MRS Instructions for Use	407308-AH

Specifications

Specifications subject to change without notice.

Maximum Recommended Patient Weight*	500 lbs (227 kg)
Weight (based on 35" x 80" x 7" MRS, other sizes will vary)	30 lbs (14 kg)
Dimensions (varies depending on model)**	
Minimum Length.....	75 in (190.5 cm)
Maximum Length.....	77 in (195.6 cm)
Minimum Width.....	25 in (63.5 cm)
Maximum Width.....	30 in (76.2 cm)
Height	5 in (12.7 cm)

* Patient weight capacity may vary depending on frame use. Please consult frame manufacturer.

** Additional sizes are available. Contact Arjo for more information.

Symbols Used

Important Operational Information

Foot End

Warning of possible hazard
to system, patient or staff

Consult User Guide

Manufacturer

= Safe Working Load

Wash at 60°C, 95°C
Max for 15 minutes

Serial Number

The operator must read this
document (this instruction for use)
before use. Note: the symbol is blue
on the product label.

Do Not Iron

Use Solution diluted to 1000 ppm
of available chlorine

Attention - See User Guide

CE marking indicating conformity
with European Community
harmonised legislation

Wipe down Only

Tumble Dry at 60°C, 80°C Max

Do not use Phenol-based
cleaning solutions

Indicates the product is a Medical
Device according to EU Medical
Device Regulation 2017/745

Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative, or visit www.arjo.com.

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At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden

www.arjo.com

arjo

CE