Instructions for Use • Brugsvejledning
Instrucciones de uso • Mode d'emploi
Istruzioni per l'uso • Gebruiksaanwijzing • Bruksanvisning

English

Dansk

Español

Français

Italiano

Nederlands

Svensk

Contents

General Safety	iii
Introduction	1
About this Manual	1
About Alpha Active 4	1
Alpha Active 4 Pump	1
Alpha Active 4 Mattress Replacement	2
Clinical Applications	3
Indications	3
Contra-Indications	3
Installation	4
Preparing the Alpha Active 4 system for use	4
Installing the Alpha Active 4 Mattress	4
Controls, Alarms and Indicators	6
Controls	6
Alarms and Indicators	7
Operation	8
Quick Start	8
Support Setting Procedure	8
Changes in Patient Position	9
To Deflate and Store the Alpha Active 4 Mattress	10
CPR Control	11
To activate CPR	11
To reset CPR	11
Decontamination	12
Routine Maintenance	14
Alpha Active 4 System	14
Alpha Active 4 Pump	14
Alpha Active 4 Mattress Replacement	14
Serial Labels	14
Troubleshooting	15
Technical Description	16
Pump	16
Mattress	17
Mattress Size Information	17
Cover Specification	17
Cleaning Symbols	18

GENERAL SAFETY

Before you connect the system pump to a mains socket, read carefully all the installation instructions contained within this manual.

The system has been designed to comply with regulatory safety standards including:

- EN60601-1:1990/A13:1996 and IEC 60601-1:1988/A2:1995
- UL60601-1 and CAN/CSA C22.2 No. 601.1-M90

Safety Warnings

- It is the responsibility of the care giver to ensure that the user can use this product safely.
- Whilst the patient is unattended, safety sides should be used based on clinical assessment and in line with local policy.
- Alignment of the bed frame, safety sides and the mattress should leave no gap wide enough to entrap a patient's head or body, or to allow egress to occur in a hazardous manner where entanglement with the mains power cable and tubeset or air hoses may result. Care should be exercised to prevent occurrence of gaps by compression or movement of the mattress. Death or serious injury may occur.
- Make sure that the mains power cable and tubeset or air hoses are positioned to avoid causing a trip or other hazard, and are clear of moving bed mechanisms or other possible entrapment areas. Where cable management flaps are provided along the sides of the mattress, these should be used to cover the mains power cable.
- Electrical equipment may be hazardous if misused. There are no user-serviceable parts inside the pump. The pump's case must only be removed by authorised technical personnel. No modification of this equipment is allowed.
- The mains power socket/plug must be accessible at all times. To disconnect the pump completely from the electricity supply, remove the plug from the mains power socket.
- The CPR control and/or the CPR indicator tag must be visible and accessible at all times.
- Disconnect the pump from the mains power socket before cleaning and inspecting.
- Keep the pump away from sources of liquids and do not immerse in water.
- Do not use the pump in the presence of uncontained flammable liquids or gasses.
- The cover of this product is vapour permeable but not air permeable and may present a suffocation risk.
- Only the pump and mattress combination as indicated by ArjoHuntleigh should be used. The correct function of the product cannot be guaranteed if incorrect pump and mattress combinations are used.

Precautions

For your own safety and the safety of the equipment, always take the following precautions:

- Placing extra layers between the patient and the mattress potentially reduces the benefits provided by the mattress and should be avoided or kept to a minimum. As part of sensible pressure area care, it is advisable to avoid wearing clothing which may cause areas of localised high pressure due to creases, seams, etc. Placing objects in pockets should be avoided for the same reason.
- Do not expose the system, especially the mattress, to naked flames, such as cigarettes, etc.
- Do not store the system in direct sunlight.
- Do not use phenol-based solutions to clean the system.
- Make sure the system is clean and dry prior to use or storage.
- Never use sharp objects or electrically heated under blankets on or under the system.
- Store the pump and mattress in the protective bags supplied.

Electromagnetic Compatibility (EMC)

This product complies with the requirements of applicable EMC Standards. Medical electrical equipment needs special precautions regarding EMC and needs to be installed in accordance with the following instructions:

- The use of accessories not specified by the manufacturer may result in increased emissions by, or decreased immunity of, the equipment, affecting its performance.
- Portable and mobile radio frequency (RF) communications equipment (e.g. mobile/cell phones) can affect medical electrical equipment.
- If this equipment needs to be used adjacent to other electrical equipment, normal operation must be checked before use.
- For detailed EMC information contact ArjoHuntleigh service personnel.

Environmental Protection

Incorrect disposal of this equipment and its component parts, particularly batteries or other electrical components, may produce substances that are hazardous to the environment. To minimise these hazards, contact ArjoHuntleigh for information on correct disposal.

Service Information

ArjoHuntleigh recommend that this system should be serviced every 12 calendar months or, where applicable, when the service indicator is illuminated.

Design Policy and Copyright

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies. As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

© ArjoHuntleigh 2009.

1. Introduction

About this Manual This manual is your introduction to the **Alpha Active 4**. Use it to initially set up the mattress and keep it as a reference for day-to-day routines and as a guide to maintenance.

About Alpha Active 4 The **Alpha Active 4** systems comprise of a mattress replacement and pump. The support system can be used on hospital and domestic beds.

Alpha Active 4 Pump The **Alpha Active 4** pump comprises of a moulded case with non slip feet on the base and integral hanging brackets.

The controls are situated on the front of the pump and a sophisticated alarm system differentiates between normal operation and genuine system faults. If an alarm situation is detected an indicator will illuminate on the front of the pump and an audible warning will sound.

Alpha Active 4 Mattress Replacement

The **Alpha Active 4** mattress replacement comprises the following components:

Detachable Cover The standard cover comprises of a 2-way stretch PU (polyurethane) coated knitted fabric zipped to a durable nylon base. The zips are protected by flaps to prevent ingress of contaminants and allow easy removal of the cover for cleaning.


Cells The mattress comprises of 19 PU cells, 16 providing support to the user in either Alternating (Active) or Static (Reactive) mode and 3 Static head cells.

CPR function A CPR (Cardio-Pulmonary Resuscitation) control is positioned at the head end of the mattress to allow rapid deflation of the mattress.

Tube-set The tube-set has a 3-way pneumatic connection which incorporates a flexible, compact anti-kink tube that is resistant to crushing and any subsequent obstruction of air flow.

When disconnecting the tube-set, place the attached cover over the end to place the mattress in transport mode.

Mattress Replacement Base Cover The base cover for the mattress replacement is PVC (Polyvinyl Chloride) coated nylon on the underside. There are six straps located on the underside of the mattress replacement to secure the mattress to the bedframe.

2. Clinical Applications

Indications The **Alpha Active 4** system is indicated for the prevention and management of pressure ulcers when combined with an individualised monitoring, repositioning and wound care programme.

The **Alpha Active 4** mattress is designed for patients weighting up to 200 kg (440lb).

Contra-Indications The **Alpha Active 4** system should not be used for patients with unstable spinal fractures.

In the case of patients with other unstable fractures, where a moving surface can be harmful, advice should be obtained from the appropriate physician before using the **Alpha Active 4** system.

The above are guidelines only and should not replace clinical judgement or experience.

3. Installation

Preparing the Alpha Active 4 system for use

The system is very simple to set up using the following guidelines:

1. Remove the system from the packaging. You should have the following items:
 - **Alpha Active 4** pump including mains power cord and hanging brackets.
 - **Alpha Active 4** mattress replacement with integral tube-set.
 - Cover.

Installing the Alpha Active 4 Mattress

The **Alpha Active 4** mattress replacement system should be installed as follows:

1. Remove the existing mattress from the bed frame and check that there are no protruding bed springs or sharp objects on the bed frame surface.
2. Unroll the mattress onto the bed frame and ensure that the tube-set is located near the foot end of the bed and the CPR at the head end. The cells of the mattress must be uppermost.
3. Secure the mattress to the bed frame using the 6 fastener straps.

 *If the bed can be profiled to any position (i.e. raised or lowered), attach the mattress to the movable parts of the bed only.*

To Complete the Mattress Installation

Complete the installation of the mattress replacement as follows:

1. If not already fitted, place the protective cover over the mattress. Ensure that the logo is uppermost and at the foot end of the mattress.
2. Zip the cover onto the mattress starting from the head end and taking care not to trap any material in the zip.
3. Ensure that the CPR unit is secured in its closed position.

 *The CPR must be accessible at all times.*

WARNING

Make sure that the mains power cable is positioned to avoid causing a hazard and is clear of moving bed mechanisms or other possible entrapment areas.


Installing the Pump The pump should be installed as follows:

1. Position the pump, feet down, on any convenient horizontal surface or alternatively suspend from the bed foot rail by means of the integral hanging brackets.
2. Ensure that the mattress tube-set is not "kinked" or twisted and connect it to the pump until it clicks into place. Ensure that the tube-set is securely connected to the pump.
3. Insert the mains power plug into a suitable mains power socket.

System Operation The system is now ready for use. See "Operation" on page 8 for day-to-day operating instructions.

4. Controls, Alarms and Indicators

Controls

On/Off Button

Pressing the **On/Off** button will activate the pump. The switch will illuminate when the pump is on.

Alarm Mute

During an alarm condition the sound of the alarm can be muted by pressing this button. The yellow indicator will then remain on but the alarm will be muted.

Mode Select

The mode selector allows the mode of operation to be selected.

Two modes of operation are available, Alternating or Static (non-alternating). The indicator next to the respective icon indicates the currently selected mode.

Auto-firm Mode

Selects the **Auto-firm**, non-alternating mode. The orange indicator will illuminate when the pump is in this mode.

Auto-firm lasts 30 minutes, after which the pump will default back to the previous setting.

Patient Weight Selector

Two buttons (+ and -) are used to select the patient weight. This sets the pressure inside the mattress replacement and therefore the amount of support the patient receives.

The pressure setting is indicated by the green indicator above the selected setting.

Cycle Time

A single button allows the cycle time to be selected for patient and/or clinician preference.

The cycle time setting is indicated by the green indicator above the selected setting.

Alarms and Indicators

Low Pressure Indicator

The red **Low Pressure** indicator is illuminated whenever the pump detects low pressure within the mattress replacement. An audible alarm will sound unless cancelled by the **mute** button.

The indicator will extinguish once normal pressure is reached.

 See “Troubleshooting” on page 15 for possible causes of Low Pressure.

Service Indicator/Pump Fault

The yellow **Service/pump fault** indicator will illuminate and remain on if the pump has detected an internal fault. A Service Engineer should be called.

Power Fail Indicator

The red **Power Fail** indicator will illuminate when a mains power failure has been detected. An audible alarm will sound until power is resumed or the pump is switched off using the on/off button.

5. Operation

These instructions cover day-to-day operation of the system. Other operations, such as maintenance and repair, should only be carried out by suitably qualified personnel.

Quick Start Before using the **Alpha Active 4** mattress replacement make sure it has been installed correctly in accordance with “Installing the Alpha Active 4 Mattress” on page 4, and ensure that the CPR unit on the mattress is clicked into the closed position.

1. Press the On/Off button on the side of the pump.
2. Allow approximately 30 minutes for the mattress replacement to inflate fully.
3. Place a bed sheet loosely over the mattress without tucking, to gain maximum pressure redistribution.

Support Setting Procedure It is important to follow the correct support setting procedure to ensure the patient receives adequate support (lift) while achieving maximum pressure redistribution and comfort.

1. Lie or sit the patient on the mattress.
2. Select appropriate cycle time.
3. Select patient weight on pump. This should serve as an approximate guide only. An independent clinical determination needs to confirm that the patient is properly supported.
4. Wait 10/20 minutes while the pump adjusts the pressures.
5. Ensure that the patient is not ‘bottoming out’. To do this, unfasten the cover and slide a hand beneath the patients sacral area to ensure at least 2.5cm clearance between the sacrum and bed base.
6. If the caregiver feels less than 2.5cm of support material, the patient has bottomed out and the support pressures should be adjusted accordingly.

Changes in Patient Position

When a patient is in the lying or supine position, their body weight is dispersed over a large area. When in the sitting position, their body weight is concentrated within a much smaller area and therefore will require more support than in the lying position.

Therefore, when the patient changes position, it may be necessary, in order to maximise the benefit of the support surface, to make adjustments to the setting of the patient weight selector.

From Lying to Sitting - Increase weight selector setting.

From Sitting to Lying - Decrease weight selector setting.

This adjustment should be in conjunction with independent clinical determination of appropriate support.

Static Provides a stable, non-moving support surface for instances where active therapy is not indicated e.g. to carry out nursing procedures or for patients unable to tolerate a moving surface. In **Static** mode the support surface remains constant (all cells are equally inflated). Additional nursing assessment must be undertaken in order to direct an individualised repositioning programme.

When operating the system in **Static** mode it may be necessary, where possible, to reduce the pressure setting to increase patient comfort and safety.

Power Fail Condition If a **Power Fail** condition arises disconnect the tubeset from the pump and place the attached cap over the end of the tubeset to put the mattress into transport mode. Once power is resumed, re-connect the tubeset to the pump to continue therapy.

To Disconnect the Tube-set

To disconnect the tube-set at any time, depress the buttons on the top and bottom of the tubeset connector pull the tube-set connector away from the pump.

To deflate the mattress Refer to “To Deflate and Store the Alpha Active 4 Mattress” on page 10.

Transport Mode To transport a patient using the **Alpha Active 4** mattress replacement, disconnect the tube-set from the pump and place the attached cap over the end of the tubeset to put the mattress into transport mode. This will automatically switch the mattress into transport mode.

The patient will remain supported by the mattress for up to 12 hours.

To resume normal operation, simply reconnect the tube-set and run the pump.

Caution

Transport mode is non-therapeutic offering support only for up to 12 hours. It is recommended that when in transport mode the patient is frequently monitored.

To Deflate and Store the Alpha Active 4 Mattress

- To deflate the mattress:**
1. Disconnect the tube-set from the pump.
 2. Activate the CPR control to deflate the mattress.

- To store the mattress** Following deflation:
1. Bring the tubeset over the mattress to lie parallel to the foot end of the mattress.
 2. Roll the mattress from the foot end toward the head end of the mattress.


CPR Control

IMPORTANT

IN THE EVENT OF CARDIAC ARREST

In the event of a patient suffering cardiac arrest and CPR needing to be administered:

To activate CPR Located at the head end of the mattress replacement (on the same side as the tubeset) is a red strap labelled CPR. In the event of a cardiac arrest pull this from the mattress to deflate.

To reset CPR To re-inflate the mattress, simply replace the stopper securely into the manifold.

6. Decontamination

The following processes are recommended, but should be adapted to comply with the local or national guidelines (Decontamination of Medical Devices) which may apply within the Healthcare Facility or the country of use. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The **Alpha Active 4** system should be routinely decontaminated between patients and at regular intervals while in use; as is good practice for all reusable medical devices.

WARNING

Remove the electrical supply to the pump by disconnecting the mains power cord from the mains power supply before cleaning. Protective clothing should always be worn when carrying out decontamination procedures.

Caution

Do not use Phenol-based solutions or abrasive compounds or pads during the decontamination process as these will damage the surface coating. Do not boil or autoclave the cover. Avoid immersing electrical parts in water during the cleaning process. Do not spray cleaning solutions directly onto the pump.

To clean Clean all exposed surfaces and remove any organic debris by wiping with a cloth moistened with a simple (neutral) detergent and water. Dry thoroughly.

Chemical Disinfection To protect the integrity of the cover we recommend a chlorine-releasing agent, such as sodium hypochlorite, at a strength of 1,000ppm available chlorine (this may vary from 250ppm to 10,000ppm depending on local policy and contamination status).

Wipe all cleaned surfaces with the solution, rinse and dry thoroughly.

Alcohol based disinfectants (maximum strength 70%) may be used as an alternative.

Ensure the product is dry before storage.

If an alternative disinfectant is selected from the wide variety available we recommend that suitability for use is confirmed with the chemical supplier prior to use.

DO NOT WRING/MANGLE, AUTOCLAVE OR USE
PHENOLIC BASED SOLUTIONS.

Thermal Disinfection For information for the mattress top cover, including
laundering guidelines, refer to “Cover Specification” on
page 17.

7. Routine Maintenance

Alpha Active 4 System

- Maintenance** The equipment has been designed to be virtually maintenance-free between service periods.
- Servicing** **ArjoHuntleigh** will make available on request service manuals, component parts lists and other information necessary for **ArjoHuntleigh** trained personnel to repair the system.

Alpha Active 4 Pump

General Care, Maintenance and Inspection

Check all electrical connections and power cord for signs of excessive wear.

In the event of the pump being subjected to abnormal treatment, e.g. immersed in water or dropped, the unit must be returned to an authorised service centre.

Alpha Active 4 Mattress Replacement

- General Care** Remove the top cover and inspect for signs of wear or any tears.
- Check all zips are secure.
- Check integrity of all connectors, including cell to manifold connections.
- Ensure all cell fasteners are correctly connected to the mattress base sheet and are not loose or damaged.

Serial Labels

- Pump** The serial number for the pump is on the label on the back of the pump case.
- Mattress** The mattress serial label can be found just inside the base cover above the tubeset. Refer to “Alpha Active 4 Mattress Replacement” on page 2.

8. Troubleshooting

The following table provides a troubleshooting guide for the **Alpha Active 4** system in the event of malfunction.

Indicator	Possible Cause	Remedy
LOW PRESSURE	The tube-set is not connected properly.	Check the tube-set connector and ensure it is securely fitted to the pump.
	CPR not fully closed.	Close CPR unit.
	There is a leak in the system.	Call service engineer.
POWER FAILURE	Power has been removed from the pump.	Re-apply power or switch the pump off.
SERVICE	Pump has detected an internal fault.	Switch the pump off and call service engineer.

9. Technical Description

PUMP							
Model:	Alpha Active 4						
Part Numbers:	648311	UK					
	648319AU	AUSTRALIA					
Supply Voltage:	230V						
Supply Frequency:	50-60Hz						
Power Input:	230VA						
Size:	(L)280mm x (W)205mm x (H)112mm						
Weight:	2.5kg						
Case Material:	ABS Plastic						
Plug Fuse Rating:	5A to BS1362 (UK ONLY)						
Pump Fuse Rating:	2 x T1AL 250V						
Degree of protection against electric shock:	Class II Type BF						
Degree of protection against liquid ingress:	IPX0						
Mode of operation:	Continuous						
Cycle Times:	10 mins	20 mins					
	Inflate - 4.5 mins	Inflate - 9.5 mins					
	Crossover - 30 secs	Crossover - 30 secs					
	Deflate - 4.5 mins	Deflate - 9.5 mins					
	Crossover - 30 secs	Crossover - 30 secs					
SYMBOLS							
	Power Disconnects from the mains supply		Do not dispose of in the domestic refuse		Alternating Current		Double Insulated
	Power Connects to the mains supply		Refer to accompanying documents		Type BF		Refer to the User Manual
	Dangerous voltage		Fuse (If shown on pump label)	Ref:	Model number	SN:	Serial Number

ENVIRONMENTAL INFORMATION			
Condition	Temperature Range	Relative Humidity	Atmospheric Pressure
<i>Operating (Pump)</i>	+10°C to +40°C	30% to 75%	700hPa to 1060 hPa
<i>Storage and Transport (Pump)</i>	-40°C to +70°C	10% to 95% (non-condensing)	500 hPa to 1060 hPa

MATTRESS					
Description		Cell Material	Base Pad Material		
ALPHA ACTIVE 4 MR 90		Polyurethane	PVC Nylon Coated		
ALPHA ACTIVE 4 MR 85		Polyurethane	PVC Nylon Coated		
ALPHA ACTIVE 4 MR 85 (NON-PVC)		Polyurethane	Nylon PU Coated		
MATTRESS SIZE INFORMATION					
Part No.	Description	Spare Cover	Length mm	Width mm	Height mm
648322	ALPHA ACTIVE 4 MR 90	648437	1950(77")	900(35")	200 (8")
648324	ALPHA ACTIVE 4 MR 85	648461	1950(77")	850(33")	200 (8")
648326	ALPHA ACTIVE 4 MR 85 (NON-PVC)	648461	1950(77")	850(33")	200 (8")

COVER SPECIFICATION	
Feature	Standard Cover (Dartex) [®]
<i>Removable Cover</i>	Yes
<i>Moisture Vapour Permeable</i>	Yes
<i>Air Permeable</i>	No
<i>Low Friction</i>	Yes
<i>Water Resistant / Repellent</i>	Yes
<i>Infection Control</i>	Material coating is Bacteriostatic, fungistatic, antimicrobial
<i>Fire Retardant</i>	BS 7175: 0,1 & 5
<i>2-Way Stretch</i>	Yes
<i>Washing Conditions</i>	MAX 95°C (203°F) for 15 mins ¹
<i>Drying Conditions</i>	Tumble Dry up to 130°C (266°F) or Air Dry
<i>Life Span</i>	50 Wash Cycles (minimum)
<i>Application Area</i>	Acute and Homecare

1. Check your local policy to determine the time/temperature ratio required to achieve thermal disinfection.

CLEANING SYMBOLS			
	Wash at 80°C (176°F)		Do not tumble dry above 130°C.
	Do not iron		Do Not Use Phenol-based cleaning Solutions
	Wipe surface with damp cloth		Use solution diluted to 1000 ppm of Available Chlorine

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ArjoHuntleigh is a branch of Arjo Ltd Med AB, Arjo Ltd Med AB Registered No: FC024121