Bari-Breeze

Instructions for Use

ARJOHUNTLEIGH
GETINGE GROUP

Bedienungsanleitung • Instrucciones de uso • Notice d'utilisation
Istruzioni per l'uso • Gebruiksaanwijzing

EN

DE

ES

FR

IT

NL

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product.

Mandatory to read the Instructions for Use

Design Policy and Copyright

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies.

© ArjoHuntleigh 2015.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of ArjoHuntleigh.

English

Deutsch

Español

Français

Italiano

Nederlands

Intentionally left blank

English

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product.

Mandatory to read the Instructions for Use

Design Policy and Copyright

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies.

© ArjoHuntleigh 2015.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of ArjoHuntleigh.

Contents

Introduction	1
About this Manual	1
About Bari-Breeze	1
Bari-Breeze Mattress	2
Clinical Applications	5
Indications	5
Contraindications	5
Precautions for Use	5
Installation and Setup	6
Warnings	6
Installation	6
CPR Facility	8
Controls, Alarms and Indicators	9
Controls and Indicators	9
Alarms	10
Operation	11
Initial Power Up	11
Comfort Controls	11
Static Mode	11
Turn Options	12
Turn Times	12
To Deflate and Store the Bari-Breeze Mattress	12
Decontamination	13
Troubleshooting Guide	15
Spare Parts	17
Warranty and Service	18
Technical Data	19

General Safety

Before you connect the system pump to a mains socket, read carefully all the installation instructions in Section 3 - Installation. The system has been designed to comply with regulatory safety standards including:

- EN60601-1:1990/A13:1996
- UL60601-1 1st Edition
- CAN/CSA C22.2 No 601.1-M90

Safety Warnings

- The cover of this product is vapour permeable but not air permeable and may present a suffocation risk. It is the responsibility of the care provider to ensure that the user can use this product safely.
- Electrical equipment may be hazardous if misused. The front and rear case of the pump unit should only be separated by authorised technical personnel.
- Do not use the pump in the presence of flammable gases such as anaesthetic agents.
- Make sure that the mains power cable is positioned to avoid causing a trip or other hazard, and is clear of moving bed mechanisms or other possible entrapment areas.

Precautions

For your own safety and the safety of the equipment, always take the following precautions:

- Keep the pump away from sources of liquids and do not immerse in water.
- Do not expose the system, especially the mattress, to naked flames, such as cigarettes etc.
- Do not store the system in direct sunlight.
- Switch off the electrical supply to the pump by disconnecting the pump from the mains socket before cleaning and inspection.
- Ensure the system is clean and dry prior to storage.
- Never use sharp objects or electrically heated under blankets on or under the system.
- Store the pump and mattress in the protective bags supplied.
- Only the pump and mattress combination as indicated by ArjoHuntleigh should be used. The correct function of the product cannot be guaranteed if the incorrect pump or mattress combinations are used.

Caution

Electromagnetic compatibility (EMC). This product complies with the requirements of applicable EMC Standards. The use of accessories not specified by the manufacturer may result in increased emissions by, or decreased immunity of, the equipment, affecting its performance.

1. Introduction

About this Manual This manual is your introduction to the **Bari-Breeze®** system. Use it to initially set up the low air loss system, keep it as a reference for day-to-day routines, and as a guide to maintenance.

About Bari-Breeze The **Bari-Breeze** is a constant low pressure support surface.

A permanently inflated sub-base and bolsters along each side of the mattress are designed to provide optimal support and comfort for patients particularly when the system is operating in **Turn** mode.


The **Bari-Breeze** pump is fully programmable and provides the following modes of operation:

- **Static Mode (Default)** - When **Turn** therapy settings are set to , the system remains flat, providing a continuous low pressure surface.
- **Turn (Rotation) Mode** - Provides an angled support surface to assist with nursing procedures.
- **Firm Mode (Maxflow)** - Provides a firm support surface (for 15 minutes) to facilitate nursing procedures and/or patient handling.

The **Bari-Breeze** cover is water resistant and vapour permeable to enhance patient comfort whilst protecting the cells from ingress of contaminants. The cover is simple to clean in situ, but may easily be removed for laundering, preventing cross contamination.

In the event of cardiac arrest the mattress replacement can be easily deflated to allow cardiac resuscitation procedure (CPR) to be performed.

Bari-Breeze Mattress

A self contained mattress replacement system with easily detachable components for cleaning.

Cells	Polyurethane cells approximately 250mm (10”) high providing support to the patient.
Tubeset	Manufactured from material to reduce accidental kinking.
Tubeset Connector/CPR	Located at the pump end of the tubeset is a specially designed twin snap lock connector. This design helps prevent accidental disconnection of the mattress from the pump. Releasing the connector allows rapid deflation of the mattress in order to facilitate CPR. For additional access or to facilitate transfer, the side bolsters and base air pad can be deflated by removing the air plug.
Cover	The mattress is totally enclosed in a double zipped cover. The top and sides of the cover are constructed from a two-way stretch and water resistant material, which is vapour permeable for maximum patient comfort. The base of the cover is constructed from a tough abrasion resistant material and is fitted with bed attachment straps.
Base Air Pad	A 50 mm (2”) thick sub-mattress made of a single air-filled cell enclosed in the base to support the patient in the event of loss of air pressure in the mattress.
Bolsters	Two side air bolsters which provide optimal support and comfort for the patient particularly when the system is operating in Turn mode.
Air Plug	The air plug is used to deflate both the side bolsters and base air pad. It is located close to the tubeset/pump connector.

Bari-Breeze Pump

Mounting Hooks Integrated mounting hooks fitted to the rear casing allow easy positioning of the pump from the bed foot rail.

Keypad The keypad is angled to provide the user with a clear visual display. An LED display and individually illuminated control buttons give precise indications of the status of the system at any point in time.

LED Display Highly visible and bright 7 segment blue LED digital display panel.

2. Clinical Applications

Indications The **Bari-Breeze** system is indicated for patients up to the highest risk of pressure ulceration and weighing up to 455Kgs (1000 lbs), and is suitable for the prevention and management of all grades of pressure ulcers.

Contraindications Patients with unstable spinal fractures should not be placed on the **Bari-Breeze** system.

Precautions for Use In the case of patients with unstable fractures advice should be obtained from the appropriate physician before using the **Bari-Breeze** system.

Care should be taken to ensure that all patient tubes and lines are suitably positioned. This will help minimise the risk of pull or disconnection whilst **Turn** mode is selected or in the event of a power failure.

Regular monitoring of patients is recommended, particularly whilst in **Turn** mode.

In any lateral rotation system there may be an increased risk of gradual patient migration into hazardous positions. This may result in entrapment and/or inadvertent bed exit.

 *A risk assessment tool combined with clinical judgement should be used when determining a patient's level of risk of developing pressure ulcers. Patient risk assessment should be an ongoing process as changes in the patient's condition may increase or decrease their risk level.*

 *The above are guidelines only and should not replace clinical judgement or experience.*

3. Installation and Setup

Warnings The following warnings should be taken note of before installing the **Bari-Breeze** system.

WARNING

The Bari-Breeze system must be installed on bed frames that are equipped with Safety Sides. Please raise safety sides on the bed and lock them in position after the patient is on the mattress.

NEVER LEAVE A PATIENT UNATTENDED ON THE MATTRESS SYSTEM WITH THE SAFETY SIDES IN THE DOWN POSITION.

WARNING

Never use sharp objects or electrically heated underblankets on or under the mattress.

WARNING

Make sure that the mains power cable is positioned to avoid causing a hazard and is clear of moving bed mechanisms or other possible entrapment areas.

Installation The following components should be found within your **Bari-Breeze** packaging:

- **Bari-Breeze** Mattress Replacement
 - Pump
 - User Manual
 - Power Cord
1. Before using the **Bari-Breeze** mattress replacement system, please remove all other mattress systems from the bed.
 2. When installing the **Bari-Breeze** mattress, care should be taken such that the mattress is placed directly on the bed frame.

Make sure that the tubeset end of the mattress is towards the foot end of the bed.

3. There are 16 nylon black straps with buckles, on the base of the mattress. Loop each strap around the bed frame and fasten it securely to the bed frame using the buckle.

If the bed can be profiled to any position (i.e. raised or lowered), attach the mattress to the movable parts of the bed only.

4. Open the hooks on the back of the pump and suspend the pump from the bed foot rail. If the bed you are using does not have a foot rail, place the pump on its base (not on its back where the filter is located) on a flat surface underneath the bed near the foot of the bed frame.

Care should be taken such that the air inlet vent on the pump is not covered, and the pump is not placed on the floor in such a manner that it is a hazard.

5. Plug the power cord into the pump and press it in place. Uncoil the power cord and plug the cord into a properly grounded AC power source.
6. Securely attach the pump connector onto the pump until it clicks, and make sure the side bolster and base air pad plug is secure.


Care should be taken such that the mattress hose is freely suspended without being pinched or kinked.

7. If not already fitted, place the protective cover over the mattress and secure in position. Make sure the **ArjoHuntleigh** logo is uppermost and at the foot end of the mattress.
8. Push the Power/Standby button to turn the pump on and enter the required therapy settings, see “Controls, Alarms and Indicators” on page 9.
9. Unzip the cover and ensure the cells, side bolsters, and the base pad are all fully inflated. Refit cover.
10. Place a bed sheet loosely over the mattress before use.

4. CPR Facility

IMPORTANT

In the event of cardiac arrest, depress the two release catches on the tubeset/pump connector, and at the same time pull the tubeset out and away from the pump.


Located on the right hand side of the pump is a CPR device marked with an arrow. In the event of cardiac arrest, press the two quick release catches on the tubeset/pump connector, and simultaneously pull the hose away from the pump.

- To reinflate the mattress**
1. Switch the pump into standby mode and reconnect the tubeset and base air pad plug if removed. Switch the pump on and reset the therapy settings.

5. Controls, Alarms and Indicators

Controls and Indicators

Standby

When the pump has power but is in Standby mode, the 7 segment LED display will show a single bar.

Power/Standby

Press the Power/Standby button to turn the pump on.

Max Flow


Press  to enable maximum air flow from the pump. A blue LED on the button will indicate that **Max Flow** is active.

If **Max Flow** 'F' is selected when the system is operating in **Turn** mode then **Turn** mode will be deselected and the indicator will show .

Unless deactivated, Max Flow will remain active for 15 minutes, during this time an audible warning will sound every three minutes. After this, the comfort level will revert to the previous setting, and remain in Static mode.

Turn Options


The  button selects the turn settings of the mattress. Each press of the button selects the next setting. For instance, to change the mattress from **Bi-lateral Turn** () to **Static** () position, press the button twice until  is selected.

In **Left Turn** () mode the right air cushions in the mattress will be held at a constant high pressure, and the left air cushions will be held at a constant low pressure. This process is reversed for the **Right Turn** () mode. When  is selected, the system will alternate between left and right turn, stopping in the central position.


During this time the display will show “C” or “d” to indicate the system is in the centralising phase.

Turn Times

The time selected for the turn is the amount of time (in minutes) the mattress will remain in each position during the selected turn sequence.

Turn Angle

Four settings are available to select the amount of turn given to the patient. These are selected using the turn angle button , where 4/4 is equal to approximately 40 degrees.

Comfort Control Level

These controls adjust the firmness of the mattress. Using the  button reduces the pressure setting, and the  button increases it.

The current comfort level (1-9) will be displayed in the LED display.

Lock Out

Pump functions (including power/standby switch) can be locked out to prevent unintentional change of modality / cell pressures. To activate, press the lock button for approximately 2-3 seconds until the blue LED illuminates. To deactivate, press and hold the lock out button until the LED turns off.

Alarms

Power Fail

In the event of a power fail situation, the pump will alert the carer by flashing the amber **Power Fail** LED and sounding the buzzer. Once the power is restored to the pump the alarm will cease and the pump will return to its previous operating settings.

Low Pressure

In the event of hose disconnection, an ‘L’ and ‘P’ will alternate on the LED display (representing "Low Pressure") while the buzzer sounds to alert the carer to the alarm condition. Once the hose is reconnected, the alarm will cease and the pump resume its previous therapy settings.

6. Operation

Before using the **Bari-Breeze** mattress system, make sure it has been installed correctly in accordance with “Installation and Setup” on page 6, and ensure that the side bolster and base plug is secure in the connector at the foot end of the mattress.

Initial Power Up Press the power/standby button to turn the pump on. An audible beep will sound. The system will go into a start up cycle for 30 seconds, during this time the pump will rapidly inflate the mattress. During this process the 7 segment display will alternate between ‘b’ and ‘f’ to indicate bolster fill.

 *These may reflect the therapy settings that were applied last time the system was in use.*

Comfort Controls Once the mattress is fully inflated and with the patient in position, adjust the comfort controls to the desired level. Adjust if necessary to ensure that the patient does not sink through to the base air pad.

Recommended Pressure Settings

For extra firm support during patient handling or nursing procedures, it is recommended to set the mattress pressure to maximum by pressing **MAX FLOW** (.

When the patient is to be placed in the gatched position, it is recommended to press the  button to *increase* the pressure/comfort level accordingly. to ensure the patient does not bottom out.

Static Mode

Max Flow Pressing the  button sets the blower to maximum flow, fully inflating the mattress overriding the comfort control. This feature is only available in static mode. If pressed during **Turn** mode, it will revert the system to **Static** mode and return the patient to a central position.

Max Flow inflates the mattress to maximum pressure for 15 minutes, and this will be indicated by a blue LED and an ‘F’ being shown in the LED display.

To cancel **Max Flow** and to enter **Turn** mode settings, press the  button once. To continue **Max Flow** for more than 15 minutes, press  once to cancel, and quickly again to enter **Max Flow** for a further 15 minutes. After 15 minutes, **Max Flow** is automatically cancelled.

 For extra firm support during patient addition/removal, wound care, patient turning, or patient hygiene procedures, it is recommended to set the mattress pressure to max by pressing  button.

Turn Options

Uni-Lateral Turn To enter uni-lateral **Turn** mode, press the  button to cycle through until  or  is indicated, this is confirmed by the illuminated blue LED.

Bi-lateral Turn To enter bi-lateral **Turn** mode, press the  button until  is indicated, this is confirmed by the illuminated blue LED.

 Selecting **Max Flow** will cancel any **Turn** options set, resulting in the system returning to the **Static** () setting.

Turn Times To set therapy turn times, press the  button until the desired time is indicated, this is confirmed by the illuminated blue LED.

 Turn times are only considered when the pump is operating in **Turn** mode.

To Deflate and Store the Bari-Breeze Mattress

To deflate the Bari-Breeze mattress:

1. Activate the CPR control to deflate the mattress.
2. Remove the side bolster and base air pad plug from the connector on the mattress hose assembly.

To store the mattress Following deflation:

Roll the mattress from the head end toward the foot end, ensuring the tubeset lies safely within the mattress.

7. Decontamination

The following processes are recommended, but should be adapted to comply with the local or national guidelines (Decontamination of Medical Devices) which may apply within the Healthcare Facility or the country of use. If you are uncertain, you should seek advice from your local Infection Control Specialist.

The **Bari-Breeze** system should be routinely decontaminated between patients and at regular intervals while in use; as is good practice for all reusable medical devices.

WARNING

Remove the power cord from the mains supply before cleaning and decontamination.

Caution

Do not use Phenol-based solutions or abrasive compounds or pads during the decontamination process as these will damage the surface coating. Do not boil or autoclave the cover.

Avoid immersing electrical parts in water during the cleaning process.

During Use *It is recommended that where practical, all parts of the system (including mattress, top cover, pump, and tubeset) should be cleaned weekly whilst in use.*

To clean Wipe all exposed surfaces with a cloth which has been moistened in simple detergent and water. Ensure to remove any organic debris that may be adherent to the surfaces. Dry.

To disinfect For disinfection we recommend a chlorine-releasing agent such as NaOCl or NaDCC diluted to 1,000 ppm (range may vary between 250 & 10,000 ppm depending on contamination status and local policy).

Alternatively use a cloth moistened with alcohol (70%) and allow to air dry.

Laundering The mattress cover can be easily unzipped for complete removal to allow laundering at 65°C for 10 minutes or 80°C for 3 minutes to achieve thermal disinfection. This complies with HSG (95) 18 Hospital Laundry arrangements for used and infected linen (UK) 1995.

DO NOT TUMBLE DRY COVERS ABOVE 30°C.

Intentionally left blank

8. Troubleshooting Guide

The following table contains fault symptoms, their possible causes and suggests steps to rectify the problem. These may occur during initial installation and power up.

Symptom	Possible Cause	Check	Action
Pump keyboard does not respond.	Lock out is on preventing set up.	Lock out button LED illuminated.	Press and hold the lock button until the LED extinguishes.
Mattress does not inflate when the system is switched on.	Pump is left in standby mode, indicated by a single dash on the LED display.	If the pump is plugged into a power socket, a dash will be visible on the LED display to indicate standby mode. Check status of LED display.	Turn on the pump unit. If the LED display is blank, check integrity of power supply at the wall socket.
	Blower failure within pump unit.	Check to ensure blower is generating airflow through the manifold connector.	If blower has failed call a Service Engineer to replace the blower/pump unit.
Mattress inflates but appears to be noisy as if a leak has appeared in the tube-set connection with the pump.	Hose-set connection not correctly fitted.	Check to ensure the hose connector clicks into place.	Re-connect hose connection to pump, ensuring firm fitment.
	Bolster/Sub-Mattress plug not fitted correctly or failing.	Check plug is fitted correctly in hose manifold. Check integrity of plug.	Ensure plug is firmly located. Replace if necessary.
Mattress inflates but doesn't appear to turn.	No turn angle set.	Therapy Settings.	Set Turn angles.
	Firm inflate activated.	Maxflow LED illuminated.	Press the  button to deactivate, and enter required therapy settings.
	Mattress remains central (check turn time).	Check status of LED display.	If LED display is illuminated with 'c' mattress is centralised.
Mattress inflates but only appears to turn to one side.	Therapy is set for left or right turn only.	Ensure bilateral therapy setting has been set.	Enter Bi-lateral therapy settings.
	Both setting may not be illuminated.	For bilateral turn ensure  is activated (LED illuminated).	If not, press  to select  and ensure the LED illuminates.
Mattress turns but then suddenly stops mid cycle.	Maxflow activated.	LED display showing 'F'.	Press the  button to deactivate, and set turn angles.
	Mattress remains central (check turn time).	Check status of LED display.	If LED display is illuminated with 'c' mattress is centralised.

Symptom	Possible Cause	Check	Action
Mattress appears to be bottoming out.	Comfort Control Settings too soft.	Assess comfort control settings.	Adjust as necessary.
	Patient Weight.	Check patient weight is within limit for system.	Refer to User Manual. Page 4 Indications.
Blower Unit Sounds noisier than usual.	Air leak from manifold.	Check connections.	Remove connection, and reconnect.
	Loose hose connector.	Check connections in tubeset connector.	Reconnect if possible, otherwise call Service Engineer.
	Torn or damaged cell.	Check.	Call service engineer to replace damaged cell.
Blower motor inside pump sounds excessively noisy during continuous use.	The blower inside the pump unit may be approaching the end of its life.	Check to ensure the mattress is fully inflated.	Call service engineer to request a replacement pump unit (if rental system).
Mattress appears to ignore therapy settings	Pump is not responding to therapy settings.	Check.	Reset Pump and re-enter therapy settings.

9. Spare Parts

To order spare items, please contact your distributor quoting the appropriate codes as shown below:

Pump	
Order Code	Description
647301	UK Bari-Breeze Pump
647303	USA Bari-Breeze Pump
647304	RoW Bari-Breeze Pump
647309AU	Australian Bari-Breeze Pump

Mattresses	
Order Code	Description
647311	Bari-Breeze Mattress (84" x 42")
647312	Bari-Breeze Mattress (80" x 39")
647313	Bari-Breeze Mattress (80" x 42")
647314	Bari-Breeze Mattress (80" x 48")

Mattress Covers	
Order Code	Description
647321	Bari-Breeze Mattress Cover (84" x 42")
647322	Bari-Breeze Mattress Cover (80" x 39")
647323	Bari-Breeze Mattress Cover (80" x 42")
647324	Bari-Breeze Mattress Cover (80" x 48")

10. Warranty and Service

ArjoHuntleigh's standard terms and conditions apply to all sales. A copy is available on request. These contain full details of warranty terms and do not limit the statutory right of the consumer.

11. Technical Data

Pump

Order Numbers:	647301 - UK Bari-Breeze Pump 647303 - USA Bari-Breeze Pump 647304 - RoW Bari-Breeze Pump 647309AU - Australian Bari-Breeze Pump
Supply Voltage:	100-240V
Electrical Rating:	50/60 Hz
Size:	Width: 33.5cm (13") Height: 28.5cm (11") Depth: 14cm (5.5")
Weight:	4.5kg (10lb)

Electrical Safety Standards

Complies with: EN 60601-1:1990/A13:1996

UL60601-1 1st Edition, CAN/CSA C22.2 No 601.1-M90

Medical Equipment

Degree of protection against electric shock: Class I, Type BF

Degree of protection of ingress of fluids: Ordinary IPx0 (not protected)

Mode of operation: Continuous

Pump Symbols

Dangerous Voltage

Refer to User Manual

Type BF

S/N:

Serial Number

Do not dispose of in a dustbin

Alternating Current

Environmental Conditions

Operating

Temperature range: +10°C to +40°C

Relative humidity: 30% to 75% (non-condensing)

Atmospheric pressure 700 hPa to 1060 hPa

Storage

Temperature range: -40°C to +70°C

Relative humidity: 10% to 100%

Atmospheric pressure: 500 hPa to 1060 hPa

Environmental Protection: Please dispose of this pump in accordance with local regulations.

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NORGE

ArjoHuntleigh Norway AS
Olaf Helsets vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail:
no.kundeservice@arjohuntleigh.com

ÖSTERREICH

ArjoHuntleigh GmbH
Dörrstrasse 85
AT-6020 INNSBRUCK
Tel: +43 (0) 512 204 160 0
Fax: +43 (0) 512 204 160 75

POLSKA

ArjoHuntleigh Polska Sp. z o.o.
ul. Ks Piotra Wawrzyniaka 2
PL-62-052 KOMORNIKI (Poznan)
Tel: +48 61 662 15 50
Fax: +48 61 662 15 90
E-mail: arjo@arjohuntleigh.com

PORTUGAL

ArjoHuntleigh em Portugal
MAQUET Portugal, Lda. (Distribuidor Exclusivo)
Rua Poeta Bocage n.º 2 - 2G
PT-1600-233 Lisboa
Tel: +351 214 189 815
Fax: +351 214 177 413
E-mail: Portugal@arjohuntleigh.com

SUISSE / SCHWEIZ

ArjoHuntleigh AG
Fabrikstrasse 8
Postfach
CH-4614 HÄGENDORF
Tél/Tel: +41 (0) 61 337 97 77
Fax: +41 (0) 61 311 97 42

SUOMI

Oy Vestek AB
Martinkuja 4
FI-02270 ESPOO
Puh: +358 9 8870 120
E-mail: info@vestek.fi

SVERIGE

ARJO Scandinavia AB
Hans Michelsensgatan 10
SE-211 20 MALMÖ
Tel: +46 (0) 10 494 7760
Fax: +46 (0) 10 494 7761
E-mail: kundservice@arjohuntleigh.com

UNITED KINGDOM

ArjoHuntleigh UK
ArjoHuntleigh House
Houghton Hall Park
Houghton Regis
UK-DUNSTABLE LU5 5XF
Tel: +44 (0) 1582 745 700
Fax: +44 (0) 1582 745 745
E-mail: sales.admin@arjohuntleigh.com

USA

ArjoHuntleigh Inc.
2349 W Lake Street Suite 250
US-Addison, IL 60101
Tel: +1 630 307 2756
Free: +1 800 323 1245 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 630 307 6195
E-mail: us.info@arjohuntleigh.com

REV 15: 02/2015

GETINGE GROUP

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **Getinge** and **Maquet**. **Getinge** provides solutions for infection control within healthcare and contamination prevention within life sciences. **Maquet** specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

ARJOHUNTLEIGH
GETINGE GROUP

ArjoHuntleigh focuses on patient handling and hygiene, disinfection, DVT prevention, medical beds, therapeutic surfaces and diagnostics.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjohuntleigh.com