ARJOHUNTLEIGH

GETINGE GROUP

Proficare Mattress Replacement System USER MANUAL

CE
0086

P/N M6252209-AH Rev C • 02/2014

...with people in mind

User Manual	3
Benutzerhandbuch	37
Manuel d'utilisation.....	77
Manuale d'uso	115
Gebruikershandleiding.....	153
Brugervejledning.....	193
Användarhandbok.....	231

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

ARJOHUNTLEIGH HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE ARJOHUNTLEIGH PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY ARJOHUNTLEIGH SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL ARJOHUNTLEIGH BE LIABLE FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND ARJOHUNTLEIGH TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in ArjoHuntleigh printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties, except as set forth in the written limited warranty included in this publication or with this product. Information in this publication may be subject to change at any time. Consult the manufacturer for updates.

Important Information for Users

In order for ArjoHuntleigh products to perform properly, ArjoHuntleigh recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by ArjoHuntleigh. Contact ArjoHuntleigh for information regarding maintenance and repair.
- Ensure the electrical installation of the room complies with the appropriate national electrical wiring standards. To avoid the risk of electric shock, this product must be connected to a grounded power receptacle.

Specific indications, contraindications, warnings, precautions and safety information exist for ArjoHuntleigh's therapeutic support surface products. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

Notice

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the product information label for specific voltage.

TABLE OF CONTENTS

Introduction.....	8
About the <i>Proficare Mattress Replacement System</i>	8
Clinical Applications	9
Indications.....	9
Contraindications.....	9
Risks and Precautions	9
Transfer.....	9
Side Rails and Restraints.....	9
Patient Migration.....	9
Use With Other Devices.....	9
Safety Information.....	10
Side Rails / Patient Restraints.....	10
Bed Frame	10
Bed Height.....	10
Brakes.....	10
Head of Bed Elevation	10
Patient Entrance / Exit	10
Skin Care.....	11
No Smoking in Bed.....	11
Patient Weight	11
Therapy Unit	11
Interference.....	11
General Protocols	11
Disposal.....	12
Protection Against Hazards.....	12
Fluids	12
Power Cable.....	12
Enclosure	12
Product Description.....	13
General	13
Therapy Unit.....	14
Mattress.....	15
Controls and Operation.....	16
Front Panel	16
Mode Selection.....	16
<i>Soft Key Comfort Control Level Adjustment</i>	17
Travel Mode.....	17
CPR.....	18
Removable Air Cells.....	18

Alarms and Mode Indicators	19
Alarms	19
Low Pressure Alarm	19
High Pressure Alarm	19
Power Fail Alarm	19
Alternation Failure Alarm	19
Mode Selection Indicators	20
Alternating Mode	20
Static Mode	20
Travel Mode	20
Soft Key Comfort Control Level Selection Indicators	21
Self-Test Function	22
Product Usage Timer	23
Power Down Procedure	23
Preparation for Use	24
Check Out Procedures	24
Attach the Quick Reference Guide (QRG)	24
Install the <i>Proficare Therapy Unit</i> and Mattress	24
Patient Placement	26
Care and Cleaning	28
Mattress Cover	28
Inflatable Air Cells	28
Foam Cover	28
Foam	28
Other Mattress Components	28
Therapy Unit	28
Storage	29
Maintenance	29
Replacement Parts	29
Specifications	30
Therapy Unit	30
Mattress	30
Dimensions (90 cm)	30
Dimensions (80 cm)	30
Transport and Storage Conditions	31
Customer Contact Information	31
Electromagnetic Emissions Information	32
Definition of Symbols Used	35

INTRODUCTION

About the *Proficare™ Mattress Replacement System*

The *Proficare Mattress Replacement System* is a fully automatic mattress replacement system (MRS) suitable for the treatment and prevention of pressure ulcers.

The *Proficare MRS* supports patients weighing from 40 kg to 185 kg.

The *Proficare MRS* is designed for use on a hospital bed frame.
Please contact the frame manufacturer regarding individual frame capacity.

The *Proficare MRS* provides both alternating and static patient support surfaces.

The *Proficare MRS* utilizes an auto-control system which self adjusts the mattress air pressure dependent on patient weight and posture. The *Soft Key Comfort Control™* feature allows patient comfort to be optimized without interrupting operation of the auto-control system in either Alternating or Static modes.

The *Proficare MRS* facilitates in situ patient transport.

The *Proficare MRS* incorporates a self test function assuring the user of correct therapy unit and mattress operation.

The *Proficare MRS* is easy to maintain and clean.

CLINICAL APPLICATIONS

It is recommended that the patient be assessed for risk of developing pressure ulcers prior to use of the *Proficare MRS* and continued monitoring of their skin condition be carried out thereafter in accordance with your established procedures.

Indications

This product is intended to help treat and reduce the incidence of pressure ulcers while optimizing patient comfort. The *Proficare MRS* is indicated for:

- Individuals at risk, or suffering from pressure ulcers.
- Pain management as prescribed by a physician.
- Any patient using a standard mattress that would benefit from an Alternating or Static Pressure Management System.

Contraindications

Patient conditions for which the pressure relieving therapy of the *Proficare MRS* are contraindicated include:

- Cervical or skeletal traction
- Unstable spinal cord injuries

Risks and Precautions

Transfer

Standard precautions should be taken during patient transfer.

Side Rails and Restraints

WARNING: Use or non-use of restraints, including side rails, can be critical to patient safety. Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. See related Safety Information.

Patient Migration

Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. Monitor patients frequently to guard against patient entrapment.

Use With Other Devices

All *Proficare MRS* components are designed to be used as single system. The *Proficare Therapy Unit* should only be used with the *Proficare MRS*. Any attempt to connect and use the therapy unit with any other brand or model of mattress or overlay may result improper operation of equipment, possibly leading to increased risk of patient injury.

Safety Information

Side Rails / Patient Restraints

Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patients frequently to guard against patient entrapment.

When selecting a standard mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

Bed Frame

Always use a standard healthcare bed frame with safeguards or protocols that may be appropriate. It is recommended that bed and side rails (if used) comply with all applicable regional and local regulations and protocols. Frame and side rails must be properly sized relative to the mattress to help minimize any gaps that might entrap a patient's head or body. Inappropriate, incorrect, poorly fitted or poorly maintained bed rails and mattress overlays or replacements can potentially produce gaps that might entrap a patient's head or body.

Bed Height

To minimize risk of falls or injury, the bed should always be in the lowest practical position when the patient is unattended.

Brakes

Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

Head of Bed Elevation

Keep head of bed as low as possible to help prevent patient migration.

Patient Entrance / Exit

Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Skin Care

Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

No Smoking in Bed

Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

Patient Weight

The maximum patient weight for this system is 180 kg.

The minimum patient weight for this system is 40 kg.

Therapy Unit

- Only use a grounded power receptacle and the power cable supplied with the *Proficare Therapy Unit*. The power cable should be positioned to avoid a tripping hazard and / or damage to the cable. ArjoHuntleigh recommends placing the cable under the bed frame and attaching it to a mains wall socket by the head of the bed. The system should never be operated with a worn or damaged power cable. Should the power cable be found to be worn or damaged, contact ArjoHuntleigh or an ArjoHuntleigh authorized representative for a replacement.
- Only use replacement fuses that have the same specified rating. Using fuses with higher ratings could result in damage or injury.
- The therapy unit is a precision electronic product. Use care when handling or transporting. Dropping or other sudden impacts may result in damage to the therapy unit .
- Do not open the therapy unit – risk of electrical shock. Do not attempt to repair or service the therapy unit . Repairs and service should be conducted only by ArjoHuntleigh or an ArjoHuntleigh authorized representative. If the therapy unit is not functioning properly or has been damaged, unplug the unit and take it out of service immediately. Contact ArjoHuntleigh or an ArjoHuntleigh authorized representative for repair and service information.
- Do not place any objects or items, such as blankets, on or over the therapy unit .
- The air intakes for the *Proficare Therapy Unit* are on the back of the unit. To avoid the risk of blocking these intakes, it is recommended that the unit be mounted on the footboard of a bed frame. Should the bed not have a suitable footboard, it is recommended that the unit be placed on a stationary surface.

Interference

Although the *Proficare Therapy Unit* conforms to the intent of IEC 60601-1-2 in relation to electromagnetic compatibility, all electrical equipment may produce interference. If interference is suspected, move equipment from sensitive devices or contact the manufacturer. (IEC 60601-1-2. Medical Electrical Equipment – Part 1: General Requirements for Safety, Amendment No. 2. Collateral Standard: Electromagnetic Compatibility Requirements and Tests). See the Electromagnetic Emissions Information section for detailed testing data.

General Protocols

Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

Disposal

Dispose of the air filter, cushions and mattress cover sheet according to local procedures and regulations. At the end of useful life, dispose of waste according to the European Community Waste Electrical and Electronic Equipment Directive (WEEE Directive) and in compliance with relevant local regulations, or contact ArjoHuntleigh or an ArjoHuntleigh authorized representative for advice.

Protection Against Hazards

Fluids

Avoid spilling fluids on any part of the therapy unit . If spills do occur:

- Disconnect the unit from the main electricity supply
- Clean fluids from the case

Ensure that there is no moisture in or near the power inlet, power switch and power plug before reconnecting the power supply.

- Check the operation of controls and other components in the area of the spill
- Perform applicable checkout procedures

Fluids remaining on the electronic controls can cause corrosion that may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and staff.

Power Cable

The system should never be operated with a worn or damaged power cable. Should the power cable be found to be worn or damaged, contact ArjoHuntleigh or an ArjoHuntleigh authorized representative for a replacement.

Enclosure

Prolonged contact with pump unit may result in hazard to patient or caregiver.

PRODUCT DESCRIPTION

General

The *Proficare MRS* is an Alternating Pressure Mattress Replacement System with an additional Static Mode feature. The system is comprised of a therapy unit and a mattress.

During Alternating mode, the three head end air cells are statically inflated while the air cells supporting the body are sequentially inflated and exhausted throughout a ten minute cycle.

During Static mode, all air cells are simultaneously inflated thereby supporting the patient over a greater surface area.

The auto-control system self adjusts air cell pressure dependent on the selected mode, patient weight and patient posture.

Therapy Unit

The therapy unit is constructed of an impact resistant molded enclosure with anti slip pads to the base and rear.

All controls and indicators are flush mounted and situated on the front panel for ease of viewing.

An integral carrying handle / mounting bracket is provided allowing the unit to be fixed to a bed footboard.

All functions of the therapy unit and air pressures within the mattress are microprocessor controlled.

A quick reference guide is attached to the side of the therapy unit.

Mattress

The mattress features seventeen individual air cells. The air cells are constructed using polyurethane-coated nylon materials and are designed to provide comfort and pressure relief.

A foam layer with protective cover underlies the air cells for additional comfort.

A pressure-sensing pad is located between the foam layer and air cells. This pad operates as part of the auto-control system.

The mattress cover is manufactured from a two-way stretch knitted fabric with a vapor permeable polyurethane coating.

The air supply is connected from the therapy unit to the mattress by means of a single quick release connector.


The mattress incorporates a CPR release situated at the head end of the mattress.

Individual air cells may be removed during patient use.

CONTROLS AND OPERATION

Front Panel

Mode Selection

Alternating, Static and Travel Modes are selected by pressing the Mode Selection key until the relevant indicator is lit.

Static Mode cannot be selected during initial inflation.

Soft Key Comfort Control Level Adjustment

Comfort levels allow the user / caregiver to optimize patient comfort by altering the sensitivity of the automatic pressure control system.

Comfort Levels 1 (extra-soft) to 4 (firm) are selected by pressing the *Soft Key Comfort Control* Key until the relevant indicator is lit.

Comfort levels cannot be selected during Travel Mode.

Travel Mode

The mattress may be kept inflated for in situ patient transport as follows:

1. Press the Mode Selection key until the Travel Mode indicator is lit. After a few seconds the Wait indicator will also be lit.
2. Wait approximately 90 seconds until the (green) Travel Move indicator flashes with an audible signal; the therapy unit may then be switched OFF and the mains cord removed.

The mattress will now remain inflated for a minimum of three hours.

Following patient transport, the therapy unit power may be restored. Alternating Mode and Comfort Level 2 will be selected by default.

CPR

Pull downward on the CPR release strip situated at the head end / patient's left side of the mattress. The mattress will be sufficiently exhausted in approximately 15 seconds.

Removable Air Cells

To further reduce tissue pressure at specific areas, up to two of the six foot-end air cells may be removed. Caps are provided within the base of the mattress to close the air connector from which the air cell was removed. Therapy may then be continued.


ALARMS AND MODE INDICATORS

Alarms

All alarms are classified as low priority alarms.

Low Pressure Alarm


Under low pressure fault conditions the Low Pressure indicator will flash and the audible alarm will sound a single intermittent beep. If the alarm is activated, check all air connections are properly made and the equipment is installed as per the instructions in this manual. To reset, switch therapy unit OFF for a few seconds and back ON. If the alarm is activated again refer to the Self Test section of this manual or contact ArjoHuntleigh.

During initial inflation the Low Pressure and Wait indicators flash alternately, with the selected mode indicator on continuously. The audible low pressure signal is automatically muted.

High Pressure Alarm


Under high pressure fault conditions the High Pressure indicator will flash and the audible alarm will sound a continuous beep. Further inflation of the mattress is automatically prevented. If this alarm is activated, check the equipment is installed as per the instructions in this manual. To reset, switch therapy unit OFF for a few seconds and back ON. If the alarm is activated again refer to the Self Test section of this manual or contact ArjoHuntleigh.

Power Fail Alarm

If the therapy unit is switched OFF, or the power supply is removed for any reason, a continuous audible alarm will sound. Press the Mute Key to cancel the alarm.


Alternation Failure Alarm

If the alternating pressure therapy is not fully operational during Alternating Mode, the Alternating Mode indicator will flash and the audible alarm will sound a single intermittent beep.


Mode Selection Indicators

Alternating Mode


When Alternating Mode is selected the indicator illustrated above is lit. The therapy unit defaults to this mode when power is first applied.

Static Mode


When Static Mode is selected the indicator illustrated above is lit.

Travel Mode


When the Travel Mode has been selected the Travel Mode and Wait indicators illustrated above will light until the mattress achieves sufficient pressure for patient transport.


When the mattress has achieved sufficient pressure for patient transport the Wait indicator will go out, the Move indicator will flash and an intermittent double beep will sound.


Soft Key Comfort Control Level Selection Indicators


When Comfort Level 1 (extra-soft) is selected the indicator illustrated above is lit.

When Comfort Level 2 is selected the indicators illustrated above are lit. The therapy unit defaults to this setting when first powered up.

When Comfort Level 3 is selected the indicators illustrated above are lit.

When Comfort Level 4 (firm) is selected the indicators illustrated above are lit.


Adjusting comfort levels varies only the sensitivity of the Automatic Pressure Controller. The mattress pressure continues to be automatically controlled at all times.

Self-Test Function

1. SelfTest should be performed with no patient on the mattress. Ensure the mattress and therapy unit are connected in accordance with the Preparation for Use chapter of this manual.
2. To activate the SelfTest Function, switch on the therapy unit and press and hold down the Mode and *Soft Key Comfort Control* Key for five seconds.

All indicators will flash throughout the duration of the self test. The test will continue for approximately 12 - 15 minutes. (If the self test was requested during initial inflation, the remainder of the initial inflation time will be in addition to the self test time). During the test, do not disturb or interfere with the therapy unit or mattress as this may effect the test results.

3. A test result Pass will be indicated by a series of short intermittent beeps. The therapy unit may now be switched OFF and will revert to normal operation when power is next applied.
4. A test result Failure will be indicated by a long beep followed by one or more short beeps. If this alarm is activated, switch OFF the therapy unit and repeat the self test to confirm the result. If the fault persists contact ArjoHuntleigh.


One long beep followed by one short beep indicates a minor leakage. The Travel Mode holding pressure may be reduced to less than three hours, but all other modes will operate normally. A full list of alarm codes and diagnosis is available on request from ArjoHuntleigh.

Product Usage Timer

The total number of hours the *Proficare MRS* has been in use may be viewed as follows:

1. Hold the Travel key down while switching ON the therapy unit.

Comfort Level indicators 1 to 4 will sequentially flash accompanied by a series of beeps corresponding to the total hours of product usage.

Example:	X 1000 Hours	12 beeps and flashes
	X 100 Hours	5 beeps and flashes
	X 10 Hours	0 beeps or flashes
	X 1 Hour	2 beeps and flashes
		Total = 12,502 Hours

When the sequence has been completed a continuous beep will sound.

2. The therapy unit may now be switched OFF and will revert to normal operation when power is next applied, or will automatically revert to normal operation after 30 seconds.

Power Down Procedure

1. Press the power switch located at the side of the therapy unit to the OFF position.
2. Unplug the mains cord from the power supply.

PREPARATION FOR USE

Check Out Procedures

The unit is shipped from ArjoHuntleigh ready to be put into service. Before putting the unit into service, the following steps should be performed.

1. Remove the mattress and the therapy unit from their packaging.

Do not use sharp instruments to open packaging. Damage could result to mattress, cover sheet, therapy unit or mattress hose.

2. Inventory and visually inspect that the following items are present and in good condition:
 - One *Proficare Therapy Unit*
 - One *Proficare Mains Cord*
 - One *Proficare Mattress*
 - One *Proficare Hold-All*
 - One *Proficare User Manual*
 - One *Proficare Quick Reference Guide*

Contact carrier and ArjoHuntleigh immediately if damage is observed.

Attach the Quick Reference Guide (QRG)

Position the hole in the QRG over the pump-mattress connector socket and snap the locating tabs into the groove. The QRG may be removed for cleaning if required.

Install the *Proficare Therapy Unit* and Mattress

1. Remove the existing mattress and place the *Proficare Mattress* directly onto the bed frame (or base) and hang the therapy unit at the foot of the bed. (Alternatively place the therapy unit on the floor).

2. Use the quick release connector to attach the air supply hoses from the mattress onto the therapy unit.

3. Ensure CPR release is in the closed position as shown at right. Keep the CPR airways clear of obstruction and the pull strip accessible at all times.

4. Insert the socket of the mains cord into the appliance inlet on the therapy unit.
5. Connect to a suitable power supply.
6. Press the power switch located at the side of the therapy unit to the ON position. The Alternating Mode indicator will be lit, with the Wait and Low Pressure indicators flashing alternately. During initial inflation, Alternating Mode is entered by default although the Alternating Cycle will not commence until initial inflation is complete.

7. Allow approximately 30 minutes for initial inflation to be completed. The *Proficare MRS* is ready for use when the Wait and Low Pressure indicators are no longer lit.

When the patient is initially placed on the inflated mattress, the mattress pressure may at first appear firm. The automatic pressure controller will reduce this pressure in approximately 5 - 10 minutes.

PATIENT PLACEMENT


Before placing a patient on the *Proficare MRS*, it is recommended that personnel familiarize themselves with the information outlined in the Safety Information section and Controls and Operation chapter of this manual as well as the procedures outlined below.

1. Allow the mattress overlay or mattress replacement to fully inflate before placing the patient. Failure to do so could cause patient injury.
2. Once the mattress is fully inflated (as described in the previous section), bedding can be replaced. Secure sheets loosely enough to ensure there is no interference with air cell alternation. Use incontinence pads (such as *Dri-Flo™* or similar) to avoid gross contamination.
3. Select desired mode:
 - The Alternating Mode is used for normal therapeutic function. The air cells are inflated (to the set pressure) and deflated following a fixed cycle time. This provides regular periods of pressure reduction to aid blood and lymphatic flow to vulnerable tissue.
 - The Static Mode is used to support the patient but without alternating inflation and deflation. During normal use, Alternating Mode should be active.
4. To ensure patient safety, assist the patient onto the mattress. Ensure patient is centered side-to-side and patient's head is properly placed on the head cushions.
5. Use the Pressure Buttons (up and down arrows) to reach the desired setting for effective patient comfort and support: 1 = lowest pressure setting, 4 = highest pressure setting.

Pressure Level 1 (lowest pressure setting) is applicable for light weight patients.
Pressure Level 4 (highest pressure setting) is applicable for heavy weight patients.

6. When the desired comfort level is chosen in relation to patient weight, ensure the mattress is supporting the patient effectively. With the patient in a supine position, slide a hand under the patient's sacral area beneath the air cells. The patient should be supported on the inflated cells with a distance of 3 - 5 cm from the base mattress. If the distance is less than 3 - 5 cm or if the patient is touching the base mattress, readjust the pressure level to a firmer setting and check again.

When lowering the pressure, ensure the mattress is still supporting the patient effectively. Failure to do so could cause patient injury.

When the patient is initially placed on the inflated mattress, the mattress pressure may at first appear firm. Wait at least ten minutes between any pressure adjustments or patient assessment, as the air cells may take necessary cycling time to adjust to the new pressure setting.

7. Re-check after approximately 20 minutes that the patient is comfortable and that the orange LED on pressure setting and the and green light on the therapy setting are illuminated to indicate the *Proficare MRS* is functioning correctly.
8. If the patient is moved to a more upright sitting position, the pressure may need to be increased by one or two comfort settings to provide effective support. Once pressure is increased, slide a hand under the patient's sacral area beneath the air cells. Patient should still be supported on the inflated cells with a distance of 3 - 5 cm from the base mattress.

When the patient returns to the supine position, it is important to return to the original pressure setting and ensure the patient is properly supported.

CARE AND CLEANING

Always disconnect the *Proficare Therapy Unit* from mains wall socket before cleaning. Failure to do so could result in equipment damage and / or electric shock.

Mattress Cover

Wipe down the *Proficare MRS* using a well wrung cloth dipped in warm soapy water. The *Proficare MRS* may then be disinfected using approved hospital disinfectant diluted according to the manufacturer's instructions. Rinse with plain water and allow to air dry before returning to service.

If the cover sheet has become excessively soiled, it may be laundered in a washing machine at up to 95°C using normal domestic washing powder. Rinse thoroughly using plain water. Dry on low temperature setting and ensure cover sheet is completely dry before returning to use.

Inflatable Air Cells

1. Remove the mattress cover and detach the air cells from their fasteners at both sides of the mattress.
2. Detach the air cells from the air connectors at the patient left side of the mattress and withdraw each air cell from the loops of the foam cover.
3. Wipe down each air cell with water and mild detergent.
4. If required, disinfect the air cells using approved hospital disinfectant diluted according to the manufacturer's instructions.

Foam Cover

1. Remove the air cells as above.
2. Detach the foam cover fasteners from the mattress base assembly fasteners and slip the foam out of each end of the cover.
3. Wipe down the cover with water and mild detergent.
4. If required, disinfect the foam cover using approved hospital disinfectant diluted according to the manufacturer's instructions.

Foam

The foam should be replaced when soiled or contaminated.

Other Mattress Components

All other mattress components may be wiped down with water and mild detergent. If required they may then be disinfected using approved hospital disinfectant diluted according to the manufacturer's instructions.

Therapy Unit

CAUTION! Do not immerse or soak therapy unit.

1. Wipe down with a cloth dampened in water and mild detergent.
2. If required, disinfect the therapy unit using approved hospital disinfectant diluted according to the manufacturer's instructions.

Storage

1. Pull downwards on the CPR Release Strip to release the air from the mattress.
2. Roll the mattress from the head end towards the foot end expelling air as it is rolled and refit the CPR Release for next use.
3. Fasten the attached storage straps to prevent unrolling and store or transport the mattress in the hold-all provided.

Maintenance

The maintenance tasks listed below should be carried out before each use.

1. Ensure mattress cover and air cells are securely fastened.
2. Ensure the Pressure-Sensing Pad is securely fastened and free of creases.
3. Check the mains supply cord for wear or fraying etc. and replace if necessary.
4. Run Self Test Function to confirm correct operation of the therapy unit / mattress combination.

The maintenance tasks listed below should be carried out at 12 month intervals.

1. Check the mattress cover for signs of wear or damage.
2. Check the air supply tubes and pump / mattress connector O rings for signs of wear or damage.
3. Lightly lubricate the pump / mattress connector O rings with silicon grease.
4. Run Self Test Function to confirm correct operation of the therapy unit / mattress combination.

Replacement Parts

<i>Proficare Mattress (90 cm)</i>	M6265524
<i>Proficare Mattress (80 cm)</i>	M6265525
<i>Proficare Therapy Unit (230V)</i>	M6257759
<i>Proficare Therapy Unit (230V) - 60601 3rd edition</i>	413411
<i>Proficare Therapy Unit (110V)</i>	M6257762
<i>Proficare Hold-All</i>	M6248433
<i>Proficare User Manual</i>	M6252209-AH
<i>Proficare Quick Reference Guide, English</i>	M3252561-AH
<i>Proficare Quick Reference Guide, German</i>	M3252559-AH
<i>Proficare Quick Reference Guide, French</i>	M3252563-AH
<i>Proficare Quick Reference Guide, Italian</i>	M3252562-AH
<i>Proficare Mains Cord, UK</i>	M4268743
<i>Proficare Mains Cord, Europe</i>	M4268742
<i>Proficare Mains Cord, Canada</i>	M4268840
<i>Proficare Mains Cord, Australia</i>	M4268893

SPECIFICATIONS

Specifications subject to change without notice.

Therapy Unit

Width 29 cm
Height 24 cm
Depth 12 cm
Weight 3.5 kg

Enclosure Material ABS / Polycarbonate blend
Alternating Pressure Cycle 10 minutes

	EUROPEAN	CANADIAN
Supply Voltage	220 - 240 V~ 50 / 60 Hz	100 - 120 V~ 60 Hz
Current	137.5mA (220V) / 150mA (240V)	275 mA (100V) / 330 mA (120V)
Internal Fuse Rating	F315 mA	F500 mA
Mains Plug Fuse Rating	3 A (UK Only)	N/A

Degree of protection against electric shock Class II Type B
Degree of protection against ingress of fluids Ordinary (not protected)
Mode of operation Continuous
Range of intended use -10°C to +40°C

Mattress

Dimensions (90 cm)

Length 195 cm
Width 90 cm
Depth 20 cm (inflated)
Weight 11 kg

Dimensions (80 cm)

Length 195 cm
Width 80 cm
Depth 20 cm (inflated)
Weight 10 kg

Mattress Cover Materials Polyurethane coated knitted nylon
Base Materials Polyurethane coated woven nylon
Inflatable Air Cells Polyurethane coated woven nylon
Foam Cover Polyurethane
Air Supply Tubes Silicon
Pump / Mattress Connector Acetal

Transport and Storage Conditions

Ambient temperature range-40°C to +70°C
Relative humidity range10% to 100% including condensation
Atmospheric pressure range500 hpa to 1060 hpa.

Specifications subject to change without notice.

The *Proficare Mattress Replacement System* conforms to the following International Standards: EN60601-1:1990 including A13:1996, CAN/CSA-C22.2 No 601.1M90, UL2601-1:1994 1st Edition-Amended 1996, EN60601-1-2:1993, 60601-1:2005, 3rd Ed.

The *Proficare Mattress Replacement System* is classified as a Type B applied part under IEC 60601-1:2005.

CUSTOMER CONTACT INFORMATION

For questions regarding this product, supplies, maintenance or additional information about ArjoHuntleigh products and service, please contact ArjoHuntleigh or an ArjoHuntleigh authorized representative or visit www.ArjoHuntleigh.com.

ELECTROMAGNETIC EMISSIONS INFORMATION

Electromagnetic Interference - Although this equipment conforms with the intent of the directive 89/336/EEC in relation to Electromagnetic Compatibility (EMC), all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.

Portable and mobile RF communications equipment can effect medical electrical equipment. Radios, cell phones and similar devices may affect this equipment and should be kept at least 2 m (6.5 ft) away from the equipment.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the following tables.

The following tables document compliance levels and guidance from the IEC 60601-1-2 2007 Standard, for the electromagnetic environment in which the *Proficare MRS* should be used in a clinical environment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The <i>Proficare MRS</i> is intended for use in the electromagnetic environment specified below. The customer or user of the <i>Proficare MRS</i> should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	The <i>Proficare MRS</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations / flicker emissions IEC 61000-3-3	yes	

Guidance and manufacturer's declaration - electromagnetic immunity			
The <i>Proficare MRS</i> is intended for use in the electromagnetic environment specified below. The customer or user of the <i>Proficare MRS</i> should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV Contact ±8 kV Air	±6 kV Contact ±8 kV Air	In accordance with IEC 60601-1-2:2007, floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	±1 kV	±1 kV	
Surge IEC 61000-4-5	1 kV line(s) to Line(s) 2 kV line(s) to Earth	1 kV line(s) to Line(s) 2 kV line(s) to Earth	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% half cycle 40% for 5 cycle 70% for 25 cycle 5% for 5 seconds	5% half cycle 40% 5 cycles 70% 25 cycles 5% for 5 seconds	
Power frequency (50Hz / 60Hz) magnetic field IEC 61000-4-8	3 A/M	3 A/M	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U _i is the a.c. mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the <i>Proficare MRS</i>			
The <i>Proficare MRS</i> is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <i>Proficare MRS</i> can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <i>Proficare MRS</i> as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter meters		
	150 kHz to 80 MHz not applicable	80 MHz to 800 MHz $d = \left[\frac{3.5}{E1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{2}{E1} \right] \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.37	0.74
1	N/A	1.2	2.3
10	N/A	3.7	7.4
100	N/A	12	23
For transmitters rated at a maximum output power not listed above, the recommended separate distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1, At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects and people.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The <i>Proficare MRS</i> is intended for use in an electromagnetic environment specified below. The customer or user of the <i>Proficare MRS</i> should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3Vrms 150K - 80 MHz 3 Vrms 80 MHz - 2.5 GHz	3 Vrms 150K - 80 MHz 3 Vrms 80 MHz - 2.5 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the <i>Proficare MRS</i>, including cables, than the recommended separation distance calculated from the equation application to the frequency of the transmitter. Recommended Separation Distance Battery Operated Device $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (see note a) should be less than the compliance level in each frequency range (see note b). Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1, At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the <i>Proficare MRS</i> is used exceeds the applicable RF compliance level above, the <i>Proficare MRS</i> should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the <i>Proficare MRS</i>. b) Over the frequency range 150kHz, field strengths should be less than [V1] V/m.			

DEFINITION OF SYMBOLS USED

	Warning of possible hazard to system, patient or staff		ETL Listed, Conforms to AAMI ES60601-1 1st Edition, CSA C22.2#60601-1 3rd Edition and IEC 60601-1 3rd Edition
	Important operational information	IPX0	Not protected against harmful effects of water
	Tumble dry at low heat		Do not iron
	May be dry-cleaned with solvent classified with the letter P		Washable up to 95° C
	Do not bleach		Foot End
I	On		Type B Applied Part
	Off		The <i>Proficare MRS</i> conforms to the Medical Devices Directive (93/42/EEC).
	Date of Manufacture		Refer to Instruction Manual
	Manufacturer		This product is designated for separate collection at an appropriate collection point.
	Static Mode		Move
	Alternating Mode		Low Pressure
	Travel Mode		High Pressure
	Wait		Audible Signal
	Muted Audible Signal	REF	Reorder Number
	Class II Device		Temperature Limitations
SN	Serial Number		
	Alternating Current		

ARJOHUNTLEIGH

GETINGE GROUP

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **GETINGE** and **MAQUET**. **ArjoHuntleigh** focuses on patient mobility and wound management solutions, **GETINGE** provides solutions for infection control within healthcare and contamination prevention within life sciences. **MAQUET** specializes in solutions, therapies and products for surgical interventions and intensive care.

www.ArjoHuntleigh.com

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden

* and ™ are trademarks belonging to the ArjoHuntleigh group of companies
© ArjoHuntleigh 2014. P/N M6252209-AH Rev C 02/2014

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ARJOHUNTLEIGH

GETINGE GROUP

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **GETINGE** and **MAQUET**. **ArjoHuntleigh** focuses on patient mobility and wound management solutions. **GETINGE** provides solutions for infection control within healthcare and contamination prevention within life sciences. **MAQUET** specializes in solutions, therapies and products for surgical interventions and intensive care.

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