

UPGRADED DESIGN ON THE PATIENT-SPECIFIC DISPOSABLE SLINGS

Standard Loop Flites®

Always ready to use with no reprocessing requirements our new upgraded patient-specific disposable Standard Loop Flites® are designed to support infection control, streamline safe patient handling tasks, and optimise patient comfort.

New design features

A - New head support design, featuring extra soft integrated padding and minimal pressure points for patient comfort

B - Larger surface at the **sitting area** and **around the shoulders** for increased patient support

C - New strap design

- Longer for better patient head positioning in relation to the spreader bar
- Three loops on the shoulder straps, for even more patient positioning possibilities
- Slimmer and size colour-coded for ease of use

An additional XXL size to cater to a wider range of patients

- Size range from S - XXL and colour-coded



New slimmer strap design. Extra loops on shoulder straps, for even more patient positioning possibilities



All slings are size colour-coded for ease of use. Now an extra XXL size is added to the range (S, M, L, XL and XXL)

10.KX.02.1GB-INT.1.ARJO