

## MULTIPURPOSE HYGIENE CHAIR


# Carendo®



DSDC Product Accreditation  
Rating Class 1B

The Carendo® multipurpose hygiene chair facilitates an efficient hygiene and showering routine within a single transfer while allowing a comfortable working posture for the caregiver throughout the process.

Designed with Carl and Doris in mind



### Dignified care

Care Raiser™ function raises the patient or resident's lower body to enable access for performing discreet hygiene tasks while maintaining the patient or resident's dignity.

### Streamlined hygiene routine

Designed to enable a single caregiver to perform the full hygiene routine, including dressing and undressing, toileting, showering and other hygiene tasks, in a single transfer.

### Enables better working postures

Carendo has been shown to result in caregivers adopting 77% neutral back positions during hygiene tasks, compared to only 48% when washing a person in the bed.\*

\*Source: Knibbe JJ, Knibbe NE & Heitink DEBL (2016).

## Key features

### HEIGHT-ADJUSTABLE

Allows carer to reach preferred working height throughout the hygiene routine.

### FULL LENGTH BACKREST

### HANDHELD CONTROL UNIT

### EMERGENCY STOP

### EMERGENCY LOWERING


### STURDY ARMRESTS

### CARE RAISER™

Enables you to gently lift your patient or resident's lower body with the touch of a button, allowing easy access for performing discreet hygiene tasks.



## SPACE REQUIREMENT



Required working area to facilitate hygiene and transfer activities while providing adequate access for the resident, equipment and assisting carer.



The Dementia Services Development Centre at the University of Stirling have undertaken a review of this product and have rated its design in relation to dementia design principles and its usability within a dementia-inclusive environment. For ease of identification, the suitability of the product is classified by number with corresponding explanation of use, within the ratings key provided. The DSDC Accredited Product logo is awarded on the merit of the individual product and colour-way. This accreditation does not merit the manufacture or final application of the product. Whilst every endeavour has made to ensure the ratings are reflective of the products suitability and applied use at the time of print, DSDC cannot be held responsible for the application of the final product, its performance or its interface with other products or finishes. For an environment to be considered 'dementia-inclusive' or 'dementia-friendly' careful consideration must be given to the specification of adjoining finishes, their performance specification (for example light reflectance value (LRV), slip resistance and use of pattern) and their suitability for their intended use.

## PRODUCT INFORMATION

Safe working load (SWL)	136 kg (300 lb)
Weight	64 kg (141 lb)
Min length, lowest position upright	875 mm (34 ½ in)
Max length, highest position reclined	1510 mm (59 ½ in)
Total width	640 mm (25 ¼ in)
Chassis, width	640 mm (25 ¼ in)
Chassis, width between legs	508 mm (20 in)
Chassis, length	865 mm (34 in)
Min height backrest, lowest position reclined	1020 mm (40 ¾ in)
Max height backrest, top position upright	1550 mm (61 in)
Min height armrests, lowest position reclined	775 mm (30 ½ in)
Max height armrests, top position upright	1320 mm (52 in)
Seat height, lowest/highest position	510 mm (20 ¼ in) / 1050 mm (41 ¾ in)
Lifting stroke (seat)	550 mm (21 ¾ in)

## ACCESSORIES

- Back pillow
- Neck pillow
- Seat cushion
- Bed pan and holder
- Safety belt
- Surface disinfectant
- Battery and battery charger

For more detailed information visit [www.arjo.com](http://www.arjo.com)