


Maxi Twin Compact










Maxi Twin Compact

- 1 DPS spreader bar
- 2 Jib
- 3 Handset
- 4 Emergency lowering
- 5 Emergency stop
- 6 Battery
- 7 Adjustable legs
- 8 Non-braked castors
- 9 Braked castors




Handset

- A Raise jib
- B Lower jib
- C Reclined position
- D Sitting position
- E Legs close
- F Legs open




Sling profiles that can be used with Maxi Twin Compact Size S - LL


Four point toileting
sling (with headrest)




Four point amputee sling



Four point padded sling



Four point unpadded sling



Four point mesh sling

