

Wellness Nordic Relax® Chair



Helping people with cognitive disorders relax through soothing music and tactile stimulation can have a beneficial effect on their wellbeing and quality of life.

The fully-automated Wellness Nordic Relax Chair combines these functions with a rocking motion to offer a unique and calming multi-sensory experience.

- Built-in audio system with MusiCure® music
- 3 pre-set automatic programmes offering multi-sensory stimulation
- Easy-to-clean, soft imitation leather
- Swing-out armrests
- Removable incontinence cover
- Wheeled pedestal for easy integration to Arjo patient lifters
- Optional: Heavy blanket and U-shaped cushion

Colours and options



SYSTEM SPECIFICATIONS

Power supply	AC 230V / 50Hz
Power consumption	90W
Operating environment	0 to 40 °C – 32 to 104 °F
Flame rating	EN1021-1 & 1021-2

MUSIC AND VIBRATION

Type	1-way, 1 speaker system
Impedance	8Ω
Speaker cone size	Full range, 700 mm x 2 – 27 ½" x 2
Speaker amplifier	14W (7W + 7W), 1 kHz
Vibration amplifier	14W, 70 Hz
Shielded design	No
External input terminal	Stereo mini jack, 3.5 mm

REMOTE CONTROL

Operation range	Approx. 1000 mm – 39 ½"
IR receiver	Both sides of chair
Batteries	2 x AAA

PRE-SET PROGRAMMES

Duration 20 min
 - Back rest lowered to semi-reclined position at the start of the programme
 - Semi-reclined position kept for the duration of the programme
 - Returning to seated position at the end of the programme
 - Automatic Musicure relax music
 - Tactile stimulation in the back

Relax

Duration 20 min
 - Backrest lowered in 5 steps
 - Returning to seated position at the end of the programme
 - Varying swing speed
 - Automatic Musicure stimulating music
 - Tactile stimulation in the back

Refresh

Duration 120 min
 - Back rest in upright position during the entire programme
 - Constant swing speed
 - No music

Comfort

OVERALL DIMENSIONS (INCLUDING PEDESTAL)

Length, upright position	1435 mm – 56 ½"
Height, upright position	1225 mm – 48 ¼"
Length, fully articulated	1830 mm – 72"
Height, fully articulated	835 mm – 32 ¾" (limited by armrest)
Max width (pedestal)	850 mm – 33 ¾"
Length, back rest	850 mm – 33 ¾"
Length, seat section	530 mm – 20 ¾"
Length, leg rest	500 mm – 19 ¾"
Width, seat section	665 mm – 26 ¼"
Seat height, back end	530 mm – 20 ¾"
Seat height, front end	585 mm – 23"
Product weight, incl pedestal	74.4 kg – 164 lbs
Product weight, chair only (incl. mattress, pillow, remote control)	46 – 102 lbs
Castor type	Swivelling, individually lockable
Castor size	100 mm – 4"
Floor clearance	130 mm – 5"
Safe Working Load	130 kg – 286 lbs

MATERIAL

Coated cloth	Polyvinyl chloride leather (synthetic leather)
Buffer material	Urethane
Armrest and incontinence cover	Stamskin TOP®, treated against bacteria, odour and fungi with Sanitized® PL 2432 (PES)

OPTIONAL ACCESSORIES AVAILABLE IN OTHER SIZES AND WEIGHTS

Weighted fibre blanket	3.7 kg – 8 lbs, 1250 mm x 1000 mm – 49 ¼" x 39 ½"
Hygiene cover	1250 mm x 1000 mm – 49 ¼" x 39 ½"
U-shaped cushion	Available in red

February 2018 Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. ® and ™ are trademarks belonging to the Arjo group of companies.
 © Arjo, July 2020

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure injury and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

Arjo AB · Hans Michelsengatan 10 · 211 20 Malmö · Sweden · +46 10 335 4500

www.arjo.com