

# SIZEWISE™: Sapphire Comfort Turn®

- 3 modes of therapy, static low air loss, pulsating low air loss & continuous lateral rotation with low air loss
- Full body lateral rotation
- Lockout to eliminate unintentional setting change
- Variable rotation cycle times from 3 to 90 minutes
- Fowler feature for seat inflation when head is elevated
- Domestic mattress featuring high strength air cells
- Turn and Hold feature
- Low pressure/power failure audible and visual alarms

## REAL THERAPY for REAL PEOPLE™

**The Sapphire Comfort Turn®** combines true low air loss with full body lateral rotation capabilities in a complete mattress replacement system.

The microprocessor controlled blower unit features up to a 40° “Full Turn” and “Right Turn” modes of rotation. These rotation modes have a minimum cycle time of 3 minutes and a maximum of 90 minutes. The “Turn & Hold” feature allows the caregiver to rotate the patient to the left or right to help administer the patient. Additional features include: “Static” for stationary low air loss therapy, “Fowler” to boost seat inflation when the head is elevated, “Auto Firm” for rapid inflation and easy patient transfer, “Lockout” to maintain caregiver settings without accidental interruptions, “Low Pressure” and “Power Failure” audible / visual alarms and “Self Diagnostic Check” to quickly assess all switches and sensors for proper calibration and operation.

The Sapphire Comfort Turn® is a domestic mattress with 20 interchangeable air cells that are incorporated with a side support bolster that cradles the patient when rotating. The entire system is contained in a durable woven base, and the polyurethane coated top cover is waterproof yet highly vapor permeable. The entire system is portable and fits most standard medical bed frames.

Air Cells



Control Unit



800.814.9389 | [sizewise.net](http://sizewise.net)

# SIZEWISE™: Sapphire Comfort Turn®

## CONTROL UNIT

Dimensions: ..... 53cm L x 28cm W x 40cm H  
Weight: ..... 5kg  
Cycle Time: ..... Variable (3-20 minutes)  
Electrical Rating: ..... 240V / 50Hz  
Power Cord: ..... Harmonized Hospital Grade Plug  
Warranty: ..... 1 Year

## SUPPORT SURFACE

Weight Capacity:  
- 90cm x200cm - Static 204kg - Therapy Mode 155kg (rotating)  
- 100cm x200cm - Static 450kg - Therapy Mode 295kg (rotating)  
- 120cm x200cm - Static 450kg - Therapy Mode 295kg (rotating)  
Warranty: ..... 2 Year Limited  
Top Cover Warranty: .....90 Days

**MADE  
IN THE USA**



800.814.9389 | [sizewise.net](http://sizewise.net)