

# Seba: Supine to Seated Edge of Bed **positioning solution**

ARJO HUNTLEIGH  
GETINGE GROUP



reddot design award



## Now a **safe way to transfer** from a **supine to seated position** at the edge of the bed

Only Seba™ empowers caregivers with the ability to safely get patients from a **Supine to Seated Edge of Bed (SSEB)** position consistently and ergonomically, to reduce the lift load and enable a healthy working posture – all while promoting patient comfort and safety.

### **Ergonomically designed**

- Reduces the lift load
- Enhances caregiver performance
- Designed to be used by one person

### **Consistency**

- Delivers predictability and consistency in day-to-day SSEB tasks

### **Compliance**

- Promotes patient safety and participation
- Helps maintain activities of daily living
- Facilitates patient mobility



...with people in mind

# Seba: seamless efficiency in SSEB positioning



#### Positioning Grip

- Primary control for lifting and lowering
- Maximized patient stability
- Reduction of primary load during lifting and lowering



#### Shoulder Support

- Thin profile for easy positioning
- Promotes patient comfort via optimized weight distribution
- Designed to help keep patient's spine straight throughout lift



#### Comfort Handle

- Delivers patient reassurance
- Promotes patient participation and security



#### Power Grip

- Ergonomic efficiency during patient lift
- Rotates for use on both sides of bed



#### Pivot Point

- Easy and controlled lifting and maneuvering



#### Leg Support

- Keeps lower body free of stress and strain
- Designed to help keep the patient's spine straight throughout lift



Maximum patient weight	Weight	Dimensions H x W x D
408 lbs (185 kg)	9.9 lbs (4.5 kg)	48.5" x 19.25" x 6" (1238 mm x 490 mm x 151 mm)

**ARJOHUNTLEIGH**  
GETINGE GROUP

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

© and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2015

ArjoHuntleigh AB, Hans Michelsengatan 10, SE-211 20 Malmö, Sweden  
Phone: +46 (0) 10 335 45 00

[www.ArjoHuntleigh.com](http://www.ArjoHuntleigh.com)

ArjoHuntleigh Inc.  
2349 West Lake Street, Addison, Illinois 60101  
Phone: 800-323-1245, Fax: 888-594-2756

ArjoHuntleigh Canada Inc.  
90 Matheson Blvd West, Suite 300, Mississauga  
Ontario L5R 3R3, CANADA  
Phone: 800-665-4831, [Info.Canada@ArjoHuntleigh.com](mailto:Info.Canada@ArjoHuntleigh.com)